


MERCK

**pharmaceutical
secondary standards**

Impurities & Related Compounds

To meet your needs in Pharmaceutical testing



Supelco®
Analytical Products

As a worldwide leading supplier to the global Life Science Industry, we are a key provider of standards for pharmaceutical quality testing. We maintain a comprehensive inventory to support your laboratory, and can provide a stress-free, time-saving and efficient route to procuring the standards you need.

- Secondary Certified Reference Materials – a convenient alternative to preparing in-house working standards
- TraceCERT® Elemental impurity standards for testing in accordance with testing limits per ICHQ3D Guidelines
- Certipur® pH Buffer Certified Reference Materials -ready to use and in convenient formats
- Extractables & Leachables – CRM Mixes and individual reference materials
- Custom services – customized certified reference materials & packaging services

[SigmaAldrich.com/standards-pharma](https://sigmaaldrich.com/standards-pharma)



Pharma Secondary CRMs

Pharmaceutical Secondary Certified Reference Materials (CRMs) provide pharmaceutical laboratories and manufacturers with a convenient and cost-effective alternative to making in-house working standards.

Using these commercially available Secondary CRMs will reduce the need to spend time and resources preparing in-house working standards.

Important Features of Secondary CRMs:

- **Traceability** to United States Pharmacopeia (USP), European Pharmacopoeia (EP) and British Pharmacopoeia (BP) Primary Reference Standards if available.
- **Certified** purity value according to ISO 17034 utilizing mass balance and/or qNMR approaches
- **Comprehensive** certificate according to ISO Guide 31

Our Secondary Standards are Certified Reference Materials (CRMs) manufactured under the scope of ISO 17034 and tested in an ISO/IEC 17025 accredited testing lab. Comprehensive Certificates of Analysis demonstrate traceability to USP, EP and/or BP primary standards (where available). In addition, an independent certified purity value is provided. This allows these Secondary CRMs to be used as reference materials for quantitative uses in a variety of applications.

The US Food and Drug Administration (FDA) and the European Pharmacopoeia (EP) both recognize the use of secondary standards or working standards which are established with reference to the corresponding primary standard (references available on request).

Comprehensive Certificate

The certificate of analysis for Pharmaceutical Secondary CRMs typically contains:

- Certified purity value by mass balance or qNMR (according to ISO 17034)
- Traceability results versus Pharmacopeial primary Reference Standards (see the example certificate extract shown on the next page)
- Qualification data along with comprehensive details on the analytical methods used
- Detailed handling and storage instructions

You can view and download an example COA for any pharmaceutical secondary standard at **[SigmaAldrich.com/pharmastandards](https://www.sigmaaldrich.com/pharmastandards)**. Select your standard of interest, scroll down to the "Documentation" section, and click "View Sample COA."

Download the latest version of the certificate for your standard

The values on the certificate are always traceable to the current Pharmacopeial lots (where available). If a valid Pharmacopeial lot changes, our corresponding Pharmaceutical Secondary CRM will be re-certified with traceability to the new lot. Then, a new certificate will be made available on our website. For this reason, the valid certificate should always be downloaded from the website prior to use of the material. To do this, find your product at **[SigmaAldrich.com](https://www.sigmaaldrich.com)** using the part number or via our standards explorer (located at **[SigmaAldrich.com/reference-materials](https://www.sigmaaldrich.com/reference-materials)**; bottom of the page) and scroll down to the Documentation section. Here you can enter the lot number of your standard.

Our portfolio of pharmaceutical secondary standard CRMs is continually growing. See what's new at:

[SigmaAldrich.com/newsecondaries](https://www.sigmaaldrich.com/newsecondaries)

and, for the full portfolio, visit

[SigmaAldrich.com/standards-pharma](https://www.sigmaaldrich.com/standards-pharma)

Comprehensive Certificate

Example COA: Ibuprofen, PHR1004

Traceability Assay:

Comparative assay demonstrates direct traceability to Pharmacopeial Standards

ASSAY vs. BP CRS (as is basis)

ASSAY VALUE

100.08 %

vs. BP BATCH

4105

Labeled Content = 99.9 %

Method: HPLC (ref.: Ibuprofen, Current Compendial Monographs)

Column: Ascentis C18, 250 mm x 4.6mm, 5µm particle size

Mobile Phase A: (10 g/L Chloroacetic Acid in Water (pH 3.0)): (Acetonitrile)

Mobile Phase Ratio: 40:60

Flow Rate: 1.5 mL/min

Column Temperature: 30 °C

Injection Volume: 10 µL

Detector: DAD, Wavelength: 254 nm

Internal Standard: Valerophenone

ASSAY vs. USP REFERENCE STANDARD (1335508) (as is basis)

ASSAY VALUE

100.24 %

vs. USP LOT

R13060

Labeled Content = 0.998 mg/mg

ASSAY vs. EP CRS (10020000) (as is basis)

ASSAY VALUE

100.24 %

vs. EP BATCH

6.0

Labeled Content = None

Assigned Content = 99.7% *

*The assigned content of the EP CRS was determined by assay against the USP Reference Standard.

Method: HPLC (ref.: Ibuprofen, Current Compendial Monographs)

Column: Ascentis C18 4.6 x 250 mm, 5 µm particle size

Mobile Phase: 10 g/L Chloroacetic Acid in Water (pH 3.0), Acetonitrile (40:60)

Flow Rate: 1.5 mL/min

Column Temperature: 30 °C


Injection Volume: 10 µL

Detector: DAD, Wavelength: 254 nm

Internal Standard: Valerophenone

Secondary Standards for Pharmaceutical Quality Control

Save time, money, and effort with Supelco® Pharmaceutical Secondary CRMs as an alternative to in-house working standards.



Our up-to-date portfolio contains >2500 pharmaceutical secondary CRMs and continues to expand. Visit our website often to view new products.

SigmaAldrich.com/pharmaimpurities

Monthly COA Updates

Did you know that Supelco® Pharmaceutical Secondary CRMs are continually updated, always providing traceability to the current Primary Reference Standard Batch?

You can be notified in advance of any COA changes that may occur due to re-qualification etc.

This service helps ensure that you always have the most current version of a product certificate.

Sign up for these email updates on our website:

SigmaAldrich.com/standards-pharma

Supelco www.sigmaaldrich.com

Certificate of Analysis – Certified Reference Material
METFORMIN HYDROCHLORIDE

Product no.: PHR1004-510102
Lot no.: LRAD2234
Description of CRM: White Crystals
Expiry date: 31 December 2026
Storage: Refrigerator, keep away from alkali.
Certificate version: LRAD2234.01 (Note: Certificate may be updated due to Pharmaceutical Q.C. changes or the availability of new lots. Check for latest changes at the website of our user center www.sigmaaldrich.com for details.)
Chemical formula: $C_4H_8N_2O_2 \cdot HCl$
Molecular mass: 146.6
CAS No.: 1115-70-4

Analyte | **Certified Purity ± associated uncertainty U, $u(k=2)$ (Mass Balance/Basis)**

| Metformin HCl | 100.0 % UCRM ± 0.1 %, k = 2.0 (as in basis) |

Metrological traceability: Traceable to the SI and higher order standards from NIST through an unbroken chain of comparisons. Additional traceability to Primary Standards is established through comparative assay determinations. See 'Details on metrological traceability' on page 2.

Measurement method: Where applicable, the certified value is based on a purity determination by mass balance. See 'Certification process details' on page 3.

Intended use: Intended for R&D and analytical use only. Not for drug, household or other uses.

Minimum sample size: 10 mg

Instructions for handling and correct use: Do not dry, use on the as is basis. The internal pressure of the container may be slightly different from the atmospheric pressure at the user's location. Clean, steady and carefully to avoid dispersion of the material. Attachment of a 20 mm aluminum ring and recommended for sealed portions.

Health and safety information: All chemical reference materials should be considered potentially hazardous and should be used only by qualified laboratory personnel. Please refer to the Safety Data Sheet for detailed information about the nature of any hazard and appropriate precautions to be taken.

Accreditation: Sigma-Aldrich RTC is accredited by the US accreditation authority ANAB as a registered reference material producer AB-1470 in accordance with ISO 17034.

Certificate issue date: 16 November 2022

ANAB
REG. 17034
AB-1470

[Signature]
[Andy O'Brien, Quality Control]

[Signature]
[Steven Miller, Quality Assurance]

Sigma-Aldrich RTC, 1933 Tuller-Savage Rd., Lenexa, KS 66215, USA.
Tel.: +1 800 441 2420 / Fax: +1 913 441 2421, www.sigmaaldrich.com
Sigma-Aldrich RTC is a subsidiary of Merck KGaA, Darmstadt, Germany.

Certificate Page: 1 of 3 Certificate version 01

Testing for impurities?

Our growing portfolio of Pharmaceutical Secondary CRMs includes both APIs and associated impurities.

SigmaAldrich.com/pharmaimpurities

**NOW
AVAILABLE!**

Human Insulin Secondary CRM

Worldwide, diabetes is one of the most common health conditions requiring treatment. The discovery of a therapy using insulin extracted from bovine pancreas was awarded the Nobel Prize in medicine in 1923. Today, insulin and several analogs are produced recombinantly, and save lives globally every day.

The new human insulin Pharmaceutical Secondary CRM, PHR8925, is traceable to the current USP reference standard and can be used as a working standard to determine identity or potency.



Cat. No.	Description	Pack Size	CAS No.
ABACAVIR			
PHR1256-500MG	ABACAVIR SULFATE	500MG	188062-50-2
PHR2073-50MG	ABACAVIR IMPURITY E	50MG	208762-35-0
PHR2074-50MG	ABACAVIR IMPURITY F	50MG	1443421-68-8
PHR2064-50MG	ABACAVIR RELATED COMPOUND A	50MG	124752-25-6
PHR2065-50MG	ABACAVIR RELATED COMPOUND B	50MG	124752-25-6
PHR2066-50MG	ABACAVIR RELATED COMPOUND C	50MG	172015-79-1
PHR2067-50MG	ABACAVIR RELATED COMPOUND D	50MG	1443421-69-9
ABIRATERONE ACETATE			
PHR2410-500MG	ABIRATERONE ACETATE	500MG	154229-18-2
ACAMPROSATE			
PHR2821-400MG	ACAMPROSATE CALCIUM	400MG	77337-73-6
PHR2822-50MG	ACAMPROSATE RELATED COMPOUND A	50MG	3687-18-1
ACEBUTOLOL			
PHR1907-250MG	ACEBUTOLOL HYDROCHLORIDE	250MG	34381-68-5
PHR2824-50MG	ACEBUTOLOL RELATED COMPOUND A	50MG	40188-45-2
ACETAMINOPHEN			
PHR1005-1G	ACETAMINOPHEN	1G	103-90-2
PHR1005-10G	ACETAMINOPHEN	10G	103-90-2
PHR1612-50MG	ACETAMINOPHEN RELATED COMPOUND A	50MG	2623-33-8
PHR1642-50MG	ACETAMINOPHEN RELATED COMPOUND B	50MG	1693-37-4
PHR1531-200MG	ACETAMINOPHEN RELATED COMPOUND C	200MG	614-80-2
PHR1530-200MG	ACETAMINOPHEN RELATED COMPOUND D	200MG	103-84-4
PHR1150-1G	ACETAMINOPHEN RELATED COMPOUND F	1G	100-02-7
PHR1972-100MG	ACETAMINOPHEN RELATED COMPOUND G	100MG	
PHR1973-100MG	ACETAMINOPHEN RELATED COMPOUND I	100MG	118-93-4
PHR1149-1G	ACETAMINOPHEN RELATED COMPOUND J	1G	539-03-7
PHR1148-1G	ACETAMINOPHEN RELATED COMPOUND K	1G	123-30-8
PHR1970-100MG	ACETAMINOPHEN IMPURITY C	100MG	
PHR1971-100MG	ACETAMINOPHEN IMPURITY E	100MG	99-93-4
ACETAZOLAMIDE			
PHR1908-5G	ACETAZOLAMIDE	5G	59-66-5
PHR2826-50MG	ACETAZOLAMIDE RELATED COMPOUND D	50MG	14949-00-9
PHR2919-100MG	ACITRETIN RELATED COMPOUND B	100MG	54350-48-0
ACYCLOVIR			
PHR1254-1G	ACYCLOVIR	1G	59277-89-3
PHR3451-100MG	ACYCLOVIR RELATED COMPOUND A	100MG	102728-64-3
PHR3452-100MG	ACYCLOVIR RELATED COMPOUND F	100MG	110104-37-5
PHR3453-100MG	ACYCLOVIR RELATED COMPOUND G	100MG	75128-73-3
ADAPALENE			
PHR1388-500MG	ADAPALENE	500MG	106685-40-9
PHR2827-100MG	ADAPALENE RELATED COMPOUND A	100MG	33626-98-1
PHR2828-100MG	ADAPALENE RELATED COMPOUND B	100MG	106685-41-0
PHR2829-50MG	ADAPALENE RELATED COMPOUND C	50MG	43109-77-9
PHR2830-50MG	ADAPALENE RELATED COMPOUND D	50MG	932033-57-3
PHR2935-50MG	ADAPALENE RELATED COMPOUND E	50MG	932033-58-4
ADEMETIONINE			
PHR3168-1G	ADEMETIONINE DISULFATE TOSYLATE	1G	97540-22-2
ADENOSINE			
PHR1138-1G	ADENOSINE	1G	58-61-7
PHR2244-100MG	ADENOSINE IMPURITY F	100MG	58-96-8
PHR2412-100MG	ADENOSINE IMPURITY G	100MG	58-63-9
PHR2411-100MG	ADENOSINE IMPURITY H, (GUANOSINE)	100MG	118-00-3
ALBUTEROL			
PHR3454-500MG	ALBUTEROL	500MG	18559-94-9
PHR1967-30MG	ALBUTEROL RELATED COMPOUND A	30MG	

Cat. No.	Description	Pack Size	CAS No.
ALBUTEROL			
PHR2040-50MG	ALBUTEROL RELATED COMPOUND B	50MG	41489-89-8
PHR1053-1G	ALBUTEROL SULFATE	1G	51022-70-9
ALLOPURINOL			
PHR1377-1G	ALLOPURINOL	1G	315-30-0
PHR2005-50MG	ALLOPURINOL RELATED COMPOUND A	50MG	27511-79-1
PHR2006-50MG	ALLOPURINOL RELATED COMPOUND B	50MG	22407-20-1
PHR2007-50MG	ALLOPURINOL RELATED COMPOUND C	50MG	1346604-13-4
PHR2008-50MG	ALLOPURINOL RELATED COMPOUND D	50MG	6994-25-8
PHR2009-50MG	ALLOPURINOL RELATED COMPOUND E	50MG	31055-19-3
PHR1981-50MG	ALLOPURINOL RELATED COMPOUND F	50MG	321571-07-7
ALMOTRIPTAN			
PHR3515-50MG	ALMOTRIPTAN RELATED COMPOUND A	50MG	1018676-02-2
PHR3514-50MG	ALMOTRIPTAN RELATED COMPOUND D	50MG	603137-43-5
ALOSETRON			
PHR2990-500MG	ALOSETRON HYDROCHLORIDE	500MG	122852-69-1
PHR2991-50MG	ALOSETRON RELATED COMPOUND A	50MG	122852-75-9
AMANTADINE			
PHR1711-500MG	AMANTADINE HYDROCHLORIDE	500MG	665-66-7
PHR2831-50MG	AMANTADINE RELATED COMPOUND A	50MG	935-56-8
PHR2832-50MG	AMANTADINE RELATED COMPOUND B	50MG	880-52-4
AMIFOSTINE			
PHR2940-300MG	AMIFOSTINE	300MG	112901-68-5
AMIKACIN			
PHR1654-1G	AMIKACIN	1G	37517-28-5
PHR1860-1G	AMIKACIN SULFATE	1G	39831-55-5
AMILORIDE			
PHR1839-1G	AMILORIDE HYDROCHLORIDE	1G	17440-83-4
PHR2061-50MG	AMILORIDE RELATED COMPOUND A	50MG	1458-01-1
AMINOLEVULINIC ACID			
PHR2993-500MG	AMINOLEVULINIC ACID HYDROCHLORIDE	500MG	
PHR2995-50MG	AMINOLEVULINIC ACID RELATED COMPOUND A	50MG	77479-02-8
PHR2994-50MG	AMINOLEVULINIC ACID RELATED COMPOUND B	50MG	109258-71-1
AMIODARONE			
PHR1164-1G	AMIODARONE HYDROCHLORIDE	1G	19774-82-4
PHR2415-100MG	AMIODARONE RELATED COMPOUND D	100MG	1951-26-4
PHR2416-100MG	AMIODARONE RELATED COMPOUND E	100MG	52490-15-0
AMITRAZ			
PHR2996-100MG	AMITRAZ RELATED COMPOUND A	100MG	60397-77-5
PHR2997-100MG	AMITRAZ RELATED COMPOUND B	100MG	33089-74-6
PHR2998-100MG	AMITRAZ RELATED COMPOUND C	100MG	16596-04-6
AMITRIPTYLINE			
PHR1384-1G	AMITRIPTYLINE HYDROCHLORIDE	1G	549-18-8
PHR2836-50MG	AMITRIPTYLINE RELATED COMPOUND A	50MG	1210-35-1
AMLODIPINE			
PHR1185-1G	AMLODIPINE BESYLATE	1G	111470-99-6
PHR2051-50MG	AMLODIPINE IMPURITY B	50MG	721958-72-1
PHR2050-50MG	AMLODIPINE RELATED COMPOUND A	50MG	2138811-33-1
PHR2052-50MG	AMLODIPINE RELATED COMPOUND C	50MG	43067-01-2
PHR2053-50MG	AMLODIPINE RELATED COMPOUND D	50MG	88150-62-3
PHR2056-50MG	AMLODIPINE RELATED COMPOUND K	50MG	1621516-91-3
AMOXAPINE			
PHR1841-500MG	AMOXAPINE	500MG	14028-44-5
PHR2837-100MG	AMOXAPINE RELATED COMPOUND D	100MG	3158-91-6
PHR3597-100MG	AMOXAPINE RELATED COMPOUND G	100MG	

Cat. No.	Description	Pack Size	CAS No.
AMOXICILLIN			
PHR1127-1G	AMOXICILLIN TRIHYDRATE	1G	61336-70-7
PHR2191-50MG	AMOXICILLIN IMPURITY D	50MG	1204829-51-5
PHR2195-100MG	AMOXICILLIN IMPURITY H	100MG	205826-86-4
PHR2188-50MG	AMOXICILLIN RELATED COMPOUND A	50MG	
PHR2190-50MG	AMOXICILLIN RELATED COMPOUND C	50MG	
PHR2192-50MG	AMOXICILLIN RELATED COMPOUND E	50MG	178738-48-2
PHR2193-100MG	AMOXICILLIN RELATED COMPOUND F	100MG	126247-63-0
PHR2196-50MG	AMOXICILLIN RELATED COMPOUND I	50MG	
AMPICILLIN			
PHR2838-500MG	AMPICILLIN	500MG	69-53-4
PHR1424-1G	AMPICILLIN SODIUM	1G	69-52-3
PHR1393-1G	AMPICILLIN TRIHYDRATE	1G	7177-48-2
AMRINONE			
PHR3466-1G	AMRINONE	1G	60719-84-8
PHR3468-200MG	AMRINONE RELATED COMPOUND A	200MG	62749-46-6
PHR3469-200MG	AMRINONE RELATED COMPOUND B	200MG	
PHR3467-100MG	AMRINONE RELATED COMPOUND C	100MG	62749-26-2
ANAGRELIDE			
PHR3545-50MG	ANAGRELIDE RELATED COMPOUND A	50MG	85325-12-8
PHR3547-50MG	ANAGRELIDE RELATED COMPOUND C	50MG	70381-75-8
ANASTROZOLE			
PHR1783-500MG	ANASTROZOLE	500MG	120511-73-1
ANTAZOLINE			
PHR3455-500MG	ANTAZOLINE PHOSPHATE	500MG	154-68-7
PHR2999-40MG	ANTAZOLINE RELATED COMPOUND A	40MG	26953-37-7
ANTIPYRINE			
PHR2347-500MG	ANTIPYRINE	500MG	60-80-0
PHR2348-50MG	ANTIPYRINE RELATED COMPOUND A	50MG	
ARIPIRAZOLE			
PHR1784-500MG	ARIPIRAZOLE	500MG	129722-12-9
PHR2320-100MG	ARIPIRAZOLE RELATED COMPOUND B	100MG	889443-20-3
PHR2321-100MG	ARIPIRAZOLE RELATED COMPOUND C	100MG	119532-26-2
PHR2322-100MG	ARIPIRAZOLE RELATED COMPOUND F	100MG	573691-09-5
PHR2323-100MG	ARIPIRAZOLE RELATED COMPOUND G	100MG	129722-25-4
PHR2324-50MG	ARIPIRAZOLE RELATED COMPOUND H	50MG	
ARTEMETHER			
PHR3423-200MG	ARTEMETHER	200MG	71963-77-4
PHR3000-40MG	ARTEMETHER RELATED COMPOUND B	40MG	71939-51-0
ARTICAINE			
PHR1882-100MG	ARTICAINE	100MG	23964-58-1
PHR1881-500MG	ARTICAINE HYDROCHLORIDE	500MG	23964-57-0
PHR2114-100MG	ARTICAINE RELATED COMPOUND A	100MG	1712677-79-6
PHR2115-50MG	ARTICAINE RELATED COMPOUND B	50MG	114176-52-2
PHR2116-100MG	ARTICAINE RELATED COMPOUND E	100MG	
ASCORBIC ACID			
PHR1008-10G	ASCORBIC ACID	10G	50-81-7
PHR1008-2G	ASCORBIC ACID	2G	50-81-7
PHR1708-30MG	ASCORBIC ACID IMP C	30MG	21675-47-8
PHR1709-30MG	ASCORBIC ACID IMP D	30MG	
PHR2418-100MG	ASCORBIC ACID IMPURITY F	100MG	89-65-6
PHR1433-1G	ASCORBYL GLUCOSIDE	1G	129499-78-1
PHR1455-2G	ASCORBYL PALMITATE	2G	137-66-6
ASPARAGINE			
PHR2349-300MG	ASPARAGINE ANHYDROUS	300MG	70-47-3
PHR2350-300MG	ASPARAGINE MONOHYDRATE	300MG	

Cat. No.	Description	Pack Size	CAS No.
ASPARTAME			
PHR1381-1G	ASPARTAME	1G	22839-47-0
PHR2419-300MG	ASPARTAME ACESULFAME	300MG	106372-55-8
PHR9041-50MG	ASPARTAME ACID (ASP-PHE)	50MG	13433-09-5
PHR2420-50MG	ASPARTAME RELATED COMPOUND A	50MG	5262-10-2
ASPIRIN			
PHR1003-1G	ASPIRIN (ACETYL SALICYLIC ACID)	1G	50-78-2
PHR1942-50MG	ACETYLSALICYLIC ANHYDRIDE(ASPIRIN IMPURITY F)	50MG	1466-82-6
PHR1941-200MG	SALSALATE (ASPIRIN IMPURITY E)	200MG	552-94-3
ATAZANAVIR			
PHR3002-500MG	ATAZANAVIR SULFATE	500MG	229975-97-7
PHR3001-50MG	ATAZANAVIR RELATED COMPOUND A	50MG	162537-11-3
A-TOCOPHEROL			
PHR1031-500MG	A-TOCOPHEROL	500MG	10191-41-0
PHR1030-500MG	A-TOCOPHERYL ACETATE	500MG	7695-91-2
PHR1029-500MG	A-TOCOPHERYL ACID SUCCINATE	500MG	4345-03-3
ATOMOXETINE			
PHR1679-500MG	ATOMOXETINE HYDROCHLORIDE	500MG	82248-59-7
PHR2840-25MG	ATOMOXETINE RELATED COMPOUND A	25MG	42142-52-9
PHR3288-30MG	ATOMOXETINE RELATED COMPOUND B	30MG	873310-28-2
PHR3460-40MG	ATOMOXETINE RELATED COMPOUND C	40MG	873310-31-7
PHR3400-50MG	ATOMOXETINE S-ISOMER	50MG	82857-39-4
ATORVASTATIN			
PHR1422-1G	ATORVASTATIN CALCIUM	1G	344423-98-9
PHR3858-30MG	ATORVASTATIN RELATED COMPOUND A	30MG	433289-83-9
PHR1869-30MG	ATORVASTATIN RELATED COMPOUND B	30MG	887196-25-0
PHR1870-30MG	ATORVASTATIN RELATED COMPOUND C	30MG	693793-53-2
PHR1871-30MG	ATORVASTATIN RELATED COMPOUND D	30MG	148146-51-4
PHR1872-30MG	ATORVASTATIN RELATED COMPOUND E	30MG	1105067-88-6
PHR1873-30MG	ATORVASTATIN RELATED COMPOUND H	30MG	125995-03-1
PHR1874-30MG	ATORVASTATIN RELATED COMPOUND I	30MG	125971-95-1
ATOVAQUONE			
PHR1591-500MG	ATOVAQUONE	500MG	95233-18-4
PHR3792-50MG	ATOVAQUONE RELATED COMPOUND A	50MG	137732-39-9
AZACITIDINE			
PHR1911-300MG	AZACITIDINE	300MG	320-67-2
PHR3473-50MG	AZACITIDINE RELATED COMPOUND A	50MG	931-86-2
PHR3472-30MG	AZACITIDINE RELATED COMPOUND B	30MG	6974-32-9
AZATHIOPRINE			
PHR1282-1G	AZATHIOPRINE	1G	446-86-6
PHR2841-30MG	AZATHIOPRINE RELATED COMPOUND A	30MG	4531-54-8
AZELASTINE			
PHR1636-1G	AZELASTINE HYDROCHLORIDE	1G	79307-93-0
PHR2842-50MG	AZELASTINE RELATED COMPOUND D	50MG	53242-88-9
PHR2843-50MG	AZELASTINE RELATED COMPOUND E	50MG	20526-97-0
AZITHROMYCIN			
PHR1088-1G	AZITHROMYCIN	1G	117772-70-0
PHR2782-30MG	AZITHROMYCIN N-OXIDE	30MG	90503-06-3
PHR2781-50MG	AZITHROMYCIN RELATED COMPOUND F	50MG	612069-28-0
PHR3403-30MG	AZITHROMYCIN RELATED COMPOUND H	30MG	612069-30-4
PHR2625-30MG	DESOSAMINYLAZITHROMYCIN	30MG	117693-41-1
PHR3759-50MG	N-DEMETHYLAZITHROMYCIN	50MG	172617-84-4
BACLOFEN			
PHR1682-1G	BACLOFEN	1G	1134-47-0
PHR8947-100MG	BACLOFEN RELATED COMPOUND A	100MG	22518-27-0

Cat. No.	Description	Pack Size	CAS No.
BALSALAZIDE			
PHR2929-100MG	BALSALAZIDE RELATED COMPOUND A	100MG	
PHR3272-50MG	BALSALAZIDE RELATED COMPOUND B	50MG	1798395-96-6
BENAZEPRIL			
PHR1912-250MG	BENAZEPRIL HYDROCHLORIDE	250MG	86541-74-4
PHR8936-50MG	BENAZEPRIL RELATED COMPOUND A	50MG	215447-89-5
PHR8937-30MG	BENAZEPRIL RELATED COMPOUND B	30MG	86541-77-7
PHR3283-75MG	BENAZEPRIL RELATED COMPOUND C	75MG	86541-78-8
PHR2922-50MG	BENAZEPRIL RELATED COMPOUND D	50MG	
PHR2422-50MG	BENAZEPRIL RELATED COMPOUND E	50MG	88372-47-8
PHR2423-100MG	BENAZEPRIL RELATED COMPOUND F	100MG	109010-60-8
PHR2930-75MG	BENAZEPRIL RELATED COMPOUND G	75MG	103129-58-4
BENDAMUSTINE			
PHR2333-500MG	BENDAMUSTINE HYDROCHLORIDE MONOHYDRATE	500MG	
PHR2944-50MG	BENDAMUSTINE RELATED COMPOUND A	50MG	109882-29-3
PHR2945-40MG	BENDAMUSTINE RELATED COMPOUND B	40MG	
PHR3292-40MG	BENDAMUSTINE RELATED COMPOUND C	40MG	3543-74-6
PHR2948-40MG	BENDAMUSTINE RELATED COMPOUND F	40MG	
PHR2951-50MG	BENDAMUSTINE RELATED COMPOUND I	50MG	87475-54-5
BENZOTHIADIAZINE			
PHR1192-500MG	BENZOTHIADIAZINE RELATED COMPOUND A	500MG	121-30-2
BENZOTROPINE			
PHR3517-100MG	BENZTROPINE RELATED COMPOUND A	100MG	25471-67-4
BETAMETHASONE			
PHR1398-1G	BETAMETHASONE	1G	378-44-9
PHR1779-1G	BETAMETHASONE ACETATE	1G	987-24-6
PHR1399-1G	BETAMETHASONE DIPROPIONATE	1G	5593-20-4
PHR3475-50MG	BETAMETHASONE DIPROPIONATE RELATED COMPOUND B	50MG	5534-13-4
PHR3476-50MG	BETAMETHASONE DIPROPIONATE RELATED COMPOUND C	50MG	75883-07-7
PHR2424-1G	BETAMETHASONE SODIUM PHOSPHATE	1G	
PHR1780-1G	BETAMETHASONE VALERATE	1G	2152-44-5
PHR3474-100MG	BETAMETHASONE VALERATE RELATED COMPOUND A	100MG	2240-28-0
BETAXOLOL			
PHR3386-50MG	BETAXOLOL RELATED COMPOUND A	100MG	67193-95-7
BICALUTAMIDE			
PHR1678-1G	BICALUTAMIDE	1G	90357-06-5
PHR3394-50MG	BICALUTAMIDE RELATED COMPOUND A	50MG	945419-64-7
PHR3392-40MG	BICALUTAMIDE RELATED COMPOUND B	40MG	1166228-30-3
BISACODYL			
PHR2794-250MG	BISACODYL	250MG	603-50-9
PHR3388-50MG	BISACODYL RELATED COMPOUND A	50MG	603-41-8
PHR3548-50MG	BISACODYL RELATED COMPOUND B	50MG	16985-05-0
PHR3385-50MG	BISACODYL RELATED COMPOUND C	50MG	72901-16-7
PHR3390-40MG	BISACODYL RELATED COMPOUND E	40MG	111664-35-8
BISOCTRIZOLE			
PHR1913-500MG	BISOCTRIZOLE	500MG	103597-45-1
PHR2844-50MG	BISOCTRIZOLE RELATED COMPOUND A	50MG	3147-75-9
BRIMONIDINE			
PHR2795-500MG	BRIMONIDINE TARTRATE	500MG	70359-46-5
PHR3480-100MG	BRIMONIDINE RELATED COMPOUND A	100MG	50358-63-9
PHR3481-100MG	BRIMONIDINE RELATED COMPOUND B	100MG	842138-74-3
PHR3479-30MG	BRIMONIDINE RELATED COMPOUND E	30MG	168329-48-4
BUDESONIDE			
PHR1178-500MG	BUDESONIDE	500MG	51333-22-3
PHR1822-50MG	BUDESONIDE RELATED COMPOUND G	50MG	137174-25-5
PHR1823-50MG	BUDESONIDE RELATED COMPOUND L	50MG	216453-74-6
PHR2574-300MG	BUDESONIDE ACETONIDE	300MG	638-94-8

Cat. No.	Description	Pack Size	CAS No.
BUMETANIDE			
PHR3393-40MG	BUMETANIDE RELATED COMPOUND A	40MG	28328-54-3
PHR3391-50MG	BUMETANIDE RELATED COMPOUND B	50MG	28328-53-2
BUPIVACAINE			
PHR1128-1G	BUPIVACAINE HYDROCHLORIDE	1G	73360-54-0
PHR1738-40MG	BUPIVACAINE RELATED COMPOUND A	40MG	
PHR1756-20MG	BUPIVACAINE RELATED COMPOUND B	20MG	15883-20-2
PHR9084-100MG	(S)-BUPIVACAINE	100MG	27262-47-1
BUPROPION			
PHR3178-500MG	BUPROPION HYDROBROMIDE	500MG	905818-69-1
PHR1730-1G	BUPROPION HYDROCHLORIDE	1G	31677-93-7
PHR1752-30MG	BUPROPION HYDROCHLORIDE RELATED COMPOUND A	30MG	1049718-72-0
PHR1753-30MG	BUPROPION HYDROCHLORIDE RELATED COMPOUND B	30MG	1049718-43-5
PHR1754-30MG	BUPROPION RELATED COMPOUND C	30MG	152943-33-4
PHR1755-30MG	BUPROPION RELATED COMPOUND F	30MG	857233-13-7
PHR3482-100MG	BUPROPION RELATED COMPOUND G	100MG	34911-51-8
BUSPIRONE			
PHR1917-500MG	BUSPIRONE HYDROCHLORIDE	500MG	33386-08-2
PHR2845-2X100MG	BUSPIRONE RELATED COMPOUND A	2X100MG	20980-22-7
PHR2846-100MG	BUSPIRONE RELATED COMPOUND G	100MG	84746-24-7
PHR2847-100MG	BUSPIRONE RELATED COMPOUND K	100MG	1075-89-4
PHR2848-100MG	BUSPIRONE RELATED COMPOUND N	100MG	257877-44-4
PHR3483-150MG	BUSPIRONE RELATED COMPOUND L	150MG	21098-11-3
CAFFEINE			
PHR1009-1G	CAFFEINE	1G	58-08-2
PHR2038-50MG	CAFFEINE IMPURITY B	50MG	7597-60-6
PHR1999-50MG	CAFFEINE IMPURITY C	50MG	519-32-4
PHR1951-30MG	CAFFEINE IMPURITY E NITRATE	30MG	
PHR2000-50MG	CAFFEINE IMPURITY F	50MG	611-59-6
PHR1095-1G	CAFFEINE MELTING POINT STANDARD	1G	58-08-2
CANDESARTAN			
PHR1854-200MG	CANDESARTAN CILEXETIL	200MG	145040-37-5
PHR2035-20MG	CANDESARTAN CILEXETIL IMPURITY C	20MG	
PHR2036-20MG	CANDESARTAN CILEXETIL IMPURITY E	20MG	914613-35-7
PHR2022-20MG	CANDESARTAN CILEXETIL IMPURITY I	20MG	139481-69-9
PHR2017-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND A	20MG	139481-58-6
PHR2018-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND B	20MG	869631-11-8
PHR2019-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND D	20MG	
PHR2020-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND F	20MG	
PHR2021-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND G	20MG	139481-59-7
CAPECITABINE			
PHR1405-1G	CAPECITABINE	1G	154361-50-9
PHR2849-50MG	CAPECITABINE RELATED COMPOUND A	50MG	66335-38-4
PHR2850-50MG	CAPECITABINE RELATED COMPOUND B	50MG	3094-09-5
CARBAMAZEPINE			
PHR1067-1G	CARBAMAZEPINE	1G	298-46-4
PHR2237-50MG	CARBAMAZEPINE IMPURITY B	50MG	611-64-3
PHR2330-50MG	CARBAMAZEPINE IMPURITY C	50MG	1219170-51-0
PHR2433-100MG	CARBAMAZEPINE IMPURITY F	100MG	33948-22-0
PHR1918-50MG	CARBAMAZEPINE RELATED COMPOUND A	50MG	3564-73-6
PHR1197-500MG	CARBAMAZEPINE RELATED COMPOUND B	500MG	256-96-2
PHR3293-50MG	10-METHOXYCARBAMAZEPINE	50MG	28721-09-7
CARBINOXAMINE			
PHR2802-400MG	CARBINOXAMINE MALEATE	400MG	3505-38-2
PHR3593-50MG	CARBINOXAMINE RELATED COMPOUND B	50MG	27652-89-7
PHR3387-50MG	CARBINOXAMINE RELATED COMPOUND C	50MG	

Cat. No.	Description	Pack Size	CAS No.
CARPROFEN			
PHR1452-1G	CARPROFEN	1G	53716-49-7
PHR1564-200MG	CARPROFEN RELATED COMPOUND A	200MG	86-74-8
CARVEDIOL			
PHR1265-1G	CARVEDIOL	1G	72956-09-3
PHR3762-50MG	CARVEDIOL RELATED COMPOUND A	50MG	1198090-73-1
PHR2952-40MG	CARVEDIOL RELATED COMPOUND B	40MG	918903-20-5
PHR2953-40MG	CARVEDIOL RELATED COMPOUND C	40MG	72955-94-3
PHR2434-50MG	CARVEDIOL RELATED COMPOUND D	50MG	51997-51-4
PHR2435-50MG	CARVEDIOL RELATED COMPOUND E	50MG	
CEFADROXIL			
PHR3579-400MG	CEFADROXIL	400MG	66592-87-8
PHR3004-50MG	CEFADROXIL RELATED COMPOUND B	50MG	22252-43-3
CEFAZOLIN			
PHR1291-500MG	CEFAZOLIN	500MG	25953-19-9
PHR2437-100MG	CEFAZOLIN IMPURITY E	100MG	29490-19-5
PHR2436-100MG	CEFAZOLIN IMPURITY H	100MG	957-68-6
PHR2681-100MG	CEFAZOLIN SODIUM SALT	100MG	27164-46-1
CEFDINIR			
PHR2665-500MG	CEFDINIR	500MG	91832-40-5
PHR2667-20MG	CEFDINIR RELATED COMPOUND B	20MG	79350-10-0
CEFEPIME			
PHR1763-1G	CEFEPIME HYDROCHLORIDE	1G	123171-59-5
PHR2851-50MG	CEFEPIME RELATED COMPOUND D	50MG	65872-41-5
CEFUROXIME			
PHR3296-700MG	CEFUROXIME AXETIL	700MG	64544-07-6
PHR2683-300MG	CEFUROXIME SODIUM	300MG	56238-63-2
CELECOXIB			
PHR1683-1G	CELECOXIB	1G	169590-42-5
PHR3485-50MG	CELECOXIB RELATED COMPOUND C	50MG	170570-09-9
PHR3484-50MG	CELECOXIB RELATED COMPOUND D	50MG	1061214-06-9
CETIRIZINE			
PHR1656-1G	CETIRIZINE HYDROCHLORIDE	1G	83881-52-1
PHR1710-40MG	CETIRIZINE HYDROCHLORIDE RELATED COMPOUND A	40MG	246870-46-2
PHR3280-40MG	CETIRIZINE RELATED COMPOUND B	40MG	108983-83-1
PHR3642-50MG	CETIRIZINE RELATED COMPOUND C	50MG	83881-37-2
PHR2852-100MG	CETIRIZINE RELATED COMPOUND G	100MG	164726-80-1
PHR3754-150MG	LEVOCETIRIZINE AMIDE	150MG	909779-33-5
PHR3195-500MG	LEVOCETIRIZINE DIHYDROCHLORIDE	500MG	130018-87-0
CHLORHEXIDINE			
PHR1421-1G	CHLORHEXIDINE	1G	55-56-1
PHR1222-500MG	CHLORHEXIDINE ACETATE	500MG	56-95-1
PHR1294-1G	CHLORHEXIDINE DIGLUCONATE SOLUTION	1G	18472-51-0
PHR1200-1G	CHLOROHEXIDINE IMP P (p-CHLOROANLINE)	1G	203-401-0
CHLOROQUINE			
PHR1258-1G	CHLOROQUINE PHOSPHATE	1G	50-63-5
PHR2438-50MG	CHLOROQUINE RELATED COMPOUND A	50MG	86-98-6
PHR3277-40MG	CHLOROQUINE RELATED COMPOUND D	40MG	1476-52-4
PHR3278-40MG	CHLOROQUINE RELATED COMPOUND E	40MG	
PHR2933-40MG	CHLOROQUINE RELATED COMPOUND G	40MG	
PHR1782-1G	HYDROXYCHLOROQUINE SULFATE	1G	747-36-4
CHLORPHENIRAMINE			
PHR1016-500MG	CHLORPHENIRAMINE MALEATE	500MG	113-92-8
PHR2105-100MG	CHLORPHENAMINE N-OXIDE MALEATE	100MG	120244-82-8
PHR2004-50MG	CHLORPHENIRAMINE IMPURITY A	50MG	1246816-57-8
PHR2103-100MG	CHLORPHENAMINE IMPURITY D	100MG	

Cat. No.	Description	Pack Size	CAS No.
CHLORPHENIRAMINE			
PHR1647-200MG	CHLORPHENIRAMINE RELATED COMPOUND B	200MG	1202-34-2
PHR1688-50MG	CHLORPHENIRAMINE RELATED COMPOUND C	50MG	22630-25-7
PHR2803-500MG	DEXCHLORPHENIRAMINE MALEATE	500MG	2438-32-6
CHLORTHALIDONE			
PHR3404-400MG	CHLORTHALIDONE	400MG	77-36-1
PHR2915-50MG	CHLORTHALIDONE RELATED COMPOUND A	50MG	5270-74-6
CHLORZOXAZONE			
PHR1919-500MG	CHLORZOXAZONE	500MG	95-25-0
PHR2699-100MG	CHLORZOXAZONE RELATED COMPOUND A	100MG	95-85-2
CHOLECALCIFEROL			
PHR1237-500MG	CHOLECALCIFEROL	500MG	67-97-0
PHR2439-100MG	CHOLECALCIFEROL IMPURITY B	100MG	434-16-2
PHR2663-100MG	VITAMIN D3 ISOMER MIX	100MG	
CICLOPIROX			
PHR1920-500MG	CICLOPIROX	500MG	29342-05-0
PHR2854-50MG	CICLOPIROX RELATED COMPOUND B	50MG	14818-35-0
PHR3646-50MG	CICLOPIROX RELATED COMPOUND C	50MG	67587-24-0
CIDOFOVIR			
PHR3583-50MG	CIDOFOVIR RELATED COMPOUND A HYDROCHLORIDE	50MG	1312776-50-3
PHR3616-50MG	CIDOFOVIR RELATED COMPOUND B	50MG	120362-35-8
CILOSTAZOL			
PHR1503-1G	CILOSTAZOL	1G	73963-72-1
PHR2855-100MG	CILOSTAZOL RELATED COMPOUND A	100MG	54197-66-9
PHR1676-100MG	CILOSTAZOL RELATED COMPOUND C	100MG	865792-18-3
CIMETIDINE			
PHR1075-1G	CIMETIDINE	1G	51481-61-9
PHR2250-50MG	CIMETIDINE IMPURITY B	50MG	73785-35-0
PHR2251-100MG	CIMETIDINE IMPURITY C	100MG	52568-80-6
PHR2242-25MG	CIMETIDINE IMPURITY D	25MG	59660-24-1
PHR2243-25MG	CIMETIDINE IMPURITY E	25MG	54237-72-8
PHR2258-50MG	CIMETIDINE IMPURITY F	50MG	
PHR2259-50MG	CIMETIDINE IMPURITY G	50MG	31857-31-5
PHR2255-50MG	CIMETIDINE IMPURITY H	50MG	74886-59-2
PHR2260-50MG	CIMETIDINE IMPURITY I	50MG	29636-87-1
PHR2252-100MG	CIMETIDINE IMPURITY J	100MG	38603-72-4
CINACALCET			
PHR3709-100MG	CINACALCET RELATED COMPOUND D	100MG	1217809-88-5
CIPROFLOXACIN			
PHR1167-1G	CIPROFLOXACIN	1G	85721-33-1
PHR2440-200MG	CIPROFLOXACIN FORMAMIDE	200MG	93594-39-9
PHR1044-1G	CIPROFLOXACIN HYDROCHLORIDE	1G	86393-32-0
PHR1174-500MG	CIPROFLOXACIN IMPURITY A (PHEUR	500MG	86393-33-1
PHR3731-50MG	CIPROFLOXACIN IMPURITY B (PH. EUR.)	50MG	93107-11-0
PHR1910-50MG	CIPROFLOXACIN RELATED COMPOUND A	50MG	526204-10-4
CITALOPRAM			
PHR1640-1G	CITALOPRAM HYDROBROMIDE	1G	59729-32-7
PHR1901-50MG	CITALOPRAM RELATED COMPOUND A	50MG	64372-56-1
PHR8899-50MG	CITALOPRAM RELATED COMPOUND B	50MG	1332724-03-4
PHR1902-50MG	CITALOPRAM RELATED COMPOUND C	50MG	
PHR1903-30MG	CITALOPRAM RELATED COMPOUND D	30MG	1188264-72-3
PHR1904-50MG	CITALOPRAM RELATED COMPOUND E	50MG	62498-71-9
PHR1905-25MG	CITALOPRAM RELATED COMPOUND G	25MG	1332724-08-9
CLARITHROMYCIN			
PHR1038-500MG	CLARITHROMYCIN	500MG	81103-11-9
PHR3705-40MG	CLARITHROMYCIN IMPURITY D	40MG	101666-68-6

Cat. No.	Description	Pack Size	CAS No.
CLARITHROMYCIN			
PHR1966-10MG	CLARITHROMYCIN IMPURITIES G AND O	10MG	
PHR2059-25MG	CLARITHROMYCIN IMPURITY K	25MG	127157-35-1
PHR1980-20MG	CLARITHROMYCIN IMPURITY N	20MG	
PHR1968-20MG	CLARITHROMYCIN IMPURITY Q	20MG	118074-07-0
CLENBUTEROL			
PHR3488-200MG	CLENBUTEROL HYDROCHLORIDE	200MG	21898-19-1
PHR3489-30MG	CLENBUTEROL RELATED COMPOUND B	30MG	37148-49-5
CLIDINIUM BROMIDE			
PHR3564-400MG	CLIDINIUM BROMIDE RELATED COMPOUND A	400MG	
CLINDAMYCIN			
PHR1159-1G	CLINDAMYCIN HYDROCHLORIDE	1G	58207-19-5
PHR8882-30MG	CLINDAMYCIN IMPURITY C	30MG	16684-06-3
PHR8883-50MG	CLINDAMYCIN IMPURITY D	50MG	26389-84-4
PHR1021-1G	CLINDAMYCIN PHOSPHATE	1G	24729-96-2
PHR8886-100MG	TRIACETYLCINDAMYCIN	100MG	
CLOMIPHENE			
PHR3742-150MG	CLOMIPHENE RELATED COMPOUND A	150MG	74056-26-1
CLOMIPRAMINE			
PHR1720-1G	CLOMIPRAMINE HYDROCHLORIDE	1G	17321-77-6
PHR3720-50MG	CLOMIPRAMINE RELATED COMPOUND A	50MG	2305824-51-3
PHR3721-50MG	CLOMIPRAMINE RELATED COMPOUND C	50MG	2095-42-3
CLONIDINE			
PHR3626-50MG	CLONIDINE RELATED COMPOUND A	50MG	54707-71-0
PHR3177-200MG	APRACLONIDINE HYDROCHLORIDE	200MG	73218-79-8
CLOPIDOGREL			
PHR1431-1G	CLOPIDOGREL BISULFATE	1G	120202-66-6
PHR1990-40MG	CLOPIDOGREL RELATED COMPOUND A	40MG	144750-42-5
PHR1991-40MG	CLOPIDOGREL RELATED COMPOUND B	40MG	2137493-39-9
PHR1992-40MG	CLOPIDOGREL RELATED COMPOUND C	40MG	120202-71-3
CLOTRIMAZOLE			
PHR1058-1G	CLOTRIMAZOLE	1G	23593-75-1
PHR3612-50MG	CLOTRIMAZOLE IMPURITY B	50MG	23593-71-7
PHR3005-50MG	CLOTRIMAZOLE IMPURITY F	50MG	15469-97-3
PHR2700-50MG	CLOTRIMAZOLE RELATED COMPOUND A	50MG	66774-02-5
CLOZAPINE			
PHR3655-200MG	CLOZAPINE	200MG	5786-21-0
CROMOLYN			
PHR2727-1G	CROMOLYN SODIUM	1G	15826-37-6
PHR3297-50MG	CROMOLYN RELATED COMPOUND A	50MG	16150-44-0
PHR3298-50MG	CROMOLYN RELATED COMPOUND B	50MG	16150-45-1
CYCLOBENZAPRINE			
PHR2797-400MG	CYCLOBENZAPRINE HYDROCHLORIDE	400MG	6202-23-9
PHR3299-30MG	CYCLOBENZAPRINE RELATED COMPOUND A	30MG	18029-54-4
CYCLOPHOSPHAMIDE			
PHR1404-1G	CYCLOPHOSPHAMIDE	1G	6055-19-2
PHR3619-50MG	CYCLOPHOSPHAMIDE RELATED COMPOUND B	50MG	158401-52-6
CYCLOMETHICONE			
PHR1565-500MG	D4 CYCLOMETHICONE	500MG	556-67-2
PHR1566-500MG	D5 CYCLOMETHICONE	500MG	541-02-6
CYPROHEPTADINE			
PHR2684-500MG	CYPROHEPTADINE HYDROCHLORIDE	500MG	41354-29-4
PHR2685-50MG	CYPROHEPTADINE RELATED COMPOUND A	50MG	256-81-5
CAMPHOR			
PHR1119-1G	D-CAMPHOR	1G	464-49-3
PHR2427-100MG	D-CAMPHOR IMPURITY A	100MG	80-56-8

Cat. No.	Description	Pack Size	CAS No.
CAMPHOR			
PHR2426-1G	CAMPHOR IMPURITY B, CAMPHENE	1G	79-92-5
PHR9058-100MG	D-CAMPHOR IMPURITY C (PHEUR)	100MG	19902-08-0
PHR2425-100MG	D-CAMPHOR IMPURITY D	100MG	7787-20-4
PHR2430-100MG	D-CAMPHOR IMPURITY F	100MG	2217-02-9
PHR2431-100MG	D-CAMPHOR IMPURITY I	100MG	124-76-5
PHR2429-100MG	D-CAMPHOR IMPURITY J	100MG	464-45-9
PHR2805-400MG	METHYL BENZYLIDENE CAMPHOR	400MG	36861-47-9
DANTROLENE			
PHR2291-500MG	DANTROLENE	500MG	7261-97-4
PHR2819-200MG	DANTROLENE SODIUM	200MG	24868-20-0
PHR2292-100MG	DANTROLENE RELATED COMPOUND A	100MG	
PHR2293-100MG	DANTROLENE RELATED COMPOUND B	100MG	
PHR2294-100MG	DANTROLENE RELATED COMPOUND C	100MG	
DACARBAZINE			
PHR1924-250MG	DACARBAZINE	250MG	4342-03-4
PHR2304-100MG	DACARBAZINE RELATED COMPOUND A	100MG	
DES Loratadine			
PHR1680-1G	DES Loratadine	1G	100643-71-8
PHR2314-50MG	DES Loratadine RELATED COMPOUND A	50MG	117796-50-6
PHR2315-50MG	DES Loratadine RELATED COMPOUND B	50MG	432543-89-0
PHR2316-100MG	DES Loratadine RELATED COMPOUND F	100MG	
DEXAMETHASONE			
PHR1526-1G	DEXAMETHASONE	1G	50-02-2
PHR1527-1G	DEXAMETHASONE ACETATE	1G	1177-87-3
PHR1628-500MG	DEXAMETHASONE PHOSPHATE	500MG	
PHR1768-1G	DEXAMETHASONE SODIUM PHOSPHATE	1G	2392-39-4
DEXTROMETHORPHAN			
PHR1018-500MG	DEXTROMETHORPHAN HYDROBROMIDE	500MG	6700-34-1
PHR3300-50MG	DEXTROMETHORPHAN IMPURITY A	50MG	1087-69-0
PHR1976-30MG	DEXTROMETHORPHAN IMPURITY D	30MG	1453167-99-1
PHR2206-100MG	DEXTROMETHORPHAN N-OXIDE HYDROCHLORIDE	100MG	1177419-85-0
DEXTROSE			
PHR1000-1G	DEXTROSE	1G	50-99-7
PHR2085-100MG	DEXTROSE IMPURITY B	100MG	499-40-1
PHR2086-1G	DEXTROSE IMPURITY C	1G	1109-28-0
DIATRIZOIC ACID			
PHR3303-200MG	DIATRIZOIC ACID DIHYDRATE	200MG	50978-11-5
PHR3302-100MG	DIATRIZOIC ACID RELATED COMPOUND A	100MG	1713-07-1
DICLOFENAC			
PHR2444-100MG	DICLOFENAC IMPURITY B	100MG	22121-58-0
PHR2445-100MG	DICLOFENAC IMPURITY C	100MG	27204-57-5
PHR2447-100MG	DICLOFENAC IMPURITY D	100MG	127792-23-8
PHR2443-100MG	DICLOFENAC IMPURITY E	100MG	59-48-3
PHR2446-100MG	DICLOFENAC IMPURITY F	100MG	560075-65-2
PHR2441-500MG	DICLOFENAC POTASSIUM	500MG	15307-81-0
PHR2442-100MG	DICLOFENAC RELATED COMPOUND A	100MG	15362-40-0
PHR1144-1G	DICLOFENAC SODIUM	1G	15307-79-6
DIETHYLTOLUAMIDE			
PHR3304-3G	DIETHYLTOLUAMIDE	3G	134-62-3
PHR3305-50MG	DIETHYLTOLUAMIDE RELATED COMPOUND A	50MG	2728-05-4
DILTIAZEM			
PHR2856-500MG	DILTIAZEM HYDROCHLORIDE	500MG	33286-22-5
PHR8701-100MG	DESACETYL DILTIAZEM HCL	100MG	75472-91-2

Cat. No.	Description	Pack Size	CAS No.
DIPHENHYDRAMINE			
PHR2268-1G	DIPHENHYDRAMINE CITRATE	1G	88637-37-0
PHR1015-1G	DIPHENHYDRAMINE HYDROCHLORIDE	1G	147-24-0
PHR1897-50MG	DIPHENHYDRAMINE IMPURITY B	50MG	4024-34-4
PHR2071-100MG	DIPHENHYDRAMINE IMPURITY D	100MG	91-01-0
PHR1993-50MG	DIPHENHYDRAMINE RELATED COMPOUND A	50MG	53499-40-4
PHR1899-50MG	DIPHENHYDRAMINE RELATED COMPOUND B	50MG	19574-66-4
PHR1898-50MG	DIPHENHYDRAMINE N-OXIDE HYDROCHLORIDE	50MG	13168-00-8
DOPAMINE			
PHR1090-1G	DOPAMINE HYDROCHLORIDE	1G	62-31-7
DORZOLAMIDE			
PHR1925-500MG	DORZOLAMIDE HYDROCHLORIDE	500MG	130693-82-2
PHR8704-30MG	DORZOLAMIDE RELATED COMPOUND B	30MG	
PHR8705-50MG	DORZOLAMIDE RELATED COMPOUND D	50MG	164455-27-0
DONEPEZIL			
PHR1584-1G	DONEPEZIL HYDROCHLORIDE	1G	120011-70-3
PHR2023-50MG	DONEPEZIL RELATED COMPOUND A	50MG	145546-80-1
PHR2025-50MG	DESBENZYL DONEPEZIL HYDROCHLORIDE	50MG	120013-39-0
DOXAZOSIN			
PHR3310-500MG	DOXAZOSIN MESYLATE	500MG	77883-43-3
PHR3313-50MG	DOXAZOSIN RELATED COMPOUND A	50MG	70918-74-0
PHR3311-50MG	DOXAZOSIN RELATED COMPOUND B	50MG	28888-44-0
PHR3308-50MG	DOXAZOSIN RELATED COMPOUND C	50MG	23680-84-4
PHR3312-50MG	DOXAZOSIN RELATED COMPOUND D	50MG	3663-80-7
PHR3309-50MG	DOXAZOSIN RELATED COMPOUND E	50MG	27631-29-4
PHR2954-40MG	DOXAZOSIN RELATED COMPOUND F	40MG	617677-53-9
DRONEDARONE			
PHR3615-50MG	DRONEDARONE RELATED COMPOUND C	50MG	1278585-68-4
DOXEPIN			
PHR8924-100MG	DOXEPIN RELATED COMPOUND B	100MG	4504-88-5
PHR8923-50MG	DOXEPIN RELATED COMPOUND C	50MG	2887-91-4
DULOXETINE			
PHR1865-500MG	DULOXETINE HYDROCHLORIDE	500MG	136434-34-9
PHR2226-25MG	DULOXETINE RELATED COMPOUND A	25MG	910138-96-4
PHR2167-50MG	DULOXETINE RELATED COMPOUND H	50MG	199191-66-7
ECONAZOLE			
PHR2791	ECONAZOLE		27220-47-9
PHR2232-500MG	ECONAZOLE IMPURITY D	500MG	24168-96-5
PHR2231-100MG	ECONAZOLE RELATED COMPOUND A	100MG	24155-42-8
PHR3637-50MG	ECONAZOLE RELATED COMPOUND C	50MG	
EDETATE DISODIUM			
PHR1068-1G	EDETATE DISODIUM	1G	6381-92-6
PHR1510-1G	EDETATE CALCIUM DISODIUM	1G	23411-34-9
EFAVIRENZ			
PHR2072-500MG	EFAVIRENZ	500MG	154598-52-4
PHR2068-50MG	EFAVIRENZ RELATED COMPOUND A	50MG	209414-27-7
PHR2069-50MG	EFAVIRENZ RELATED COMPOUND B	50MG	440124-96-9
PHR2070-50MG	EFAVIRENZ RELATED COMPOUND C	50MG	391860-73-4
ERYTHROMYCIN			
PHR1039-1G	ERYTHROMYCIN	1G	114-07-8
PHR1610-250MG	ERYTHROMYCIN ETHYLSUCCINATE	250MG	1264-62-6
PHR2780-200MG	AZAERYTHROMYCIN A	200MG	76801-85-9
ESCITALOPRAM			
PHR1733-1G	ESCITALOPRAM OXALATE	1G	219861-08-2
PHR3580-50MG	ESCITALOPRAM RELATED COMPOUND A	50MG	128173-53-5

Cat. No.	Description	Pack Size	CAS No.
ESTRADIOL			
PHR1353-1G	ESTRADIOL	1G	50-28-2
PHR2228-50MG	ESTRADIOL 3-METHYL ETHER	50MG	1035-77-4
PHR2858-500MG	ESTRADIOL BENZOATE	500MG	50-50-0
PHR2839-500MG	ESTRADIOL CYPIONATE	500MG	313-06-4
PHR2227-50MG	ESTRADIOL IMPURITY L	50MG	57-91-0
PHR2859-25MG	ESTRADIOL RELATED COMPOUND C	25MG	571-92-6
PHR2857-200MG	ESTRADIOL VALERATE	200MG	979-32-8
PHR1480-200MG	ETHINYL ESTRADIOL	200MG	57-63-6
ESTRIOL			
PHR2800-200MG	ESTRIOL	200MG	50-27-1
ETHAMBUTOL			
PHR1930-500MG	ETHAMBUTOL HYDROCHLORIDE	500MG	1070-11-7
PHR3710-50MG	ETHAMBUTOL RELATED COMPOUND B	50MG	134566-79-3
ETHOSUXIMIDE			
PHR1413-1G	ETHOSUXIMIDE	1MG	77-67-8
PHR2452-100MG	ETHOSUXIMIDE IMPURITY A	100MG	631-31-2
EXEMESTANE			
PHR1634-1G	EXEMESTANE	1G	107868-30-4
PHR1894-30MG	EXEMESTANE RELATED COMPOUND A	30MG	19457-55-7
PHR1895-30MG	EXEMESTANE RELATED COMPOUND B	30MG	184972-09-6
PHR8973-30MG	EXEMESTANE RELATED COMPOUND C	30MG	897-06-3
PHR1896-30MG	EXEMESTANE RELATED COMPOUND D	30MG	
FAMCICLOVIR			
PHR1696-1G	FAMCICLOVIR	1G	104227-87-4
PHR3743-50MG	FAMCICLOVIR RELATED COMPOUND A	50MG	246021-75-0
PHR3744-50MG	FAMCICLOVIR RELATED COMPOUND B	50MG	
PHR2860-50MG	FAMCICLOVIR RELATED COMPOUND E	50MG	452-06-2
PHR2861-50MG	FAMCICLOVIR RELATED COMPOUND F	50MG	10310-21-1
FAMOTIDINE			
PHR1055-1G	FAMOTIDINE	1G	76824-35-6
PHR2909-50MG	FAMOTIDINE RELATED COMPOUND B	50MG	89268-62-2
PHR2910-50MG	FAMOTIDINE RELATED COMPOUND C	50MG	76824-17-4
PHR2911-50MG	FAMOTIDINE RELATED COMPOUND D	50MG	76824-16-3
PHR2912-50MG	FAMOTIDINE RELATED COMPOUND E	50MG	129083-44-9
PHR2908-50MG	FAMOTIDINE RELATED COMPOUND F	50MG	107880-74-0
FELBAMATE			
PHR3639-50MG	FELBAMATE RELATED COMPOUND B	50MG	6326-19-8
FENBENDAZOLE			
PHR1832-200MG	FENBENDAZOLE	200MG	43210-67-9
PHR2862-50MG	FENBENDAZOLE RELATED COMPOUND A	50MG	10605-21-7
PHR2863-50MG	FENBENDAZOLE RELATED COMPOUND B	50MG	20367-38-8
FENOFIBRATE			
PHR1246-500MG	FENOFIBRATE	500MG	49562-28-9
PHR2864-50MG	FENOFIBRATE RELATED COMPOUND A	50MG	42019-78-3
PHR3315-50MG	FENOFIBRATE RELATED COMPOUND B	50MG	"42017-89-0"
FEXOFENADINE			
PHR1685-1G	FEXOFENADINE HYDROCHLORIDE	1G	153439-40-8
PHR9103-50MG	FEXOFENADINE METHYL ESTER	50MG	154825-96-4
PHR2454-50MG	FEXOFENADINE RELATED COMPOUND A	50MG	76811-98-8
FLAVOXATE			
PHR3316-500MG	FLAVOXATE HYDROCHLORIDE	500MG	3717-88-2
PHR3317-50MG	FLAVOXATE RELATED COMPOUND A	50MG	3468-01-7
PHR3008-50MG	FLAVOXATE RELATED COMPOUND B	50MG	90101-87-4
PHR3592-40MG	FLAVOXATE RELATED COMPOUND C	40MG	35888-94-9

Cat. No.	Description	Pack Size	CAS No.
FLUCONAZOLE			
PHR1160-1G	FLUCONAZOLE	1G	86386-73-4
PHR1900-50MG	FLUCONAZOLE IMPURITY A	50MG	89429-59-4
PHR2046-30MG	FLUCONAZOLE IMPURITY E	30MG	86404-63-9
PHR2047-30MG	FLUCONAZOLE IMPURITY F	30MG	118689-07-9
PHR2048-30MG	FLUCONAZOLE IMPURITY G	30MG	86386-76-7
PHR2037-30MG	FLUCONAZOLE IMPURITY H	30MG	
PHR1862-20MG	FLUCONAZOLE RELATED COMPOUND B	20MG	81886-51-3
PHR1826-20MG	FLUCONAZOLE RELATED COMPOUND C	20MG	514222-44-7
FLUMAZENIL			
PHR2702-500MG	FLUMAZENIL	500MG	78755-81-4
PHR3318-1G	FLUMAZENIL RELATED COMPOUND C	1MG	1188-33-6
FLUNIXIN			
PHR1442-1G	FLUNIXIN MEGLUMINE	1G	42461-84-7
PHR2305-50MG	FLUNIXIN RELATED COMPOUND A	50MG	
PHR2823-50MG	FLUNIXIN RELATED COMPOUND B	50MG	54396-44-0
PHR2413-50MG	FLUNIXIN RELATED COMPOUND C	50MG	1452-94-4
FLUOROURACIL			
PHR1227-500MG	FLUOROURACIL	500MG	51-21-8
PHR2866-50MG	FLUOROURACIL RELATED COMPOUND B	50MG	496-76-4
PHR2867-50MG	FLUOROURACIL RELATED COMPOUND E	50MG	1820-81-1
PHR2868-50MG	FLUOROURACIL RELATED COMPOUND F	50MG	56177-80-1
FLUOXETINE			
PHR1394-1G	FLUOXETINE HYDROCHLORIDE	1G	56296-78-7
PHR1952-30MG	FLUOXETINE RELATED COMPOUND A	30MG	79088-29-2
PHR2041-5ML	FLUOXETINE RELATED COMPOUND B	5ML	23580-89-4
PHR1953-30MG	FLUOXETINE RELATED COMPOUND C	30MG	1026723-45-4
FLURBIPROFEN			
PHR1499-500MG	FLURBIPROFEN	500MG	5104-49-4
PHR2710-100MG	FLURBIPROFEN RELATED COMPOUND A	100MG	6341-72-6
PHR2711-500MG	FLURBIPROFEN SODIUM	500MG	56767-76-1
FLUTICASONE			
PHR1702-200MG	FLUTICASONE PROPIONATE	200MG	80474-14-2
PHR2270-50MG	FLUTICASONE PROPIONATE IMPURITY A	50MG	65429-42-7
PHR8880-50MG	FLUTICASONE PROPIONATE IMPURITY E	50MG	
PHR1842-50MG	FLUTICASONE RELATED COMPOUND A	50MG	948566-12-9
PHR1843-10MG	FLUTICASONE RELATED COMPOUND C	10MG	80474-24-4
PHR1844-30MG	FLUTICASONE RELATED COMPOUND D	30MG	73205-13-7
PHR1845-50MG	FLUTICASONE RELATED COMPOUND F	50MG	1219174-94-3
PHR1846-30MG	FLUTICASONE RELATED COMPOUND J	30MG	28416-82-2
FOLIC ACID			
PHR1035-1G	FOLIC ACID	1G	59-30-3
PHR1334-500MG	FOLIC ACID IMPURITY A	500MG	4271-30-1
PHR2001-20MG	FOLIC ACID IMPURITY B	20MG	35011-47-3
PHR2002-100MG	FOLIC ACID IMPURITY D	100MG	119-24-4
PHR2003-20MG	FOLIC ACID IMPURITY F	20MG	1391194-56-1
PHR2334-50MG	10-FORMYLFOLIC ACID	50MG	
FOSCARNET			
PHR1436-1G	FOSCARNET SODIUM	1G	34156-56-4
PHR2869-50MG	FOSCARNET RELATED COMPOUND B	50MG	55920-24-6
PHR2870-0.5ML	FOSCARNET RELATED COMPOUND D	0.5ML	1474-78-8
FOSINOPRIL			
PHR8709-50MG	FOSINOPRIL RELATED COMPOUND A	50MG	95399-71-6
PHR3356-500MG	FOSINOPRIL SODIUM	500MG	88889-14-9

Cat. No.	Description	Pack Size	CAS No.
FUROSEMIDE			
PHR1057-1G	FUROSEMIDE	1G	54-31-9
PHR2457-100MG	FUROSEMIDE IMPURITY E	100MG	50-84-0
PHR1816-100MG	FUROSEMIDE RELATED COMPOUND B	100MG	3086-91-7
GABAPENTIN			
PHR1049-1G	GABAPENTIN	1G	60142-96-3
PHR1322-100MG	GABAPENTIN RELATED COMPOUND A	100MG	64744-50-9
PHR2011-100MG	GABAPENTIN RELATED COMPOUND B	100MG	133481-09-1
PHR2012-20MG	GABAPENTIN RELATED COMPOUND D	20MG	1076198-17-8
PHR2013-50MG	GABAPENTIN RELATED COMPOUND E	50MG	67950-95-2
GATIFLOXACIN			
PHR3849-500MG	GATIFLOXACIN	500MG	
PHR3711-50MG	GATIFLOXACIN RELATED COMPOUND A	50MG	616205-76-6
PHR3712-50MG	GATIFLOXACIN RELATED COMPOUND D	50MG	112811-57-1
GLIMEPIRIDE			
PHR1617-1G	GLIMEPIRIDE	1G	93479-97-1
PHR1726-40MG	GLIMEPIRIDE SULFONAMIDE	40MG	
PHR1725-40MG	GLIMEPIRIDE RELATED COMPOUND A	40MG	684286-46-2
PHR1727-40MG	GLIMEPIRIDE RELATED COMPOUND C	40MG	119018-30-3
GLIPIZIDE			
PHR1818-1G	GLIPIZIDE	1G	29094-61-9
PHR1813-50MG	GLIPIZIDE RELATED COMPOUND A	50MG	33288-71-0
GLUCOSAMINE			
PHR1199-500MG	GLUCOSAMINE HYDROCHLORIDE	500MG	66-84-2
PHR2290-50MG	GLUCOSAMINE IMPURITY B	50MG	13185-73-4
PHR2460-100MG	GLUCOSAMINE IMPURITY E	100MG	67-47-0
PHR1432-1G	N-ACETYLGUCOSAMINE	1G	7512-17-6
GLUTATHIONE			
PHR1359-500MG	GLUTATHIONE	500MG	70-18-8
PHR2462-100MG	GLUTATHIONE IMPURITY A	100MG	19246-18-5
PHR2461-100MG	GLUTATHIONE IMPURITY C	100MG	27025-41-8
PHR2463-100MG	GLUTATHIONE IMPURITY D	100MG	636-58-8
GLYBURIDE			
PHR1287-1G	GLYBURIDE (GLIBENCLAMIDE)	1G	10238-21-8
PHR2414-50MG	GLYBURIDE RELATED COMPOUND A	50MG	
PHR3594-50MG	GLYBURIDE RELATED COMPOUND B	50MG	21165-77-5
GLYCOPYRROLATE			
PHR2712-50MG	GLYCOPYRROLATE ERYTHRO ISOMER	50MG	
PHR2465-100MG	GLYCOPYRROLATE RELATED COMPOUND A	100MG	
PHR2466-100MG	GLYCOPYRROLATE RELATED COMPOUND B	100MG	13118-11-1
PHR2467-50MG	GLYCOPYRROLATE RELATED COMPOUND C	50MG	427-49-6
PHR2713-30MG	GLYCOPYRROLATE RELATED COMPOUND I	30MG	
PHR2468-100MG	GLYCOPYRROLATE RELATED COMPOUND L	100MG	19833-96-6
GRANISETRON			
PHR1616-500MG	GRANISETRON HYDROCHLORIDE	500MG	107007-99-8
PHR2871-50MG	GRANISETRON RELATED COMPOUND D	50MG	50890-83-0
GUAIACOL			
PHR1136-1.5G	GUAIACOL	1.5G	90-05-1
PHR2474-100MG	GUAIACOL IMPURITY A	100MG	120-80-9
PHR2469-100MG	GUAIACOL IMPURITY C	100MG	91-16-7
PHR9076	GUAIACOL IMPURITY D		18113-18-3
PHR2475-100MG	GUAIACOL IMPURITY E	100MG	93-58-3
PHR2471-100MG	GUAIACOL IMPURITY F	100MG	
PHR2473-100MG	GUAIACOL IMPURITY G	100MG	150-76-5
PHR2472-100MG	GUAIACOL IMPURITY H	100MG	150-19-6

Cat. No.	Description	Pack Size	CAS No.
GUAIFENESIN			
PHR1027-1G	GUAIFENESIN	1G	93-14-1
PHR2325-50MG	GUAIFENESIN IMPURITY C	50MG	
PHR2326-50MG	GUAIFENESIN IMPURITY D	50MG	16929-60-5
PHR1543-50MG	GUAIFENESIN IMPURITY B	50MG	14007-09-1
HALOPERIDOL			
PHR1724-400MG	HALOPERIDOL	400MG	52-86-8
PHR8951-25MG	HALOPERIDOL RELATED COMPOUND A	25MG	67987-08-0
HYDROCORTISONE			
PHR1014-500MG	HYDROCORTISONE	500MG	50-23-7
PHR1606-1G	HYDROCORTISONE ACETATE	1G	50-03-3
PHR1926-500MG	HYDROCORTISONE HEMISUCCINATE	500MG	83784-20-7
PHR2156-100MG	HYDROCORTISONE IMPURITY D	100MG	53-35-0
PHR2158-100MG	HYDROCORTISONE IMPURITY F	100MG	152-58-9
PHR2159-50MG	HYDROCORTISONE IMPURITY G	50MG	14760-49-7
PHR2160-50MG	HYDROCORTISONE IMPURITY J	50MG	
PHR2161-50MG	HYDROCORTISONE IMPURITY K	50MG	640-87-9
PHR2162-50MG	HYDROCORTISONE IMPURITY L	50MG	641-77-0
PHR2163-50MG	HYDROCORTISONE IMPURITY M	50MG	566-35-8
PHR3185-500MG	HYDROCORTISONE VALERATE	500MG	57524-89-7
PHR2164-50MG	6A-HYDROXYHYDROCORTISONE	50MG	2242-98-0
HYDROXYZINE			
PHR1607-1G	HYDROXYZINE DIHYDROCHLORIDE	1G	2192-20-3
PHR2714-500MG	HYDROXYZINE PAMOATE	500MG	10246-75-0
PHR2872-50MG	HYDROXYZINE RELATED COMPOUND A	50MG	303-26-4
HYOSCYAMINE			
PHR3012-25MG	HYOSCYAMINE RELATED COMPOUND A	25MG	
IBANDRONATE			
PHR2477-50MG	IBANDRONATE RELATED COMPOUND A	50MG	625120-81-2
PHR2715-500MG	IBANDRONATE SODIUM	500MG	138926-19-9
IBUPROFEN			
PHR1004-1G	IBUPROFEN	1G	15687-27-1
PHR3625-50MG	IBUPROFEN EP IMPURITY M	50MG	60057-62-7
PHR1146-500MG	IBUPROFEN RELATED COMPOUND C	500MG	38861-78-8
PHR1933-20MG	IBUPROFEN IMPURITY A	20MG	66622-47-7
PHR1934-20MG	IBUPROFEN IMPURITY B	20MG	3585-49-7
PHR3326-50MG	IBUPROFEN IMPURITY C	50MG	59512-17-3
PHR1964-30MG	IBUPROFEN IMPURITY F	30MG	65322-85-2
PHR1978-50MG	IBUPROFEN IMPURITY J	50MG	65813-55-0
PHR1965-50MG	IBUPROFEN IMPURITY G	50MG	1391054-15-1
PHR2028-50MG	IBUPROFEN IMPURITY H	50MG	
PHR1984-50MG	IBUPROFEN IMPURITY I	50MG	
PHR3324-50MG	IBUPROFEN IMPURITY K	50MG	43153-07-7
PHR3325-50MG	IBUPROFEN IMPURITY L	50MG	53949-53-4
PHR1935-20MG	IBUPROFEN IMPURITY N	20MG	3585-52-2
PHR1983-30MG	IBUPROFEN IMPURITY O	30MG	64451-76-9
PHR1982-50MG	IBUPROFEN IMPURITY R	50MG	
PHR2145-50MG	IBUPROFEN 1,2-PROPANETRIOL ESTERS	50MG	
PHR2144-50MG	IBUPROFEN 1,2-PROPYLENE GLYCOL ESTERS	50MG	
PHR2142-100MG	IBUPROFEN 1,3-BUTYLENE GLYCOL ESTERS	100MG	
PHR2143-50MG	IBUPROFEN 2,3-BUTYLENE GLYCOL ESTER	50MG	
PHR2138-50MG	IBUPROFEN ETHYL ESTER	50MG	41283-72-1
PHR2133-50MG	IBUPROFEN INTERMEDIATE	50MG	121734-80-3
PHR2137-1ML	IBUPROFEN METHYL ESTER	1ML	61566-34-5

Cat. No.	Description	Pack Size	CAS No.
IMIQUIMOD			
PHR2716-500MG	IMIQUIMOD	500MG	99011-02-6
PHR2489-50MG	IMIQUIMOD RELATED COMPOUND A	50MG	99010-24-9
PHR3014-50MG	IMIQUIMOD RELATED COMPOUND B	50MG	99010-63-6
PHR2490-50MG	IMIQUIMOD RELATED COMPOUND C	50MG	99010-64-7
PHR3013-50MG	IMIQUIMOD RELATED COMPOUND D	50MG	853792-81-1
INDAPAMIDE			
PHR1660-1G	INDAPAMIDE	1G	26807-65-8
PHR3665-100MG	INDAPAMIDE RELATED COMPOUND A	100MG	63968-75-2
INDOMETHACIN			
PHR1247-500MG	INDOMETHACIN	500MG	53-86-1
PHR2770-50MG	INDOMETHACIN RELATED COMPOUND A	50MG	2882-15-7
PHR2581-50MG	INDOMETHACIN RELATED COMPOUND B	50MG	74-11-3
IOPAMIDOL			
PHR1985-500MG	IOPAMIDOL	500MG	60166-93-0
PHR3753-100MG	IOPAMIDOL RELATED COMPOUND A	100MG	60166-98-5
IPRATROPIUM			
PHR2796-400MG	IPRATROPIUM BROMIDE	400MG	66985-17-9
PHR3017-25MG	IPRATROPIUM BROMIDE RELATED COMPOUND A	25MG	58005-18-8
PHR2873-100MG	IPRATROPIUM BROMIDE RELATED COMPOUND C	100MG	552-63-6
IRBESARTAN			
PHR1443-1G	IRBESARTAN	1G	138402-11-6
PHR1958-50MG	IRBESARTAN RELATED COMPOUND A	50MG	748812-53-5
IRINOTECAN			
PHR2717-200MG	IRINOTECAN HYDROCHLORIDE	200MG	136572-09-3
PHR2491-25MG	IRINOTECAN RELATED COMPOUND A	25MG	19685-09-7
PHR2718-20MG	IRINOTECAN RELATED COMPOUND B	20MG	86639-52-3
PHR2492-30MG	IRINOTECAN RELATED COMPOUND E	30MG	78287-27-1
ISOFLURANE			
PHR2874-1.2ML	ISOFLURANE	1.2ML	26675-46-7
PHR2875-0.5ML	ISOFLURANE RELATED COMPOUND A	0.5ML	32778-08-8
PHR2876-0.5ML	ISOFLURANE RELATED COMPOUND B	0.5ML	1885-48-9
KETOPROFEN			
PHR2626-500MG	DEKKETOPROFEN TROMETHAMINE	500MG	156604-79-4
PHR1375-500MG	KETOPROFEN	500MG	22071-15-4
PHR3018-50MG	KETOPROFEN RELATED COMPOUND C	50MG	68432-95-1
PHR2955-40MG	KETOPROFEN RELATED COMPOUND D	40MG	66067-44-5
PHR2877-50MG	KETOPROFEN RELATED COMPOUND F	50MG	42872-30-0
KETOROLAC			
PHR8721-100MG	KETOROLAC IMPURITY J	100MG	108061-03-6
PHR1140-500MG	KETOROLAC TROMETHAMINE	500MG	74103-07-4
PHR2493-50MG	KETOROLAC RELATED COMPOUND A	50MG	167105-80-8
PHR2494-50MG	KETOROLAC RELATED COMPOUND C	50MG	113502-52-6
PHR3605-50MG	KETOROLAC RELATED COMPOUND D	50MG	113502-55-9
LAMIVUDINE			
PHR1365-1G	LAMIVUDINE	1G	134678-17-4
PHR8731-50MG	LAMIVUDINE IMPURITY D	50MG	
PHR3357-50MG	LAMIVUDINE IMPURITY J	50MG	145986-07-8
PHR8938-50MG	LAMIVUDINE (R)-SULFOXIDE	50MG	
LAMOTRIGINE			
PHR1392-1G	LAMOTRIGINE	1G	84057-84-1
PHR8729-20MG	LAMOTRIGINE RELATED COMPOUND C	20MG	252186-78-0
PHR8730-20MG	LAMOTRIGINE RELATED COMPOUND D	20MG	252186-79-1
LANSOPRAZOLE			
PHR1390-1G	LANSOPRAZOLE	1G	103577-45-3
PHR3019-50MG	LANSOPRAZOLE RELATED COMPOUND A	50MG	131926-99-3

Cat. No.	Description	Pack Size	CAS No.
LANSOPRAZOLE			
PHR2956-100MG	LANSOPRAZOLE	1G	103577-45-3
LATANOPROST			
PHR1884-50MG	LATANOPROST	50MG	130209-82-4
PHR3022-50MG	LATANOPROST RELATED COMPOUND B	50MG	145773-22-4
LEFLUNOMIDE			
PHR1378-1G	LEFLUNOMIDE	1G	75706-12-6
PHR1568-200MG	LEFLUNOMIDE RELATED COMPOUND A	200MG	455-14-1
PHR2878-50MG	LEFLUNOMIDE RELATED COMPOUND B	50MG	108605-62-5
LEVETIRACETAM			
PHR1447-1G	LEVETIRACETAM	1G	102767-28-2
PHR2417-30MG	LEVETIRACETAM RACEMIC MIXTURE	30MG	
PHR2428-50MG	LEVETIRACETAM RELATED COMPOUND A	50MG	
PHR2432-50MG	LEVETIRACETAM RELATED COMPOUND B	50MG	7682-20-4
LEVALBUTEROL			
PHR3457-500MG	LEVALBUTEROL HYDROCHLORIDE	500MG	50293-90-8
PHR3622-50MG	LEVALBUTEROL RELATED COMPOUND C HYDROCHLORIDE	50MG	27566-30-9
PHR3623-50MG	LEVALBUTEROL RELATED COMPOUND E	50MG	754926-25-5
PHR3624-50MG	LEVALBUTEROL RELATED COMPOUND F	50MG	174607-68-2
LEVOCARNITINE			
PHR1569-1G	LEVOCARNITINE	1G	541-15-1
PHR2879-100MG	LEVOCARNITINE RELATED COMPOUND A	100MG	"6538-82-5"
LEVODOPA			
PHR1271-500MG	LEVODOPA	500MG	59-92-7
PHR2095-100MG	LEVODOPA RELATED COMPOUND A	100MG	27244-64-0
PHR2096-100MG	LEVODOPA RELATED COMPOUND B	100MG	7636-26-2
LEVOFLOXACIN			
PHR1697-1G	LEVOFLOXACIN HEMIHYDRATE	1G	138199-71-0
PHR2229-100MG	LEVOFLOXACIN RELATED COMPOUND B	100MG	100986-89-8
PHR2230-50MG	LEVOFLOXACIN RELATED COMPOUND C	50MG	177472-30-9
PHR3023-50MG	LEVOFLOXACIN RELATED COMPOUND E	50MG	1346603-62-0
LIDOCAINE			
PHR1034-1G	LIDOCAINE	1G	137-58-6
PHR1257-500MG	LIDOCAINE HYDROCHLORIDE	500MG	6108-05-0
PHR1669-100MG	LIDOCAINE RELATED COMPOUND A	100MG	201-758-7
PHR1670-100MG	LIDOCAINE RELATED COMPOUND B	100MG	
PHR1672-100MG	LIDOCAINE RELATED COMPOUND E	100MG	
PHR1673-100MG	LIDOCAINE RELATED COMPOUND H	100MG	1131-01-7
PHR1674-100MG	LIDOCAINE RELATED COMPOUND I	100MG	
PHR1675-100MG	LIDOCAINE RELATED COMPOUND J	100MG	
PHR1977-50MG	LIDOCAINE IMPURITY D	50MG	
PHR1994-50MG	LIDOCAINE IMPURITY F	50MG	857170-72-0
PHR1995-50MG	LIDOCAINE IMPURITY G	50MG	
LINEZOLID			
PHR1885-500MG	LINEZOLID	500MG	165800-03-3
PHR3027-50MG	LINEZOLID RELATED COMPOUND A	50MG	168828-84-0
PHR3028-50MG	LINEZOLID RELATED COMPOUND B	50MG	216868-57-4
PHR3024-50MG	LINEZOLID RELATED COMPOUND C	50MG	168828-90-8
PHR3025-50MG	LINEZOLID RELATED COMPOUND D	50MG	174649-09-3
PHR3026-50MG	LINEZOLID R-ISOMER	50MG	872992-20-6
LISINAPRIL			
PHR1143-1G	LISINAPRIL	1G	83915-83-7
PHR2239-50MG	LISINAPRIL IMPURITY A	50MG	1012-05-1
PHR2495-100MG	LISINAPRIL IMPURITY B	100MG	6192-52-5
PHR8894-50MG	LISINAPRIL R,S,S-DIKETOPIPERAZINE	50MG	219677-82-4
PHR8931-50MG	LISINAPRIL S,S,S-DIKETOPIPERAZINE	50MG	328385-86-0

Cat. No.	Description	Pack Size	CAS No.
LOMUSTINE			
PHR3029-500MG	LOMUSTINE	500MG	13010-47-4
PHR3030-50MG	LOMUSTINE RELATED COMPOUND B	50MG	13908-11-7
PHR2924-500MG	LOMUSTINE RELATED COMPOUND C	500MG	2387-23-7
LOPERAMIDE			
PHR1162-1G	LOPERAMIDE HYDROCHLORIDE	1G	34552-83-5
PHR2916-50MG	LOPERAMIDE RELATED COMPOUND F	50MG	106900-12-3
LORATADINE			
PHR1376-1G	LORATADINE	1G	79794-75-5
PHR2310-50MG	LORATADINE IMPURITY A	50MG	133284-74-9
PHR2311-50MG	LORATADINE IMPURITY E	50MG	170727-59-0
PHR2313-100MG	LORATADINE IMPURITY H	100MG	29976-53-2
PHR2308-50MG	LORATADINE RELATED COMPOUND B	50MG	38092-89-6
PHR2309-50MG	LORATADINE RELATED COMPOUND C	50MG	31251-41-9
LOSARTAN			
PHR1602-1G	LOSARTAN POTASSIUM	1MG	124750-99-8
PHR8605-50MG	LOSARTAN RELATED COMPOUND C	50MG	114798-36-6
LOVASTATIN			
PHR1285-1G	LOVASTATIN	1G	75330-75-5
PHR2240-50MG	LOVASTATIN IMPURITY A	50MG	73573-88-3
PHR3741-50MG	LOVASTATIN EP IMPURITY B	50MG	75225-51-3
PHR2278-50MG	LOVASTATIN IMPURITY C	50MG	109273-98-5
PHR2279-20MG	LOVASTATIN IMPURITY D	20MG	149756-62-7
PHR2277-100MG	LOVASTATIN RELATED COMPOUND A	100MG	77517-29-4
LUFENURON			
PHR3031-500MG	LUFENURON	500MG	103055-07-8
MAFENIDE			
PHR3032-100MG	MAFENIDE RELATED COMPOUND A	100MG	3240-35-5
MECLIZINE			
PHR3571-1G	MECLIZINE HYDROCHLORIDE	1G	31884-77-2
PHR3034-30MG	MECLIZINE RELATED COMPOUND A	30MG	119-56-2
PHR3728-50MG	MECLIZINE RELATED COMPOUND B	50MG	
MEDROXYPROGESTERONE			
PHR1589-500MG	MEDROXYPROGESTERONE ACETATE	500MG	71-58-9
PHR3717-50MG	MEDROXYPROGESTERONE ACETATE RELATED COMPOUND A	50MG	
MELATONIN			
PHR1767-500MG	MELATONIN	500MG	73-31-4
PHR2241-50MG	MELATONIN RELATED COMPOUND A	50MG	608-07-1
MELOXICAM			
PHR1799-500MG	MELOXICAM	500MG	71125-38-7
PHR2327-50MG	MELOXICAM IMPURITY C	50MG	
PHR2131-50MG	MELOXICAM IMPURITY D	50MG	
PHR2503-50MG	MELOXICAM IMPURITY E	50MG	35511-15-0
PHR2500-50MG	MELOXICAM RELATED COMPOUND A	50MG	24683-26-9
PHR2501-50MG	MELOXICAM RELATED COMPOUND B	50MG	7305-71-7
PHR2502-50MG	MELOXICAM RELATED COMPOUND C	50MG	118854-48-1
PHR3395-50MG	MELOXICAM RELATED COMPOUND D	50MG	
MEMANTINE			
PHR1886-1G	MEMANTINE HYDROCHLORIDE	1G	41100-52-1
PHR2449-100MG	MEMANTINE RELATED COMPOUND A	100MG	702-79-4
PHR2450-50MG	MEMANTINE RELATED COMPOUND B	50MG	707-37-9
PHR2451-100MG	MEMANTINE RELATED COMPOUND C	100MG	707-36-8
PHR2583-100MG	MEMANTINE RELATED COMPOUND D	100MG	941-37-7
PHR2931-30MG	MEMANTINE RELATED COMPOUND E	30MG	351329-88-9
PHR2759-100MG	MEMANTINE RELATED COMPOUND F	100MG	356572-08-2
PHR2761-50MG	MEMANTINE RELATED COMPOUND G	50MG	33103-93-4

Cat. No.	Description	Pack Size	CAS No.
MEMANTINE			
PHR2764-50MG	MEMANTINE RELATED COMPOUND H	50MG	15210-60-3
PHR2765-50MG	MEMANTINE RELATED COMPOUND I	50MG	80121-67-1
MEROPENEM			
PHR1772-500MG	MEROPENEM	500MG	119478-56-7
MESALAMINE			
PHR1060-1G	MESALAMINE (MESALAZINE)	1G	89-57-6
PHR1739-200MG	2-AMINOPHENOL (MESALAZINE IMPURITY C)	200MG	95-55-6
MESNA			
PHR1570-500MG	MESNA	500MG	19767-45-4
PHR2818-100MG	MESNA RELATED COMPOUND A	100MG	
PHR2121-100MG	MESNA RELATED COMPOUND B	100MG	
METAMIZOLE			
PHR1800-100MG	METAMIZOLE SODIUM	100MG	5907-38-0
PHR2783-100MG	METAMIZOLE IMPURITY A	100MG	1672-58-8
PHR2784-100MG	METAMIZOLE IMPURITY B	100MG	83-07-8
PHR2785-100MG	METAMIZOLE IMPURITY C	100MG	856307-27-2
PHR2786-100MG	METAMIZOLE IMPURITY D	100MG	58-15-1
METAXALONE			
PHR3036-500MG	METAXALONE	500MG	1665-48-1
METFORMINE			
PHR1084-500MG	METFORMIN HYDROCHLORIDE	500MG	1115-70-4
PHR1331-500MG	METFORMIN RELATED COMPOUND A	500MG	461-58-5
PHR1505-50MG	METFORMIN RELATED COMPOUND B	50MG	1674-62-0
PHR1506-50MG	METFORMIN RELATED COMPOUND C	50MG	1985-46-2
PHR1532-200MG	METFORMIN IMPURITY F	200MG	506-59-2
PHR1274-1G	METFORMIN IMPURITY D	1G	108-78-1
METHIMAZOLE			
PHR3037-500MG	METHIMAZOLE	500MG	60-56-0
PHR3038-100MG	METHIMAZOLE RELATED COMPOUND A	100MG	122-07-6
PHR3039-100MG	METHIMAZOLE RELATED COMPOUND C	100MG	14486-52-3
METHOTREXATE			
PHR1396-1G	METHOTREXATE	1G	59-05-2
METHOXSALEN			
PHR3040-500MG	METHOXSALEN	500MG	298-81-7
PHR3041-50MG	METHOXSALEN RELATED COMPOUND A	50MG	484-20-8
METHYLDOPA			
PHR2672-1G	METHYLDOPA	1G	41372-08-1
PHR2673-30MG	METHYLDOPA RELATED COMPOUND B	30MG	35026-10-9
PHR2674-30MG	METHYLDOPA RELATED COMPOUND C	30MG	
METHYLPHENIDATE			
PHR3576-100MG	METHYLPHENIDATE RELATED COMPOUND A	100MG	19395-40-5
PHR8987-1ML	METHYLPHENIDATE ERYTHRO ISOMER	1ML	29419-97-4
METHYLPREDNISOLONE			
PHR1717-500MG	METHYLPREDNISOLONE	500MG	83-43-2
PHR1718-500MG	METHYLPREDNISOLONE ACETATE	500MG	53-36-1
PHR3045-100MG	METHYLPREDNISOLONE HEMISUCCINATE RELATED COMPOUND D	100MG	119657-85-1
PHR3043-50MG	METHYLPREDNISOLONE RELATED COMPOUND C	50MG	61919-52-6
METHYL SALICYLATE			
PHR1214-3.2ML	METHYL SALICYLATE	3.2ML	119-36-8
PHR3042-100MG	METHYL SALICYLATE RELATED COMPOUND A	100MG	5985-24-0
METOCLOPRAMIDE			
PHR1132-1G	METOCLOPRAMIDE HYDROCHLORIDE	1G	
PHR2507-100MG	METOCLOPRAMIDE IMPURITY C	100MG	7206-70-4
PHR2286-100MG	METOCLOPRAMIDE IMPURITY E	100MG	100-36-7
PHR2506-100MG	METOCLOPRAMIDE IMPURITY H	100MG	

Cat. No.	Description	Pack Size	CAS No.
METOCLOPRAMIDE			
PHR2283-100MG	METOCLOPRAMIDE RELATED COMPOUND A	100MG	5608-13-9
PHR2284-100MG	METOCLOPRAMIDE RELATED COMPOUND B	100MG	4093-31-6
PHR2285-100MG	METOCLOPRAMIDE RELATED COMPOUND D	100MG	4093-29-2
METOPROLOL			
PHR3046-500MG	METOPROLOL FUMARATE	500MG	119637-66-0
PHR9068-100MG	METOPROLOL EP IMPURITY H	100MG	62572-94-5
PHR9088-100MG	METOPROLOL EP IMPURITY J	100MG	163685-37-8
PHR9192-50MG	METOPROLOL EP IMPURITY N	50MG	6452-57-9
PHR9101-100MG	METOPROLOL EP IMPURITY E	100MG	163685-38-9
PHR9102-100MG	METOPROLOL EP IMPURITY D	100MG	62572-90-1
PHR9188-50MG	METOPROLOL EP IMPURITY G	50MG	501-94-0
PHR2509-50MG	METOPROLOL RELATED COMPOUND A	50MG	109632-08-8
PHR2510-30MG	METOPROLOL RELATED COMPOUND C	30MG	
PHR3641-50MG	METOPROLOL RELATED COMPOUND D	50MG	1486464-40-7
PHR2511-500MG	METOPROLOL SUCCINATE	500MG	98418-47-4
PHR1076-1G	METOPROLOL TARTRATE	1G	56392-17-7
METRONIDAZOLE			
PHR1052-1G	METRONIDAZOLE	1G	443-48-1
PHR2514-100MG	METRONIDAZOLE IMPURITY B	100MG	3034-38-6
PHR2512-100MG	METRONIDAZOLE IMPURITY E	100MG	705-19-1
PHR2165-100MG	METRONIDAZOLE IMPURITY F	100MG	16156-94-8
PHR3331-50MG	METRONIDAZOLE IMPURITY G	50MG	1010-93-1
MICONAZOLE			
PHR1618-1G	MICONAZOLE	1G	22916-47-8
PHR1163-1G	(±)-MICONAZOLE NITRATE SALT	1G	22832-87-7
PHR3561-100MG	MICONAZOLE EP IMPURITY F	100MG	27220-35-5
PHR3562-100MG	MICONAZOLE RELATED COMPOUND I	100MG	24169-00-4
MIDODRINE			
PHR3047-500MG	MIDODRINE HYDROCHLORIDE	500MG	3092-17-9
PHR3048-100MG	MIDODRINE RELATED COMPOUND A	100MG	"3600-87-1"
MILRINONE			
PHR2515-1G	MILRINONE	1G	78415-72-2
PHR3049-100MG	MILRINONE RELATED COMPOUND A	100MG	80047-24-1
MINOXIDIL			
PHR2123-200MG	MINOXIDIL	200MG	38304-91-5
PHR3050-50MG	MINOXIDIL RELATED COMPOUND E	50MG	24867-26-3
PHR3051-500MG	MINOXIDIL RELATED COMPOUND E	500MG	70476-82-3
MIRTAZAPINE			
PHR8751-50MG	MIRTAZAPINE RELATED COMPOUND A	50MG	61337-68-6
PHR8736-50MG	MIRTAZAPINE RELATED COMPOUND B	50MG	155172-12-6
PHR8737-100MG	MIRTAZAPINE RELATED COMPOUND D	100MG	191546-97-1
MOEXIPRIL			
PHR3606-50MG	MOEXIPRIL RELATED COMPOUND A	50MG	103775-14-0
PHR3719-50MG	MOEXIPRIL RELATED COMPOUND B	50MG	103733-51-3
PHR3052-50MG	MOEXIPRIL RELATED COMPOUND E	50MG	82586-62-7
PHR3053-50MG	MOEXIPRIL RELATED COMPOUND F	50MG	82717-96-2
MOXIDECTIN			
PHR1827-500MG	MOXIDECTIN	500MG	113507-06-5
MOXIFLOXACIN			
PHR1542-1G	MOXIFLOXACIN HYDROCHLORIDE	1G	186826-86-8
MYCOPHENOLATE MOFETIL			
PHR1451-1G	MYCOPHENOLATE MOFETIL	1G	128794-94-5
PHR2776-100MG	MYCOPHENOLATE MOFETIL RELATED COMPOUND B	100MG	26675-76-3
NALTREXONE			
PHR8930-25MG	NALTREXONE RELATED COMPOUND A	25G	131670-05-8

Cat. No.	Description	Pack Size	CAS No.
NAPHAZOLINE			
PHR2455-50MG	NAPHAZOLINE RELATED COMPOUND A	50MG	36321-43-4
PHR1836-500MG	NAPHAZOLINE HYDROCHLORIDE	500MG	550-99-2
NAPROXEN			
PHR1040-500MG	NAPROXEN	500MG	22204-53-1
PHR2254-50MG	NAPROXEN 1,2,3 PROPANETRIOL ESTERS	50MG	
PHR3332-50MG	NAPROXEN EP IMPURITY D	50MG	116883-62-6
PHR2521-100MG	NAPROXEN IMPURITY A	100MG	52079-10-4
PHR2520-100MG	NAPROXEN IMPURITY B	100MG	89617-86-7
PHR2519-100MG	NAPROXEN IMPURITY C	100MG	84236-26-0
PHR2529-100MG	NAPROXEN IMPURITY E	100MG	26159-35-3
PHR2528-100MG	NAPROXEN IMPURITY F	100MG	31220-35-6
PHR2526-100MG	NAPROXEN IMPURITY H	100MG	5111-66-0
PHR2522-100MG	NAPROXEN IMPURITY I	100MG	23981-47-7
PHR3333-50MG	NAPROXEN IMPURITY J	50MG	21388-17-0
PHR2525-100MG	NAPROXEN IMPURITY M	100MG	93-04-9
PHR2524-100MG	NAPROXEN IMPURITY N	100MG	5111-65-9
PHR2527-100MG	NAPROXEN RELATED COMPOUND A	100MG	
PHR2518-100MG	NAPROXEN RELATED COMPOUND G	100MG	23979-41-1
PHR2517-100MG	NAPROXEN RELATED COMPOUND K	100MG	77301-42-9
PHR2523-100MG	NAPROXEN RELATED COMPOUND L	100MG	3900-45-6
PHR1165-1G	NAPROXEN SODIUM	1G	26159-34-2
PHR2530-100MG	RAC-NAPROXEN	100MG	23981-80-8
NATEGLINIDE			
PHR3058-30MG	NATEGLINIDE RELATED COMPOUND A	30MG	7077-05-6
NEBIVOLOL			
PHR3751-150MG	NEBIVOLOL RELATED COMPOUND C	150MG	876514-31-7
PHR3730-100MG	NEBIVOLOL RELATED COMPOUND D	100MG	793669-26-8
NEVIRAPINE			
PHR1757-1G	NEVIRAPINE	1G	129618-40-2
PHR1760-50MG	NEVIRAPINE RELATED COMPOUND A	50MG	133627-17-5
PHR1761-50MG	NEVIRAPINE RELATED COMPOUND B	50MG	287980-84-1
PHR1762-50MG	NEVIRAPINE RELATED COMPOUND C	50MG	287980-85-2
NICOTINAMIDE			
PHR2087-100MG	NICOTINAMIDE IMPURITY B	100MG	100-54-9
PHR2088-100MG	NICOTINAMIDE IMPURITY C	100MG	1452-77-3
PHR2089-100MG	NICOTINAMIDE IMPURITY D	100MG	1453-82-3
PHR2090-100MG	NICOTINAMIDE IMPURITY E	100MG	1986-81-8
PHR2091-100MG	NICOTINAMIDE IMPURITY F	100MG	1005-24-9
NICOTINE			
PHR2532-500MG	NICOTINE BITARTRATE DIHYDRATE	500MG	6019-06-3
PHR3060-50MG	NICOTINE RELATED COMPOUND A/EP IMPURITY	50MG	2743-90-0
PHR3061-50MG	NICOTINE RELATED COMPOUND B	50MG	487-19-4
PHR2533-50MG	NICOTINE RELATED COMPOUND C	50MG	486-56-6
PHR3062-50MG	NICOTINE RELATED COMPOUND D	50MG	
PHR3063-50MG	NICOTINE RELATED COMPOUND F	50MG	5746-86-1
PHR3064-50MG	NICOTINE RELATED COMPOUND G	50MG	494-52-0
NILOTINIB			
PHR2833-500MG	NILOTINIB	500MG	923288-90-8
PHR2834-50MG	NILOTINIB RELATED COMPOUND A	50MG	641571-11-1
PHR2835-50MG	NILOTINIB RELATED COMPOUND C	50MG	2458-12-0
NIMODIPINE			
PHR1293-1G	NIMODIPINE	1MG	66085-59-4
PHR3560-100MG	NIMODIPINE RELATED COMPOUND B	100MG	70172-96-2
NITAZOXANIDE			
PHR2888-400MG	NITAZOXANIDE	400MG	55981-09-4

Cat. No.	Description	Pack Size	CAS No.
NITAZOXANIDE			
PHR2889-50MG	NITAZOXANIDE RELATED COMPOUND A	50MG	121-66-4
NITROFURAZONE			
PHR1196-1G	NITROFURAZONE	1G	59-87-0
PHR3067-1G	NITROFURAZONE RELATED COMPOUND A	1G	112537-96-9
NORETHINDRONE			
PHR1714-1G	NORETHINDRONE (NORETHISTERONE)	1G	68-22-4
PHR2534-200MG	NORETHINDRONE ACETATE	200MG	51-98-9
NORFLOXACIN			
PHR1614-1G	NORFLOXACIN	1G	70458-96-7
PHR2233-100MG	NORFLOXACIN RELATED COMPOUND A	100MG	68077-26-9
NORGESTREL			
PHR1850-500MG	LEVONORGESTREL	500MG	797-63-7
PHR2487-200MG	NORGESTREL	200MG	6533-00-2
PHR3068-50MG	NORGESTREL RELATED COMPOUND L	50MG	21800-83-9
OFLOXACIN			
PHR1168-1G	OFLOXACIN	1G	82419-36-1
PHR2234-100MG	OFLOXACIN IMPURITY A	100MG	82419-35-0
PHR9016-50MG	OFLOXACIN IMPURITY C (PHEUR)	50MG	95848-94-5
PHR2235-50MG	OFLOXACIN RELATED COMPOUND A	50MG	82419-52-1
OLANZAPINE			
PHR1825-1G	OLANZAPINE	1G	132539-06-1
PHR2536-30MG	OLANZAPINE RELATED COMPOUND A	30MG	138564-59-7
PHR2537-30MG	OLANZAPINE RELATED COMPOUND B	30MG	221176-49-4
PHR2538-30MG	OLANZAPINE RELATED COMPOUND C	30MG	174794-02-6
OLMESARTAN			
PHR1851-200MG	OLMESARTAN MEDOXOMIL	200MG	144689-63-4
PHR2016-40MG	OLMESARTAN MEDOXOMIL RELATED COMPOUND A	40MG	849206-43-5
OLOPATADINE			
PHR3069-300MG	OLOPATADINE HYDROCHLORIDE	300MG	140462-76-6
PHR2539-50MG	OLOPATADINE RELATED COMPOUND C	50MG	
OMEPRAZOLE			
PHR1585-1G	ESOMEPRAZOLE MAGNESIUM	1G	217087-09-7
PHR1059-1G	OMEPRAZOLE	1G	73590-58-6
PHR1987-30MG	OMEPRAZOLE IMPURITY B	30MG	110374-16-8
PHR1988-30MG	OMEPRAZOLE IMPURITY C	30MG	73590-85-9
PHR1989-30MG	OMEPRAZOLE IMPURITY H	30MG	863029-89-4
PHR1648-30MG	OMEPRAZOLE RELATED COMPOUND A	30MG	88546-55-8
PHR1986-30MG	OMEPRAZOLE RELATED COMPOUND B	30MG	37052-78-1
PHR1649-30MG	OMEPRAZOLE RELATED COMPOUND E	30MG	176219-04-8
PHR1650-30MG	OMEPRAZOLE RELATED COMPOUND I	30MG	
PHR1954-30MG	N- METHYL OMEPRAZOLE	30MG	
ONDANSETRON			
PHR1141-1G	ONDANSETRON HYDROCHLORIDE	1G	103639-04-9
PHR2223-100MG	ONDANSETRON RELATED COMPOUND A	100MG	119812-29-2
PHR2224-100MG	ONDANSETRON RELATED COMPOUND C	100MG	27387-31-1
PHR2221-100MG	ONDANSETRON RELATED COMPOUND D	100MG	99614-64-9
ORLISTAT			
PHR1445-1G	ORLISTAT	1G	96829-58-2
PHR2543-50MG	ORLISTAT RELATED COMPOUND B	50MG	19740-72-8
PHR2544-50MG	ORLISTAT RELATED COMPOUND C	50MG	791-28-6
ORPHENADRINE			
PHR1777-1G	ORPHENADRINE CITRATE	1G	4682-36-4
PHR2957-40MG	ORPHENADRINE RELATED COMPOUND B	40MG	93940-17-1
PHR2958-40MG	ORPHENADRINE RELATED COMPOUND C	40MG	10488-36-5
PHR2959-100MG	ORPHENADRINE RELATED COMPOUND E	100MG	21945-86-8

Cat. No.	Description	Pack Size	CAS No.
OXACILLIN			
PHR2496-50MG	OXACILLIN RELATED COMPOUND C	50MG	1136-45-4
PHR2488-300MG	OXACILLIN SODIUM	300MG	7240-38-2
OXALIPLATIN			
PHR1528-200MG	OXALIPLATIN	200MG	61825-94-3
PHR2545-200MG	OXALIPLATIN RELATED COMPOUND A	200MG	6153-56-6
OXCARBAZEPINE			
PHR1635-1G	OXCARBAZEPINE	1G	28721-07-5
PHR3628-50MG	OXCARBAZEPINE RELATED COMPOUND A	50MG	1346601-76-0
PHR3071-25MG	OXCARBAZEPINE RELATED COMPOUND C	25MG	578-95-0
PHR2238-50MG	OXCARBAZEPINE RELATED COMPOUND E	50MG	21737-58-6
OXICONAZOLE			
PHR3736-40MG	OXICONAZOLE RELATED COMPOUND C	40MG	98164-08-0
OXYBUTYNIN			
PHR2547-500MG	OXYBUTYNIN CHLORIDE	500MG	1508-65-2
PHR2548-200MG	OXYBUTYNIN RELATED COMPOUND A	200MG	4335-77-7
PHR2549-50MG	OXYBUTYNIN RELATED COMPOUND B	50MG	10399-13-0
OXYMETAZOLIN			
PHR1572-500MG	OXYMETAZOLINE HYDROCHLORIDE	500MG	2315-02-8
PHR2215-100MG	OXYMETAZOLINE RELATED COMPOUND A	100MG	1391053-50-1
PACLITAXEL			
PHR1803-200MG	PACLITAXEL	200MG	33069-62-4
PHR2111-50MG	PACLITAXEL RELATED COMPOUND A	50MG	71610-00-9
PHR2112-50MG	PACLITAXEL RELATED COMPOUND B	50MG	78454-17-8
PHR9150-50MG	PACLITAXEL IMPURITY K	50MG	148930-55-6
PHR9149-50MG	7-EPI-PACLITAXEL	50MG	105454-04-4
PHR9146-50MG	7-O-ACETYLPACLITAXEL	50MG	92950-39-5
PALIPERIDONE			
PHR3072-50MG	PALIPERIDONE RELATED COMPOUND B	50MG	84163-13-3
PANTOPRAZOLE			
PHR1604-1G	PANTOPRAZOLE SODIUM	1G	164579-32-2
PHR2458-50MG	PANTOPRAZOLE RELATED COMPOUND A	50MG	127780-16-9
PHR3074-50MG	PANTOPRAZOLE RELATED COMPOUND B	50MG	102625-64-9
PHR3075-50MG	PANTOPRAZOLE RELATED COMPOUND C	50MG	97963-62-7
PAROXETINE			
PHR1804-500MG	PAROXETINE HYDROCHLORIDE	500MG	110429-35-1
PHR3076-50MG	PAROXETINE RELATED COMPOUND C	50MG	
PHR3077-50MG	PAROXETINE RELATED COMPOUND F	50MG	110429-36-2
PERINDOPRIL			
PHR1887-500MG	PERINDOPRIL ERBUMINE	500MG	107133-36-8
PHR3078-50MG	PERINDOPRIL RELATED COMPOUND A	50MG	80875-98-5
PHR3859-40MG	PERINDOPRIL RELATED COMPOUND B	40MG	1418095-48-3
PERPHENAZINE			
PHR3079-500MG	PERPHENAZINE	500MG	58-39-9
PHR3718-50MG	PERPHENAZINE RELATED COMPOUND B	50MG	3533-97-9
PHR3729-50MG	PERPHENAZINE SULFOXIDE	50MG	10078-25-8
PHENIRAMINE			
PHR2719-50MG	PHENIRAMINE IMPURITY A	50MG	101-82-6
PHR2646-300MG	PHENIRAMINE MALEATE	300MG	132-20-7
PHENYLBUTYRATE			
PHR3080-50MG	PHENYLBUTYRATE RELATED COMPOUND A	50MG	2051-95-8
PHR3081-200MG	PHENYLBUTYRATE RELATED COMPOUND B	200MG	529-34-0
PHR3082-50MG	PHENYLBUTYRATE RELATED COMPOUND C	50MG	4441-63-8
PHR3237-500MG	SODIUM PHENYLBUTYRATE	500MG	1716-12-7
PHENYLEPHRINE			
PHR1017-500MG	PHENYLEPHRINE HYDROCHLORIDE	500MG	61-76-7
PHR2151-100MG	PHENYLEPHRINE IMPURITY 01	100MG	57196-61-9

Cat. No.	Description	Pack Size	CAS No.
PHENYLEPHRINE			
PHR2153-100MG	PHENYLEPHRINE IMPURITY 03	100MG	
PHR2154-100MG	PHENYLEPHRINE IMPURITY 04	100MG	
PHR1695-50MG	PHENYLEPHRINE RELATED COMPOUND C	50MG	94240-17-2
PHR2148-100MG	PHENYLEPHRINE RELATED COMPOUND D	100MG	1367567-95-0
PHR1712-50MG	PHENYLEPHRINE RELATED COMPOUND E	50MG	71786-67-9
PHR2149-100MG	PHENYLEPHRINE RELATED COMPOUND F	100MG	
PHR2150-100MG	PHENYLEPHRINE RELATED COMPOUND G	100MG	1094089-46-9
PHENYTOIN			
PHR1139-1G	PHENYTOIN	1G	57-41-0
PHR1492-1G	PHENYTOIN SODIUM	1G	630-93-3
PHR1692-100MG	PHENYTOIN RELATED COMPOUND B	100MG	6802-95-5
PHR2553-100MG	PHENYTOIN IMPURITY F	100MG	
PHR2246-100MG	PHENYTOIN INTERMEDIATE COMPOUND M01	100MG	6859-11-6
PIPERACILLIN			
PHR1805-1G	PIPERACILLIN	1G	66258-76-2
PHR3085-50MG	PIPERACILLIN RELATED COMPOUND E	50MG	59702-31-7
PIROCTONE			
PHR2187-500MG	PIROCTONE OLAMINE	500MG	68890-66-4
PHR3086-50MG	PIROCTONE RELATED COMPOUND B	50MG	50650-75-4
PLASTIC ADDITIVE			
PHR2173-1G	PLASTIC ADDITIVE 2	1G	6683-19-8
PHR3090-400MG	PLASTIC ADDITIVE 3	400MG	1709-70-2
PHR2174-1G	PLASTIC ADDITIVE 4	1G	2082-79-3
PHR3091-400MG	PLASTIC ADDITIVE 5	400MG	31570-04-4
PHR3092-400MG	PLASTIC ADDITIVE 6	400MG	27676-62-6
PHR3093-200MG	PLASTIC ADDITIVE 8	200MG	2500-88-1
PHR3094-200MG	PLASTIC ADDITIVE 9	200MG	123-28-4
PHR3095-200MG	PLASTIC ADDITIVE 10	200MG	693-36-7
PHR3096-200MG	PLASTIC ADDITIVE 11	200MG	65447-77-0
PHR3088-200MG	PLASTIC ADDITIVE 13	200MG	112-84-5
PHR3089-200MG	PLASTIC ADDITIVE 14	200MG	117-81-7
PHR2965-200MG	PLASTIC ADDITIVE 16	200MG	8016-11-3
POLYETHYLENE GLYCOL			
PHR3335-1G	POLYETHYLENE GLYCOL 200	1G	25322-68-3
PHR3336-1G	POLYETHYLENE GLYCOL 300	1G	25322-68-3
PHR2891-1G	POLYETHYLENE GLYCOL 400	1G	25322-68-3
PHR3338-1G	POLYETHYLENE GLYCOL 600	1G	25322-68-3
PHR2890-1G	POLYETHYLENE GLYCOL 1000	1G	25322-68-3
PHR3334-1G	POLYETHYLENE GLYCOL 1500	1G	25322-68-3
PHR3337-1G	POLYETHYLENE GLYCOL 3000	1G	25322-68-3
PHR2362-1G	POLYETHYLENE GLYCOL 3350	1G	25322-68-3
PHR2892-1G	POLYETHYLENE GLYCOL 4000	1G	25322-68-3
PHR2893-1G	POLYETHYLENE GLYCOL 6000	1G	25322-68-3
PHR2894-1G	POLYETHYLENE GLYCOL 8000	1G	25322-68-3
PRAMIPEXOLE			
PHR1598-500MG	PRAMIPEXOLE DIHYDROCHLORIDE	500MG	191217-81-9
PHR2010-50MG	PRAMIPEXOLE RELATED COMPOUND A	50MG	106092-09-5
PRAVASTATIN			
PHR3098-30MG	PRAVASTATIN RELATED COMPOUND A	30MG	81093-43-8
PHR3407-100MG	PRAVASTATIN SODIUM	100MG	81131-70-6
PHR3406-400MG	PRAVASTATIN 1,1,3,3-TETRAMETHYLBUTYLAMINE	400MG	151006-14-3
PRAZQUANTEL			
PHR1391-1G	PRAZQUANTEL	1G	55268-74-1
PHR2961-40MG	PRAZQUANTEL RELATED COMPOUND A	40MG	54761-87-4
PHR2962-100MG	PRAZQUANTEL RELATED COMPOUND B	100MG	125273-86-1

Cat. No.	Description	Pack Size	CAS No.
PRAZIQUANTEL			
PHR2963-100MG	PRAZIQUANTEL RELATED COMPOUND C	100MG	125273-88-3
PRAZOSIN			
PHR3099-500MG	PRAZOSIN HYDROCHLORIDE	500MG	19237-84-4
PHR3100-50MG	PRAZOSIN RELATED COMPOUND D	50MG	40172-95-0
PREDNICARBATE			
PHR3101-500MG	PREDNICARBATE	500MG	73771-04-7
PHR3102-50MG	PREDNICARBATE RELATED COMPOUND C	50MG	5740-62-5
PREDNISOLONE			
PHR1043-500MG	PREDNISOLONE	500MG	50-24-8
PHR1630-1G	PREDNISOLONE ACETATE	1G	52-21-1
PHR2328-30MG	PREDNISOLONE IMPURITY G	30MG	
PHR2556-50MG	PREDNISOLONE IMPURITY H	50MG	2427-64-7
PHR2557-50MG	PREDNISOLONE IMPURITY I	50MG	13479-38-4
PHR3227-500MG	PREDNISOLONE TEBUTATE	500MG	7681-14-3
PHR2816-500MG	PREDNISOLONE SODIUM PHOSPHATE	500MG	125-02-0
PHR3103-50MG	11-EPI-PREDNISOLONE	50MG	600-90-8
PREDNISON			
PHR1042-1G	PREDNISON	1G	53-03-2
PHR2288-20MG	PREDNISON IMPURITY C	20MG	70522-55-3
PHR2289-50MG	PREDNISON IMPURITY J	50MG	1807-14-3
PREGABALIN			
PHR3104-50MG	PREGABALIN RELATED COMPOUND C	500MG	61312-87-6
PRILOCAINE			
PHR3105-500MG	PRILOCAINE	500MG	721-50-6
PHR3106-500MG	PRILOCAINE HYDROCHLORIDE	500MG	1786-81-8
PHR3107-200MG	PRILOCAINE RELATED COMPOUND A	200MG	636-21-5
PRIMIDONE			
PHR3109-500MG	PRIMIDONE	500MG	125-33-7
PHR3110-50MG	PRIMIDONE RELATED COMPOUND A	50MG	7206-76-0
PHR3111-50MG	PRIMIDONE RELATED COMPOUND C	50MG	90-26-6
PROBUCOL			
PHR3112-500MG	PROBUCOL	500MG	23288-49-5
PHR3113-50MG	PROBUCOL RELATED COMPOUND A	50MG	23288-49-5
PROCHLORPERAZINE			
PHR3114-500MG	PROCHLORPERAZINE MALEATE	500MG	84-02-6
PHR3786-50MG	PROCHLORPERAZINE RELATED COMPOUND A	50MG	1346602-34-3
PROGESTERONE			
PHR1142-1G	PROGESTERONE	1G	57-83-0
PHR2280-50MG	PROGESTERONE IMPURITY J	50MG	1162-54-5
PHR2905-100MG	PROGESTERONE RELATED COMPOUND H	100MG	1162-56-7
PHR2563-100MG	11-KETOPROGESTERONE	100MG	516-15-4
PHR2562-100MG	11A-HYDROXYPROGESTERONE	100MG	80-75-1
PHR2628-300MG	HYDROXYPROGESTERONE CAPROATE	300MG	630-56-8
PROGUANIL			
PHR1713-1G	PROGUANIL HYDROCHLORIDE	1G	637-32-1
PHR3115-50MG	PROGUANIL RELATED COMPOUND A	50MG	44830-55-9
PHR3116-50MG	PROGUANIL RELATED COMPOUND C	50MG	13590-98-2
PHR3117-50MG	PROGUANIL RELATED COMPOUND E	50MG	1482-62-8
PROMETHAZINE			
PHR1467-1G	PROMETHAZINE HYDROCHLORIDE	1G	58-33-3
PHR2777-50MG	PROMETHAZINE RELATED COMPOUND B	50MG	5568-90-1
PROPAFENONE			
PHR3566-50MG	PROPAFENONE RELATED COMPOUND A	50MG	
PROPOFOL			
PHR1663-1ML	PROPOFOL	1ML	2078-54-8

Cat. No.	Description	Pack Size	CAS No.
PROPOFOL			
PHR2565-50MG	PROPOFOL RELATED COMPOUND A	50MG	2416-95-7
PHR2566-50MG	PROPOFOL RELATED COMPOUND B	50MG	1988-11-0
PHR9062-25MG	PROPOFOL EP IMPURITY K	25MG	14366-59-7
PROPRANOLOL			
PHR1308-1G	PROPRANOLOL HYDROCHLORIDE	1G	318-98-9
PHR3119-50MG	PROPRANOLOL RELATED COMPOUND A	50MG	36112-95-5
PROPYL GALLATE			
PHR1118-1G	PROPYL GALLATE	1G	121-79-9
PHR2318-100MG	PROPYL GALLATE IMPURITY A	100MG	149-91-7
PYRANTEL			
PHR2080-1G	PYRANTEL PAMOATE	1G	22204-24-6
PHR2082-50MG	PYRANTEL IMPURITY B	50MG	36700-39-7
PHR2083-50MG	PYRANTEL IMPURITY C	50MG	98-03-3
PHR2084-50MG	PYRANTEL IMPURITY D	50MG	4271-96-9
PHR2081-50MG	PYRANTEL RELATED COMPOUND A	50MG	
PYRIDOXINE			
PHR9291-100MG	PYRIDOXINE	100MG	65-23-6
PHR2213-100MG	BIS-PYRIDOXINE	100MG	19203-56-6
PHR1036-500MG	PYRIDOXINE HYDROCHLORIDE	500MG	58-56-0
PHR1693-50MG	PYRIDOXINE IMPURITY A	50MG	5196-20-3
PHR1694-50MG	PYRIDOXINE IMPURITY B	50MG	
QUETIAPINE			
PHR1856-1G	QUETIAPINE FUMARATE	1G	111974-72-2
PHR3121-40MG	QUETIAPINE RELATED COMPOUND B	40MG	111974-74-4
PHR3122-40MG	QUETIAPINE RELATED COMPOUND G	40MG	3159-07-7
PHR8631-50MG	QUETIAPINE N-OXIDE	50MG	1076199-40-0
RABEPRAZOLE			
PHR1525-1G	RABEPRAZOLE SODIUM	1G	117976-90-6
PHR3752-50MG	RABEPRAZOLE RELATED COMPOUND B	50MG	924663-38-7
PHR3123-50MG	RABEPRAZOLE RELATED COMPOUND C	50MG	583-39-1
PHR3677-50MG	RABEPRAZOLE RELATED COMPOUND D	50MG	117976-47-3
PHR3124-50MG	RABEPRAZOLE RELATED COMPOUND E	50MG	117977-21-6
PHR3125-30MG	RABEPRAZOLE RELATED COMPOUND F	30MG	168167-42-8
PHR3788-100MG	RABEPRAZOLE IMPURITY I	100MG	924663-37-6
RALOXIFENE			
PHR1852-500MG	RALOXIFENE HYDROCHLORIDE	500MG	82640-04-8
PHR2015-40MG	RALOXIFENE RELATED COMPOUND C	40MG	195454-31-0
RALTEGRAVIR			
PHR3126-500MG	RALTEGRAVIR POTASSIUM	500MG	871038-72-1
PHR3127-50MG	RALTEGRAVIR RELATED COMPOUND E	50MG	1193687-87-4
RAMIPRIL			
PHR1446-1G	RAMIPRIL	1G	87333-19-5
PHR8633-50MG	RAMIPRIL DIKETOPIPERAZINE	50MG	108731-95-9
PHR3128-50MG	RAMIPRIL RELATED COMPOUND A	50MG	108313-11-7
PHR3339-50MG	RAMIPRIL RELATED COMPOUND B	50MG	295328-72-2
RANITIDINE			
PHR1026-500MG	RANITIDINE HYDROCHLORIDE	500MG	66357-59-3
PHR1773-100MG	RANITIDINE RELATED COMPOUND A	100MG	91224-69-0
PHR1774-100MG	RANITIDINE RELATED COMPOUND B	100MG	72126-78-4
PHR1775-100MG	RANITIDINE RELATED COMPOUND C	100MG	73851-70-4
PHR2175-100MG	RANITIDINE IMPURITY D	100MG	112251-56-6
PHR2177-100MG	RANITIDINE IMPURITY F	100MG	
PHR2178-100MG	RANITIDINE IMPURITY G	100MG	112233-23-5
PHR2179-100MG	RANITIDINE IMPURITY H	100MG	72078-82-1
PHR2180-100MG	RANITIDINE IMPURITY I	100MG	

Cat. No.	Description	Pack Size	CAS No.
RANITIDINE			
PHR2181-100MG	RANITIDINE IMPURITY J	100MG	1331637-48-9
PHR2182-100MG	RANITIDINE IMPURITY K	100MG	61832-41-5
RANOLAZINE			
PHR2676-500MG	RANOLAZINE	500MG	95635-55-5
PHR2677-100MG	RANOLAZINE RELATED COMPOUND A	100MG	2210-74-4
PHR2678-100MG	RANOLAZINE RELATED COMPOUND B	100MG	
PHR2679-100MG	RANOLAZINE RELATED COMPOUND C	100MG	5294-61-1
PHR2680-100MG	RANOLAZINE RELATED COMPOUND D	100MG	380204-72-8
REPAGLINIDE			
PHR3408-400MG	REPAGLINIDE	400MG	135062-02-1
PHR3129-25MG	REPAGLINIDE RELATED COMPOUND A	25MG	219921-94-5
PHR3340-100MG	REPAGLINIDE RELATED COMPOUND B	50MG	99469-99-5
PHR3130-50MG	REPAGLINIDE RELATED COMPOUND C	50MG	1803084-79-8
PHR3131-50MG	REPAGLINIDE RELATED COMPOUND E	50MG	147852-26-4
RETINYL			
PHR1236-1G	RETINYL ACETATE	1G	127-47-9
PHR1235-1G	RETINYL PALMITATE	1G	79-81-2
PHR1668-1G	RETINYL PROPIONATE	1G	7069-42-3
RIBAVIRIN			
PHR1700-1G	RIBAVIRIN	1G	36791-04-5
PHR2779-30MG	RIBAVIRIN RELATED COMPOUND D	30MG	3641-08-5
RIBOFLAVIN			
PHR1054-1G	RIBOFLAVIN	1G	83-88-5
PHR1297-1G	PHOSPHATED RIBOFLAVIN (RIBOFLAVIN SODIUM PHOSPHATE)	1G	
PHR1936-10MG	RIBOFLAVIN IMPURITY A	10MG	
PHR2570-50MG	RIBOFLAVIN IMPURITY B	50MG	
RILUZOLE			
PHR2499-300MG	RILUZOLE	300MG	1744-22-5
PHR3132-50MG	RILUZOLE RELATED COMPOUND A	50MG	461-82-5
RISEDRONATE			
PHR1888-1G	RISEDRONATE SODIUM	1G	329003-65-8
PHR2078-50MG	RISEDRONATE IMPURITY C	50MG	105462-25-7
PHR2079-50MG	RISEDRONATE IMPURITY D	50MG	6419-36-9
PHR2075-50MG	RISEDRONATE RELATED COMPOUND A	50MG	105462-23-5
PHR2077-50MG	RISEDRONATE RELATED COMPOUND C	50MG	75755-10-1
RIVASTIGMINE			
PHR3343-0.2ML	RIVASTIGMINE	0.2ML	123441-03-2
PHR2029-30MG	RIVASTIGMINE RELATED COMPOUND A	30MG	32634-68-7
PHR2030-30MG	RIVASTIGMINE RELATED COMPOUND B	30MG	25081-93-0
PHR2031-30MG	RIVASTIGMINE RELATED COMPOUND C	30MG	139306-10-8
PHR2032-100MG	RIVASTIGMINE RELATED COMPOUND D	100MG	855300-09-3
PHR2033-30MG	RIVASTIGMINE RELATED COMPOUND E	30MG	
PHR2034-50MG	RIVASTIGMINE RELATED COMPOUND F	50MG	1346602-84-3
PHR1867-1G	RIVASTIGMINE TARTRATE	1G	129101-54-8
ROPINIROLE			
PHR1633-200MG	ROPINIROLE HYDROCHLORIDE	200MG	91374-20-8
PHR3345-100MG	ROPINIROLE RELATED COMPOUND A	100MG	139122-19-3
PHR3133-100MG	ROPINIROLE RELATED COMPOUND B	100MG	221264-21-7
ROPIVACAINE			
PHR2571-500MG	ROPIVACAINE	500MG	84057-95-4
PHR2303-500MG	ROPIVACAINE HYDROCHLORIDE	500MG	132112-35-7
PHR2572-50MG	ROPIVACAINE RELATED COMPOUND A	50MG	21436-98-6
PHR2306-100MG	ROPIVACAINE RELATED COMPOUND B/ROPIVACAI	100MG	
ROSUVASTATIN			
PHR1928-1G	ROSUVASTATIN CALCIUM	1G	147098-20-2

Cat. No.	Description	Pack Size	CAS No.
ROSUVASTATIN			
PHR3678-50MG	ROSUVASTATIN IMPURITY D	50MG	503610-43-3
PHR3346-50MG	ROSUVASTATIN RELATED COMPOUND C	50MG	
ROTIGOTINE			
PHR3725-20MG	ROTIGOTINE RELATED COMPOUND H	20MG	1148154-91-9
PHR3745-50MG	ROTIGOTINE RELATED COMPOUND K	50MG	51927-48-1
SAFINAMIDE			
PHR2379-50MG	SAFINAMIDE MESYLATE	50MG	202825-46-5
PHR2381-20MG	SAFINAMIDE DEFLUOR-DERIVATIVE	20MG	
PHR2383-20MG	SAFINAMIDE DIALANINE AMIDE ADDUCT	20MG	
PHR2380-20MG	SAFINAMIDE FREE ACID	20MG	
PHR2382-20MG	SAFINAMIDE IMIDAZOLIDINONE	20MG	
SALBUTAMOL			
PHR1959-50MG	SALBUTAMOL IMPURITY B	50MG	96948-64-0
PHR1960-50MG	SALBUTAMOL IMPURITY D	50MG	156339-88-7
PHR1961-50MG	SALBUTAMOL IMPURITY G	50MG	
PHR1962-20MG	SALBUTAMOL IMPURITY I	20MG	56796-66-8
PHR2939-50MG	SALBUTAMOL IMPURITY K	50MG	898542-80-8
SALICYLIC ACID			
PHR1940-50MG	ACETYLSALICYLSALICYLIC ACID	50MG	530-75-6
PHR9138-100MG	3-AMINOSALICYLIC ACID	100MG	570-23-0
PHR1746-500MG	4-AMINOSALICYLIC ACID	500MG	65-49-6
PHR9052-100MG	3-NITROSALICYLIC ACID	100MG	85-38-1
PHR1643-200MG	SALICYLIC ACID RELATED COMPOUND B	200MG	636-46-4
SERTRALINE			
PHR2535-200MG	SERTRALINE HYDROCHLORIDE	200MG	79559-97-0
PHR3852-50MG	SERTRALINE HYDROCHLORIDE(RACEMIC MIXTURE)	50MG	
PHR3397-50MG	SERTRALINE HYDROCHLORIDE RELATED COMPOUND A	50MG	79617-99-5
SESAME OIL			
PHR2967-1G	SESAME OIL	1G	8008-74-0
PHR2968-50MG	SESAME OIL RELATED COMPOUND B	50MG	2190-15-0
SEVOFLURANE			
PHR3348-1ML	SEVOFLURANE	1ML	28523-86-6
PHR3347-0.5ML	SEVOFLURANE RELATED COMPOUND B	0.5ML	13171-18-1
SILDENAFIL			
PHR1807-500MG	SILDENAFIL CITRATE	500MG	171599-83-0
PHR3682-50MG	SILDENAFIL RELATED COMPOUND A	50MG	1391053-95-4
SITAGLIPTIN			
PHR3270-40MG	SITAGLIPTIN IMPURITY B (PH. EUR.)	40MG	1345822-87-8
PHR3271-40MG	SITAGLIPTIN IMPURITY A	40MG	823817-58-9
PHR3269-40MG	SITAGLIPTIN IMPURITY C	40MG	1345822-86-7
PHR1857-1G	SITAGLIPTIN PHOSPHATE	1G	654671-77-9
SORAFENIB			
PHR3687-50MG	SORAFENIB RELATED COMPOUND H	50MG	1431697-81-2
SOTALOL			
PHR2798-500MG	SOTALOL HYDROCHLORIDE	500MG	959-24-0
PHR3353-100MG	SOTALOL RELATED COMPOUND B	100MG	83922-54-7
SPIRONOLACTONE			
PHR2886-200MG	SPIRONOLACTONE	200MG	52-01-7
PHR2887-150MG	SPIRONOLACTONE RELATED COMPOUND A	150MG	976-71-6
PHR3135-50MG	SPIRONOLACTONE RELATED COMPOUND C	50MG	976-70-5
SULBACTAM			
PHR2576-500MG	SULBACTAM	500MG	68373-14-8
PHR3684-40MG	SULBACTAM RELATED COMPOUND E / EP IMPURITY E	40MG	76646-91-8
PHR3685-40MG	SULBACTAM RELATED COMPOUND F / EP IMPURITY F	40MG	24158-88-1
SULCONAZOLE			
PHR3750-50MG	SULCONAZOLE RELATED COMPOUND A	50MG	77279-29-9

Cat. No.	Description	Pack Size	CAS No.
SULFAMETHOXAZOLE			
PHR1126-1G	SULFAMETHOXAZOLE	1G	723-46-6
PHR2578-30MG	SULFAMETHOXAZOLE IMPURITY A	30MG	21312-10-7
PHR3686-40MG	SULFAMETHOXAZOLE RELATED COMPOUND B	40G	135529-16-7
PHR2577-100MG	SULFAMETHOXAZOLE RELATED COMPOUND C	100MG	1072-67-9
SUMATRIPTAN			
PHR2579-100MG	SUMATRIPTAN	100MG	103628-46-2
PHR2580-500MG	SUMATRIPTAN SUCCINATE	500MG	103628-48-4
PHR2914-75MG	SUMATRIPTAN SUCCINATE RELATED COMPOUND C	75MG	
TACROLIMUS			
PHR1809-150MG	TACROLIMUS	150MG	109581-93-3
TENIPOSIDE			
PHR3362-20MG	TENIPOSIDE RELATED COMPOUND A	20MG	6559-91-7
TENOFOVIR DISOPROXIL			
PHR1957-500MG	TENOFOVIR DISOPROXIL FUMARATE	500MG	202138-50-9
PHR3693-50MG	TENOFOVIR DISOPROXIL RELATED COMPOUND A	50MG	1432630-26-6
PHR3816-50MG	TENOFOVIR DISOPROXIL RELATED COMPOUND B	50MG	1446486-33-4
TERAZOSIN			
PHR1721-1G	TERAZOSIN HYDROCHLORIDE	1G	70024-40-7
PHR3137-100MG	TERAZOSIN RELATED COMPOUND A	100MG	
PHR3138-100MG	TERAZOSIN RELATED COMPOUND B	100MG	
PHR3139-50MG	TERAZOSIN RELATED COMPOUND C	50MG	1486464-41-8
TERBINAFINE			
PHR1298-1G	TERBINAFINE HYDROCHLORIDE	1G	78628-80-5
PHR2885-50MG	TERBINAFINE RELATED COMPOUND A	50MG	65473-13-4
PHR2927-50MG	TERBINAFINE RELATED COMPOUND B	50MG	176168-78-8
PHR3363-50MG	TERBINAFINE RELATED COMPOUND C	50MG	877265-30-0
PHR3364-50MG	TERBINAFINE RELATED COMPOUND D	50MG	877265-33-3
TERBUTALINE			
PHR3246-250MG	TERBUTALINE SULFATE	250MG	23031-32-5
PHR3784-50MG	TERBUTALINE RELATED COMPOUND D		94109-61-2
TETRACAINE			
PHR1129-1G	TETRACAINE HYDROCHLORIDE	1G	136-47-0
PHR2236-50MG	TETRACAINE RELATED COMPOUND B	50MG	4740-24-3
PHR2261-100MG	TETRACAINE RELATED COMPOUND C	100MG	71839-12-8
PHR3140-100MG	TETRACAINE	100MG	94-24-6
THEOPHYLLINE			
PHR2551-50MG	THEOPHYLLINE RELATED COMPOUND D	500MG	116131-08-9
THIAMINE			
PHR2169-100MG	ETHYLTHIAMINE HYDROCHLORIDE	100MG	
PHR1037-1G	THIAMINE HYDROCHLORIDE	1G	67-03-8
PHR2216-1G	THIAMINE DISULFIDE	1G	67-16-3
PHR1996-20MG	THIAMINE IMPURITY A	20MG	
PHR1997-20MG	THIAMINE IMPURITY B	20MG	
PHR1998-20MG	THIAMINE IMPURITY C	20MG	
PHR2168-100MG	THIAMINE IMPURITY D	100MG	490-82-4
PHR1759-20MG	THIAMINE IMPURITY E	20MG	299-35-4
PHR2170-100MG	THIAMINE IMPURITY G	100MG	3419-28-1
PHR2172-50MG	THIAMINE IMPURITY I	50MG	13100-57-7
PHR1337-1G	THIAMINE NITRATE	1G	532-43-4
PHR1369-500MG	THIAMINE PYROPHOSPHATE	500MG	154-87-0
PHR2207-100MG	THIAMINE THIOCARBAMATE	100MG	89285-03-0
THIMEROSAL			
PHR1587-1G	THIMEROSAL	1G	54-64-8
PHR2590-100MG	THIMEROSAL RELATED COMPOUND A	100MG	119-80-2

Cat. No.	Description	Pack Size	CAS No.
TIAMULIN			
PHR2721-500MG	TIAMULIN FUMARATE	500MG	55297-96-6
PHR3577-100MG	TIAMULIN RELATED COMPOUND A	100MG	31716-01-5
TIGECYCLINE			
PHR2591-500MG	TIGECYCLINE	500MG	220620-09-7
PHR2592-50MG	TIGECYCLINE RELATED COMPOUND B	50MG	149934-21-4
TIMOLOL			
PHR2593-500MG	TIMOLOL MALEATE	500MG	26921-17-5
PHR3144-50MG	TIMOLOL RELATED COMPOUND A	50MG	26839-77-0
PHR3695-40MG	TIMOLOL RELATED COMPOUND D	40MG	30165-97-0
PHR3617-50MG	TIMOLOL RELATED COMPOUND B	50MG	59697-06-2
PHR2594-50MG	TIMOLOL RELATED COMPOUND F	50MG	30165-96-9
TINIDAZOLE			
PHR2807-400MG	TINIDAZOLE	400MG	19387-91-8
PHR1195-500MG	TINIDAZOLE RELATED COMPOUND A	500MG	696-23-1
PHR3696-40MG	TINIDAZOLE RELATED COMPOUND B	40MG	25459-12-5
TIOCONAZOLE			
PHR2808-400MG	TIOCONAZOLE	400MG	65899-73-2
PHR3756-50MG	TIOCONAZOLE RELATED COMPOUND B	50MG	61675-62-5
TIZANIDINE			
PHR3814-100MG	TIZANIDINE RELATED COMPOUND B	100MG	173532-15-5
PHR3694-100MG	TIZANIDINE RELATED COMPOUND C	100MG	5391-52-6
TOPOTECAN			
PHR3369-50MG	TOPOTECAN RELATED COMPOUND C	50MG	7689-03-4
TORSEMIDE			
PHR3370-500MG	TORSEMIDE	500MG	56211-40-6
PHR3145-150MG	TORSEMIDE RELATED COMPOUND A	150MG	72811-73-5
PHR3550-100MG	TORSEMIDE RELATED COMPOUND B	100MG	160972-33-8
PHR3600-50MG	TORSEMIDE RELATED COMPOUND C	50MG	58155-35-4
PHR3544-50MG	TORSEMIDE RELATED COMPOUND E	50MG	72810-61-8
TRANDOLAPRIL			
PHR3146-250MG	TRANDOLAPRIL	250MG	87679-37-6
PHR3601-50MG	TRANDOLAPRIL RELATED COMPOUND D	50MG	149881-40-3
TRANEXAMIC ACID			
PHR1812-500MG	TRANEXAMIC ACID	500MG	1197-18-8
PHR2920-30MG	TRANEXAMIC ACID RELATED COMPOUND C	30MG	330838-52-3
TRANLYCYPROMINE			
PHR3427-250MG	TRANLYCYPROMINE SULFATE	250MG	13492-01-8
PHR3147-100MG	TRANLYCYPROMINE RELATED COMPOUND A	100MG	54779-58-7
TRAVOPROST			
PHR1622-3ML	TRAVOPROST	3ML	157283-68-6
TRAZODONE			
PHR3372-500MG	TRAZODONE HYDROCHLORIDE	500MG	25332-39-2
PHR3574-50MG	TRAZODONE RELATED COMPOUND C	50MG	25332-39-2
PHR3148-50MG	TRAZODONE RELATED COMPOUND D	50MG	1263278-80-3
PHR3373-50MG	TRAZODONE RELATED COMPOUND F	50MG	52605-52-4
TRIAMCINOLONE			
PHR3253-500MG	TRIAMCINOLONE	500MG	124-94-7
PHR1701-1G	TRIAMCINOLONE ACETONIDE	1G	76-25-5
PHR3254-500MG	TRIAMCINOLONE HEXACETONIDE	500MG	5611-51-8
PHR3374-30MG	TRIAMCINOLONE ACETONIDE RELATED COMPOUND C	30G	161740-69-8
TRIAMTERENE			
PHR1722-400MG	TRIAMTERENE	400MG	396-01-1
PHR2014-100MG	TRIAMTERENE RELATED COMPOUND A	100MG	1006-23-1
PHR1820-100MG	TRIAMTERENE RELATED COMPOUND C	100MG	19152-93-3

Cat. No.	Description	Pack Size	CAS No.
TRICLOCARBAN			
PHR1303-500MG	TRICLOCARBAN	500MG	101-20-2
PHR3375-30MG	TRICLOCARBAN RELATED COMPOUND A	30MG	1219-99-4
PHR3149-40MG	TRICLOCARBAN RELATED COMPOUND C	40MG	4300-43-0
TRIHEXYPHENIDYL			
PHR3255-500MG	TRIHEXYPHENIDYL HYDROCHLORIDE	500MG	52-49-3
PHR3256-50MG	TRIHEXYPHENIDYL RELATED COMPOUND A	50MG	886-06-6
TRIMETHOPRIM			
PHR1056-1G	TRIMETHOPRIM	1G	738-70-5
PHR3603-50MG	TRIMETHOPRIM RELATED COMPOUND A	50MG	60729-91-1
PHR3602-50MG	TRIMETHOPRIM RELATED COMPOUND B	50MG	30806-86-1
TRIMIPRAMINE			
PHR3151-200MG	TRIMIPRAMINE MALEATE	200MG	521-78-8
PHR3578-50MG	TRIMIPRAMINE RELATED COMPOUND A	50MG	315-69-5
TROPICAMIDE			
PHR3152-500MG	TROPICAMIDE	500MG	1508-75-4
PHR3376-200MG	TROPICAMIDE RELATED COMPOUND A	200MG	33403-97-3
TROSPIMUM CHLORIDE			
PHR3378-200MG	TROSPIMUM CHLORIDE	200MG	10405-02-4
PHR3154-50MG	TROSPIMUM CHLORIDE RELATED COMPOUND A	50MG	76-93-7
PHR3155-50MG	TROSPIMUM CHLORIDE RELATED COMPOUND B	50MG	63516-30-3
TRYPTOPHAN			
PHR1176-1G	TRYPTOPHAN	1G	73-22-3
PHR2602-30MG	TRYPTOPHAN RELATED COMPOUND A	30MG	132685-02-0
UREA			
PHR2603-30MG	UREA RELATED COMPOUND A	30MG	"108-19-0"
VALACYCLOVIR			
PHR1601-1G	VALACYCLOVIR HYDROCHLORIDE	1G	124832-27-5
PHR3554-40MG	VALACYCLOVIR RELATED COMPOUND C	40MG	1346617-39-7
PHR3553-50MG	VALACYCLOVIR RELATED COMPOUND E	50MG	124832-31-1
PHR3618-50MG	VALACYCLOVIR RELATED COMPOUND F	50MG	86150-61-0
PHR2925-500MG	VALACYCLOVIR RELATED COMPOUND G	500MG	1122-58-3
PHR3749-50MG	D-VALACYCLOVIR	50MG	142963-60-8
VALPROIC ACID			
PHR1061-1G	VALPROIC ACID	1G	99-66-1
PHR1652-0.5ML	VALPROIC ACID RELATED COMPOUND A	0.5ML	99-67-2
PHR2281-100MG	VALPROIC ACID RELATED COMPOUND B	100MG	62391-99-5
PHR2614-100MG	VALPROIC ACID IMPURITY H	100MG	110-59-8
PHR2282-100MG	VALPROIC ACID IMPURITY K	100MG	5343-52-2
VALSARTAN			
PHR1315-1G	VALSARTAN	1G	137862-53-4
PHR1875-30MG	VALSARTAN RELATED COMPOUND A	30MG	
PHR1876-30MG	VALSARTAN RELATED COMPOUND B	30MG	952652-79-8
PHR1877-30MG	VALSARTAN RELATED COMPOUND C	30MG	
VENLAFAXINE			
PHR1736-1G	VENLAFAXINE HYDROCHLORIDE	1G	99300-78-4
PHR3697-40MG	VENLAFAXINE RELATED COMPOUND A	40MG	93413-90-2
VERAPAMIL			
PHR1131-1G	VERAPAMIL HYDROCHLORIDE	1G	152-11-4
PHR2616-100MG	VERAPAMIL IMPURITY K	100MG	20850-49-1
PHR8888-30MG	VERAPAMIL IMPURITY M		190850-50-1
PHR8890-100MG	VERAPAMIL IMPURITY P	100MG	2086275-13-8
PHR2360-100MG	VERAPAMIL RELATED COMPOUND A	100MG	67775-97-7
PHR2361-100MG	VERAPAMIL RELATED COMPOUND B	100MG	1794-55-4
VIGABATRIN			
PHR3429-400MG	VIGABATRIN	400MG	68506-86-5
PHR3158-100MG	VIGABATRIN RELATED COMPOUND A	100MG	7529-16-0

Cat. No.	Description	Pack Size	CAS No.
VINPOCETINE			
PHR3161-200MG	VINPOCETINE	200MG	42971-09-5
PHR3633-50MG	VINPOCETINE RELATED COMPOUND A	50MG	40163-56-2
PHR3634-50MG	VINPOCETINE RELATED COMPOUND B	50MG	4880-92-6
PHR3643-50MG	VINPOCETINE RELATED COMPOUND D	50MG	57327-92-1
VORICONAZOLE			
PHR1892-1G	VORICONAZOLE	1G	137234-62-9
PHR3698-40MG	VORICONAZOLE RELATED COMPOUND B	40MG	137234-63-0
ZALEPLON			
PHR3162-30MG	ZALEPLON RELATED COMPOUND A	30MG	
ZIDOVUDINE			
PHR1292-1G	ZIDOVUDINE	1G	30516-87-1
PHR1345-250MG	ZIDOVUDINE RELATED COMPOUND C	250MG	65-71-4
PHR3163-100MG	ZIDOVUDINE RELATED COMPOUND D	100MG	50-89-5
ZILEUTON			
PHR2555-200MG	ZILEUTON	200MG	111406-87-2
PHR3164-100MG	ZILEUTON RELATED COMPOUND B	100MG	97978-07-9
PHR3382-100MG	ZILEUTON RELATED COMPOUND C	100MG	22720-75-8
ZINC OXIDE			
PHR2724-500MG	ZINC OXIDE	500MG	1314-13-2
ZIPRASIDONE			
PHR3165-30MG	ZIPRASIDONE RELATED COMPOUND A	30MG	"87691-88-1"
ZOLEDRONIC ACID			
PHR1893-1G	ZOLEDRONIC ACID	1G	165800-06-6
PHR2459-50MG	ZOLEDRONIC ACID RELATED COMPOUND A	50MG	22884-10-2
ZOLMITRIPTAN			
PHR3857-50MG	ZOLMITRIPTAN RELATED COMPOUND B	50MG	152305-23-2
PHR3700-50MG	ZOLMITRIPTAN RELATED COMPOUND E	50MG	251451-30-6

Cat. No.	Description	Pack Size	CAS No.
Additional Standards			
PHR1963-50MG	(+)-CATECHIN	50MG	154-23-4
PHR2197-100MG	(2R)-2-AMINO-2-(4-HYDROXYPHENYL)-ACETAMIDE	100MG	54397-23-8
PHR2104-100MG	2-(4-CHLOROPHENYL)-4-DIMETHYLAMINO BUTYRONITRILE HYDROCHLORIDE	100MG	854863-55-1
PHR2146-100MG	2-(4-ISOBUTYLPHENYL)-N-(1-PHENYLETHYL)PROPIONAMIDE	100MG	105959-56-6
PHR2141-50MG	2-(4-ISOBUTYLPHENYL)PROP-2-ENOIC ACID	50MG	
PHR2122-50MG	2-(4-ISOBUTYRYLPHENYL)PROPANAMIDE	50MG	
PHR2147-50MG	2-(4-SEC-BUTYLPHENYL)PROP-2-ENOIC ACID	50MG	1368654-26-5
PHR1719-1G	2-AMINO-6-METHOXYPURINE	1G	20535-83-5
PHR1914-3X1.2ML	2-ETHYLHEXANOIC ACID	3X1.2ML	149-57-5
PHR2136-100MG	4-(1-HYDROXY-2-METHYLPROPYL)-ACETOPHENONE	100MG	
PHR2117-200MG	4,4'-DIFLUOROBENZOPHENONE	200MG	345-92-6
PHR2212-100MG	4-AMINO-5-AMINOMETHYL-2-METHYLPYRIMIDINE DIHYDROCHLORIDE	100MG	874-43-1
PHR1745-500MG	4-AMINOPHENOL HYDROCHLORIDE	500MG	51-78-5
PHR1916-3X1.2ML	4-HYDROXY-4-METHYL-2-PENTANONE	3X1.2ML	123-42-2
PHR2124-100MG	4-ISOBUTYRYLACETOPHENONE	100MG	103931-20-0

Cat. No.	Description	Pack Size	CAS No.
PHR1546-500MG	ACETYLCHOLINE CHLORIDE	500MG	60-31-1
PHR1599-1G	ALENDRONATE SODIUM	1G	121268-17-5
PHR2063-300MG	AMBROXOL HYDROCHLORIDE	300MG	23828-92-4
PHR1840-500MG	AMINOBENZOIC ACID	500MG	150-13-0
PHR2198-75MG	APIGENIN	75MG	520-36-5
PHR2199-75MG	APIGENIN-7-GLUCOSIDE	75MG	578-74-5
PHR2100-500MG	ARABINOSE	500MG	10323-20-3
PHR1909-500MG	ATENOLOL	500MG	29122-68-7
PHR1785-500MG	AZTREONAM	500MG	78110-38-0
PHR2097-500MG	BETAHISTINE HYDROCHLORIDE	500MG	5579-84-0
PHR1831-200MG	BROMHEXINE HYDROCHLORIDE	200MG	611-75-6
PHR1864-500MG	BROMPHENIRAMINE MALEATE	500MG	980-71-2
PHR1716-500MG	CAPREOMYCIN SULFATE	500MG	1405-37-4
PHR1684-1G	CARBIMAZOLE	1G	22232-54-8
PHR2208-20MG	CARNOSIC ACID	20MG	3650-09-7
PHR2099-1G	CEFOTAXIME SODIUM	1G	64485-93-4
PHR1847-500MG	CEFTAZIDIME PENTAHYDRATE	500MG	78439-06-2
PHR1848-1G	CEPHALEXIN	1G	15686-71-2
PHR2202-100MG	CHLOROGENIC ACID	100MG	327-97-9
PHR1786-500MG	CHONDROITIN SULFATE SODIUM	500MG	9082-07-9
PHR1921-500MG	CLOBETASOL PROPIONATE	500MG	25122-46-7
PHR1922-500MG	CLOXACILLIN SODIUM	500MG	7081-44-9
PHR1764-400MG	COLCHICINE	400MG	64-86-8
PHR1737-500MG	COLISTIMETHATE SODIUM	500MG	8068-28-8
PHR1605-1G	COLISTIN SULFATE	1G	1264-72-8
PHR1461-1G	COPPER CHLOROPHYLLIN	1G	11006-34-1
PHR2155-500MG	CORTISONE	500MG	53-06-5
PHR1923-300MG	CORTISONE ACETATE	300MG	50-04-4
PHR2209-50MG	CURCUMIN	50MG	458-37-7
PHR1765-1G	CYCLOSERINE	1G	68-41-7
PHR1787-500MG	CYTARABINE	500MG	147-94-4
PHR1858-100MG	DAIDZEIN	100MG	486-66-8
PHR1735-1G	DEMECLOCYCLINE HYDROCHLORIDE	1G	64-73-3
PHR1830-100MG	DESOGESTREL	100MG	54024-22-5
PHR1956-200MG	D-GLUCURONIC ACID	200MG	6556-12-3
PHR2094-5G	DIETHYLAMINO HYDROXYBENZOYL HEXYL BENZOATE	5G	302776-68-7
PHR1771-300MG	DIGOXIN	300MG	20830-75-5
PHR1517-500MG	DIHYDROSTREPTOMYCIN SULFATE	500MG	5490-27-7
PHR2118-500MG	DIMETHYL FUMARATE	500MG	624-49-7
PHR1690-1G	DISULFIRAM	1G	97-77-8
PHR1883-200MG	DOCETAXEL	200MG	148408-66-6
PHR1788-1G	DOCOSAHEXAENOIC ACID ETHYL ESTER	1G	84494-72-4
PHR1789-200MG	DOXORUBICIN HYDROCHLORIDE	200MG	25316-40-9
PHR1790-1G	EICOSAPENTAENOIC ACID ETHYL ESTER	1G	86227-47-6
PHR2120-500MG	EMTRICITABINE	500MG	143491-57-0
PHR1766-1G	EPRINOMECTIN	1G	123997-26-2
PHR1374-1G	ETHYLCELLULOSE	1G	
PHR1828-500MG	EUCALYPTOL	500MG	470-82-6
PHR1866-1G	EZETIMIBE	1G	163222-33-1
PHR1482-1G	FERROUS FUMARATE	1G	141-01-5
PHR1483-1G	FERROUS SULFATE	1G	7782-63-0
PHR1791-100MG	FERULIC ACID	100MG	537-98-4
PHR1792-1G	FLUPHENAZINE HYDROCHLORIDE	1G	146-56-5
PHR1793-500MG	FLUVOXAMINE MALEATE	500MG	61718-82-9

Cat. No.	Description	Pack Size	CAS No.
PHR1833-1G	GALACTITOL	1G	608-66-2
PHR2042-200MG	GALACTOSAMINE HYDROCHLORIDE	200MG	1772-03-8
PHR1623-500MG	GALANTAMINE HYDROBROMIDE	500MG	1953-04-4
PHR1593-1G	GANCICLOVIR	1G	82410-32-0
PHR1859-100MG	GENISTEIN	100MG	446-72-0
PHR1567-500MG	GIBBERELLIC ACID	500MG	77-06-5
PHR1748-3X1.5ML	GLACIAL ACETIC ACID	3X1.5ML	64-19-7
PHR1484-1G	GLYCERYL MONOOLEATE 40%	1G	25496-72-4
PHR1516-200MG	GLYCYRRHIZIC ACID	200MG	1405-86-3
PHR1794-500MG	HESPERIDIN	500MG	520-26-3
PHR1939-500MG	HYDROXYPROLINE	500MG	51-35-4
PHR1440-1G	HYDROXYPROPYL-BETA-CYCLODEXTRIN	1G	128446-35-5
PHR1829-500MG	HYPROMELLOSE	500MG	9004-65-3
PHR1795-100MG	IMINODIBENZYL	100MG	494-19-9
PHR1796-200MG	IMIPENEM	200MG	74431-23-5
PHR1797-1G	IMIPRAMINE HYDROCHLORIDE	1G	113-52-0
PHR1931-1G	IOHEXOL	1G	66108-95-0
PHR1769-1G	ISOMALT	1G	64519-82-0
PHR1932-500MG	ISONIAZID	500MG	54-85-3
PHR1834-200MG	ITRACONAZOLE	200MG	84625-61-6
PHR1798-250MG	LEVAMISOLE HYDROCHLORIDE	250MG	16595-80-5
PHR1504-500MG	LIOTHYRONINE	500MG	6893-02-3
PHR1927-1G	LOPINAVIR	1G	192725-17-0
PHR2186-500MG	LUMEFANTRINE	500MG	82186-77-4
PHR1699-3X100MG	LUTEIN	3X100MG	127-40-2
PHR1770-500MG	LYCOPENE	500MG	502-65-8
PHR1770-5X100MG	LYCOPENE	5X100MG	502-65-8
PHR1363-2.5G	MAGNESIUM STEARATE	2.5G	
PHR2044-1G	MANNOSE	1G	3458-28-4
PHR1955-500MG	MEPIVACAINE HYDROCHLORIDE	500MG	1722-62-9
PHR1943-1G	METHACHOLINE CHLORIDE	1G	62-51-1
PHR1594-1G	METHACYCLINE HYDROCHLORIDE	1G	3963-95-9
PHR2287-100MG	METHYL 4-(ACETYLAMINO)-2-HYDROXYBENZOATE	100MG	4093-28-1
PHR1835-500MG	METHYL NICOTINATE	500MG	93-60-7
PHR1571-1G	METHYL OLEATE	1G	112-62-9
PHR1414-1G	METHYL STEARATE	1G	112-61-8
PHR1801-500MG	MINOCYCLINE HYDROCHLORIDE	500MG	13614-98-7
PHR1490-1ML	MONOETHANOLAMINE	1ML	141-43-5
PHR1802-300MG	NALOXONE	300MG	51481-60-8
PHR1698-1G	NEOTAME	1G	243116-62-3
PHR1731-1G	NITISINONE	1G	104206-65-7
PHR2119-500MG	NORGESTIMATE	500MG	35189-28-7
PHR2210-50MG	NORVALINE	50MG	760-78-1
PHR1880-10MG	OCTREOTIDE ACETATE	10MG	79517-01-4
PHR1781-500MG	OSELTAMIVIR PHOSPHATE	500MG	204255-11-8
PHR1945-300MG	OUABAIN	300MG	11018-89-6
PHR1632-500MG	PIOGLITAZONE HYDROCHLORIDE	500MG	112529-15-4
PHR1588-200MG	PIROXICAM	200MG	36322-90-4
PHR1518-1G	POLYDIMETHYLSILOXANE	1G	9016-00-6
PHR2203-1G	POTASSIUM NITRATE	1G	7757-79-1
PHR2253-100MG	PROCAINE N-OXIDE	100MG	1027072-50-9
PHR1161-1G	PROCAINE HYDROCHLORIDE	1G	51-05-8"
PHR1575-500MG	PROPANTHELINE BROMIDE	500MG	50-34-0
PHR1944-1G	QUINAPRIL HYDROCHLORIDE	1G	82586-55-8

Cat. No.	Description	Pack Size	CAS No.
PHR1750-1G	RESERPINE	1G	50-55-5
PHR2101-500MG	RIBOSE	500MG	50-69-1
PHR1946-300MG	RIFAXIMIN	300MG	80621-81-4
PHR1734-1G	RITONAVIR	1G	155213-67-5
PHR1889-500MG	RIZATRIPTAN BENZOATE	500MG	145202-66-0
PHR1947-500MG	SALMETEROL XINAFOATE	500MG	94749-08-3
PHR1715-1G	SISOMICIN SULFATE	1G	53179-09-2
PHR1471-1G	SODIUM ALGINATE	1G	9005-38-3
PHR1948-2G	SODIUM CARBONATE ANHYDROUS	2G	497-19-8
PHR1645-1G	SODIUM IRON EDTA	1G	15708-41-5
PHR1949-2G	SODIUM LAURYL SULFATE	2G	151-21-3
PHR1444-1G	SODIUM STEARYL FUMARATE	1G	4070-80-8
PHR2245-1G	SODIUM TRIPOLYPHOSPHATE HEXAHYDRATE	1G	15091-98-2
PHR1950-500MG	SPIRAMYCIN	500MG	8025-81-8
PHR1001-1G	SUCROSE	1G	57-50-1
PHR1498-1G	SULFACETAMIDE	1G	144-80-9
PHR1577-1G	SULFAMERAZINE	1G	127-79-7
PHR1810-500MG	TADALAFIL	500MG	171596-29-5
PHR1811-100MG	TEICOPLANIN	100MG	61036-62-2
PHR1855-500MG	TELMISARTAN	500MG	144701-48-4
PHR1592-1G	TENOFOVIR	1G	206184-49-8
PHR1853-500MG	TETRAHYDROZOLINE HYDROCHLORIDE	500MG	522-48-5
PHR2205-1G	TETRASODIUM PYROPHOSPHATE DECAHYDRATE	1G	13472-36-1
PHR1837-100MG	THEOBROMINE	100MG	83-67-0
PHR1890-50MG	TIOTROPIUM BROMIDE MONOHYDRATE	50MG	411207-31-3
PHR1778-1G	TOLNAFTATE	1G	2398-96-1
PHR1891-500MG	TOLTERODINE TARTRATE	500MG	124937-52-6
PHR1600-1G	TOPIRAMATE	1G	97240-79-4
PHR2217-100MG	TOXOPYRIMIDINE	100MG	73-67-6
PHR2201-200MG	TRANS-RESVERATROL	200MG	501-36-0
PHR1749-1G	TRIMEPRAZINE TARTRATE	1G	4330-99-8
PHR1578-5ML	TROLAMINE	5ML	102-71-6
PHR1929-500MG	UBIDECARENONE	500MG	303-98-0
PHR1732-4X250MG	VANCOMYCIN HYDROCHLORIDE	4X250MG	1404-93-9
PHR2039-500MG	XYLOMETAZOLINE HYDROCHLORIDE	500MG	1218-35-5
PHR2102-500MG	XYLOSE	500MG	58-86-6
PHR8933-100MG	(17A)-PREGN-4-ENE-3,20-DIONE	100MG	2000-66-0
PHR8932-100MG	(20R)-20-HYDROXYPREGN-4-EN-3-ONE		145-15-3
PHR8939-100MG	(20R)-3-OXOPREGN-4-EN-20-YL ACETATE		5062-62-4
PHR8935-100MG	(20S)-20-HYDROXYPREGN-4-EN-3-ONE		145-14-2
PHR8943-100MG	(20S)-3-OXOPREGN-4-EN-20-YL ACETATE		5035-09-6
PHR9105-25MG	(2R)-2-[(2,2-DIMETHYLPROPANOYL)AMINO]-2-PHENYLACETIC ACID		40610-41-1
PHR3321-50MG	(2RS)-2-(4-METHYLPHENYL)PROPANOIC ACID		938-94-3
PHR3323-50MG	(2RS)-2-[4-(2-METHYLPROPYL)PHENYL]PROPAN-1-OL		36039-36-8
PHR3273-50MG	(2S)-2-AMINO-3-(PHENYLAMINO)PROPANOIC ACID		145545-23-9
PHR9059-100MG	(2S)-N-(2,6-DIMETHYLPHENYL)PIPERIDINE-2-CARBOXAMIDE		27262-40-4
PHR9018-100MG	(3RS)-9-FLUORO-3-METHYL-10-(4-METHYL-4-OXIDOPIPERAZIN-1-YL)-7-OXO-2,3-DIHYDRO-7H-PYRIDO[1,2,3-DE]-1,4-BENZOXAZINE-6-CARBOXYLIC ACID		104721-52-0
PHR9015-100MG	(3RS)-9-FLUORO-3-METHYL-10-(4-METHYLPIPERAZIN-1-YL)-2,3-DIHYDRO-7H-PYRIDO[1,2,3-DE]-1,4-BENZOXAZIN-7-ONE		123155-82-8
PHR8917-40MG	(5E,7E)-9,10-SECOCHOLESTA-5,7,10(19)-TRIEN-3BETA-OL		67-97-0
PHR9003-100MG	1-(2,4-DICHLOROBENZOYLMETHYL)IMIDAZOLE		46503-52-0
PHR8949-25MG	1-(3-INDOLYLMETHYL)-1,2,3,4-TETRAHYDRO-B-CARBOLINE-3-CARBOXYLIC ACID		
PHR2135-50MG	1-(4-ISOBUTYLPHENYL)ETHANOL		
PHR2421-100MG	1,1-DIMETHYLETHYL (4R-CIS)-6-AMINOETHYL-2,2-DIMETHYL-1,3-DIOXANE-4-ACETATE		

Cat. No.	Description	Pack Size	CAS No.
PHR9178-100MG	1,2,4-TRIAZOLE		288-88-0
PHR3492-1ML	1,2-DICHLOROBENZENE		95-50-1
PHR3494-1ML	1,2-DICHLOROPROPANE		78-87-5
PHR9064-100MG	1,3-DIMETHYLURIC ACID		944-73-0
PHR3167-1ML	1,3-PROPANEDIOL		504-63-2
PHR1028-1G	1,4-BENZOQUINONE		106-51-4
PHR2558-500MG	1,4-SORBITAN		27299-12-3
PHR9001-100MG	1-[(2RS)-2-(2,4-DICHLOROPHENYL)-2-[(2,5-DICHLOROPHENYL)METHYL]OXY]ETHYL]-1H-IMIDAZOLIUM NITRATE		909277-72-1
PHR9002-100MG	1-[(2RS)-2-(BENZYLOXY)-2-(2,4-DICHLOROPHENYL)ETHYL]-1H-IMIDAZOLIUM NITRATE		203264-09-9
PHR2791-100MG	1-[(2RS)-2-[[4-(4-CHLOROPHENYL)METHYL]OXY]-2-(2,4-DICHLOROPHENYL)ETHYL]-1H-IMIDAZOLE		27220-47-9
PHR9066-100MG	10-BROMO-5H-DIBENZ[B,F]AZEPINE-5-CARBOXAMIDE		59690-97-0
PHR3007-50MG	10-DESACETYL-BACCATIN III		32981-86-5
PHR9115-50MG	17,21-DIHYDROXYPREGNA-1,4,9(11)-TRIENE-3,20-DIONE		10184-69-7
PHR9004-100MG	1-CYCLOPROPYL-6-FLUORO-7-(PIPERAZIN-1-YL)QUINOLIN-4(1H)-ONE		105394-83-0
PHR9014-100MG	1-CYCLOPROPYL-7-(4-ETHYLPIPERAZIN-1-YL)-6-FLUOROQUINOLIN-4(1H)-ONE		131775-99-0
PHR2130-50MG	1-ETHYL-4-ISOBUTYLBENZENE		
PHR9010-100MG	1-ETHYL-6-FLUORO-7-(4-FORMYLPIPERAZIN-1-YL)-4-OXO-1,4-DIHYDROQUINOLINE-3-CARBOXYLIC ACID		70459-04-0
PHR9008-100MG	1-ETHYL-6-FLUORO-7-(PIPERAZIN-1-YL)QUINOLIN-4(1H)-ONE		75001-82-0
PHR9166-100MG	1-FLUORONAPHTHALENE		321-38-0
PHR2497-50MG	1-METHYL-1,2,3,4-TETRAHYDRO-B-CARBOLINE-3-CARBOXYLIC ACID		5470-37-1
PHR9089-25MG	1-METHYLETHYL 4-HYDROXY-3,5-BIS(1-METHYLETHYL)BENZOATE		2095678-97-8
PHR1740-200MG	1-VINYL-2-PYRROLIDINONE		88-12-0
PHR9031-100MG	2-(2,4-DIFLUOROPHENYL)-1,2,3-PROPANETRIOL		173837-65-5
PHR8903-50MG	2-(3-INDOLYLMETHYL)-L-TRYPTOPHAN		SIGN IN TO VIEW
PHR3715-50MG	2-(4-ISOBUTYRYLPHENYL)PROPANAMIDE		1823308-53-7
PHR2125-50MG	2-(4-ISOBUTYRYLPHENYL)PROPANE		
PHR2127-50MG	2-(4-N-PROPYLPHENYL)PROPANOIC ACID		
PHR9134-100MG	2-(4-TERT-BUTYL-3-HYDROXY-2,6-DIMETHYLPHENYL)ACETIC ACID		55699-12-2
PHR9161-50MG	2-(4-TERT-BUTYL-3-HYDROXY-2,6-DIMETHYLPHENYL)ACETONITRILE		55699-10-0
PHR9072-100MG	2-(5-NITRO-1H-IMIDAZOL-1-YL)ETHANOL		5006-68-8
PHR8743-50MG	2,2'-BISNALOXONE		51481-60-8
PHR9074-25MG	2,2-DIMETHYL-4-(1-METHYLETHYL)-1,3-BENZODIOXOLE		201166-22-5
PHR9092-100MG	2,2-DIPROPYLPENTANENITRILE		5340-48-7
PHR9093-100MG	2,2-DIPROPYLPENTANOIC ACID		52061-75-3
PHR9073-100MG	2,4-DICHLORO-5-SULFAMOYLBENZOIC ACID		2736-23-4
PHR2595-200MG	2,4-DICHLOROPHENOL		120-83-2
PHR9032-50MG	2,4-DIISOPROPYLPHENOL		2934-05-6
PHR9104-75MG	2,5-DEOXYFRUCTOSAZINE		17460-13-8
PHR3329-50MG	2,5-DIHYDROXYBENZOIC ACID		490-79-9
PHR9075-25MG	2,5-DIISOPROPYLPHENOL		35946-91-9
PHR9155-100MG	2,6-DICHLORO-N-(2,6-DIMETHYLPHENYL)HEXANAMIDE		1037184-07-8
PHR9096-100MG	2-[[4-(4-CHLORO-3,5-DIMETHYL-2-PYRIDINYL)METHYL]THIO]-6-METHOXY-1H-BENZIMIDAZOLE		220757-74-4
PHR9000-100MG	2-[1-[(2RS)-2-(2,4-DICHLOROPHENYL)-2-[(2,4-DICHLOROPHENYL)METHYL]OXY]ETHYL]-1H-IMIDAZOL-3-IUM-3-YL]-2-METHYLPROPANOATE		
PHR9044-50MG	2-[2-[4-(DIPHENYLMETHYL)-1-PIPERAZINYL]ETHOXY]ETHANOL		3733-63-9
PHR2132-50MG	2-[4-(2-HYDROXY-2-METHYLPROPYL)-PHENYL]PROPANOIC ACID		51146-55-5
PHR3322-50MG	2-[4-(2-METHYLPROPYL)PHENYL]ETHANOL		36039-35-7
PHR9169-100MG	2-[4-(BUTYLAMINO)BENZOYL]OXY-N,N-DIMETHYLETHANAMINE-N-OXIDE HYDROCHLORIDE		55750-02-2
PHR9148-50MG	2-[5-(1H-1,2,4-TRIAZOL-1-YLMETHYL)-1H-INDOL-3-YL]ETHANOL		160194-39-8
PHR3804-100MG	2'-O-METHYLCYTIDINE		2140-72-9
PHR8940-100MG	21-(CYCLOHEXYLIDENE)PREGN-4-ENE-3,20-DIONE		
PHR1507-200MG	2-AMINO-5-CHLOROBENZOPHENONE		719-59-5

Cat. No.	Description	Pack Size	CAS No.
PHR9179-100MG	2-CHLORO-2',4'-DIFLUOROACETOPHENONE	—	51336-94-8
PHR9053-100MG	2-CHLORO-3-NITROBENZOIC ACID		3970-35-2
PHR1741-200MG	2-CHLORO-5-NITROBENZOIC ACID		2516-96-3
PHR3330-50MG	2-CHLOROBENZOIC ACID		118-91-2
PHR3491-1ML	2-CHLOROETHYL VINYL ETHER		110-75-8
PHR9054-100MG	2-CHLOROMETHYL-4-METHOXY-3,5-DIMETHYLPYRIDINE HYDROCHLORIDE		86604-75-3
PHR9231-100MG	2-DEOXY- -GLUCOSE		154-17-6
PHR9061-100MG	2-ETHYLPENTANOIC ACID		20225-24-5
PHR3290-50MG	2-FURALDEHYDE		98-01-1
PHR3502-1ML	2-HEXANONE		591-78-6
PHR9219-100MG	2-HYDROXY-2,2-DI(2-THIENYL)ACETIC ACID		
PHR9056-100MG	2-HYDROXY-5-NITROBENZOIC ACID		96-97-9
PHR1742-200MG	2-HYDROXY-NITROBENZOIC ACID		619-19-2
PHR2126-50MG	2-HYDROXYPROPYL 2-(4-ISOBUTYLPHENYL)-PROPANOATE		95093-58-6
PHR1223-500MG	2-METHYLPYRAZINE		109-08-0
PHR9057-100MG	2-METHYLVALERIC ACID		97-61-0
PHR1121-1.5G	2-PHENOXYETHANOL	1.5G	122-99-6
PHR1121-10G	2-PHENOXYETHANOL	10G	122-99-6
PHR3542-1G	2-PHENYLBUTYRIC ACID	1G	90-27-7
PHR2335-5X1.2ML	2-PROPANOL IMPURITY MIXTURE	1.2ML	
PHR9077-100MG	2-PROPYLPENTANENITRILE	100MG	13310-75-3
PHR2336-500MG	2-TERT-BUTYL-4-HYDROXYANISOLE		88-32-4
PHR2619-100MG	3,4-DIMETHOXYPHENYLACETONITRILE		93-17-4
PHR1339-500MG	3-AMINO-1-PROPANOL	500MG	62-23-7
PHR1743-200MG	3-AMINOBENZOIC ACID	200MG	99-05-8
PHR1744-200MG	3-CHLOROBENZOIC ACID	200MG	535-80-8
PHR9165-50MG	3-ETHYL 5-METHYL 2-(2-AMINOETHOXYMETHYL)-4-(2-CHLOROPHENYL)-6-METHYLPYRIDINE-3,5-DICARBOXYLATE OXALATE SALT	50MG	1216406-90-4
PHR1354-500MG	3-iodo-2-propynyl N-butylcarbamate	500MG	
PHR8950-50MG	3-O-ALLYLNALOXONE	50MG	1352085-46-1
PHR8744-50MG	3-O-METHYLNALOXONE	50MG	70866-64-7
PHR9078-25MG	3-PHENYLPYRAZIN-2-OL	25MG	73200-73-4
PHR1306-500MG	3-TERT-BUTYL-4-HYDROXYANISOLE	500MG	121-00-6
PHR9195-100MG	4-(4-CHLOROPHENYL)-4-HYDROXYPIPERIDINE	100MG	39512-49-7
PHR9095-100MG	4-(DIETHYLAMINO)BUT-2-YNYL 2-HYDROXY-2,2-DIPHENYLACETATE	100MG	14943-53-4
PHR3640-50MG	4-[(1RS)-2-(TERT-BUTYLAMINO)-1-HYDROXYETHYL]-2-ETHYLPHENOL HYDROCHLORIDE	50MG	
PHR9050-50MG	4-[2-(TERT-BUTYLAMINO)-1-METHOXYETHYL]-2-(HYDROXYMETHYL)PHENOL	50MG	870076-72-5
PHR2106-50MG	4-[2-[(1,1-DIMETHYLETHYL)AMINO]ETHYL]-2-METHYLPHENOL	50MG	132183-64-3
PHR9071-50MG	4-[2-[BENZYL(TERT-BUTYL)AMINO]-1-HYDROXYETHYL]-2-(HYDROXYMETHYL)PHENOL	50MG	24085-03-8
PHR9194-50MG	4-[3-(METHYLAMINO)-1-(THIOPHEN-2-YL)PROPYL]NAPHTHALEN-1-OL HYDROCHLORIDE	50MG	953028-76-7
PHR9208-100MG	4-[4-(4-CHLOROPHENYL)-3,6-DIHYDROPYRIDIN-1(2H)-YL]-N,N-DIMETHYL-2,2-DIPHENYLBUTANAMIDE	100MG	61299-42-1
PHR9097-50MG	4-[4-[4-(DIPHENYLMETHYLENE)-1-PIPERIDINYL]-1-HYDROXYBUTYL]-A, A-DIMETHYLBENZENEACETIC ACID	50MG	1187954-57-9
PHR1747-200MG	4'-AMINOACETOPHENONE	200MG	99-92-3
PHR2483-100MG	4-ACETYLBENZOIC ACID	100MG	586-89-0
PHR8944-50MG	4-AMINO-1-[(2R,3S,5S)-2-(HYDROXYMETHYL)-1,3-OXATHIOLAN-5-YL]PYRIMIDIN-2(1H)-ONE S-OXIDE	50MG	
PHR9180-100MG	4-AMINO-4H-1,2,4-TRIAZOLE	100MG	584-13-4
PHR9189-50MG	4-BENZYLPIRIDINE	50MG	2116-65-6
PHR2620-50MG	4-CHLOROBENZOPHENONE	50MG	134-85-0
PHR1691-100MG	4-EPIANHYDROTETRACYCLINE HYDROCHLORIDE	100MG	4465-65-0
PHR8752-30MG	4-EPIOXYTETRACYCLINE HCL	30MG	

Cat. No.	Description	Pack Size	CAS No.
PHR9049-25MG	4-HYDROXY-3,5-BIS(1-METHYLETHYL)BENZOIC ACID	25MG	13423-73-9
PHR1048-1G	4-HYDROXYBENZOIC ACID	1G	99-96-7
PHR2269-100MG	4-ISOBUTYLBENZOIC ACID	100MG	38861-88-0
PHR2134-50MG	4-ISOBUTYLSTYRENE	50MG	
PHR2140-50MG	4-ISOPROPYLCINNAMIC ACID	50MG	3368-21-6
PHR9185-50MG	4-ISOPROPYLPHENOL	50MG	99-89-8
PHR3274-3X1.2ML	4-METHYL-2-PENTANOL	3X1.2ML	108-11-2
PHR2339-100MG	4-NITROBENZOIC ACID	100MG	62-23-7
PHR9223-100MG	4-NITROSO-HYDROCHLOROTHIAZIDE	100MG	63779-86-2
PHR3365-50MG	5-(2-HYDROXYETHYL)-4-METHYLTHIAZOLE	50MG	137-00-8
PHR2247-50MG	5-(4-METHYLPHENYL)HYDANTOIN	50MG	
PHR9079-100MG	5-(AMINOSULFONYL)-4-CHLORO-2-[[[(TETRAHYDRO-2-FURANYL)METHYL]AMINO]BENZOIC ACID	100MG	4793-38-8
PHR2584-100MG	5,6-DIAMINO-1,3-DIMETHYLURACIL	100MG	
PHR3831-100MG	5-[4-(AZIDOMETHYL)[1,1-BIPHENYL]-2-YL]-2H-TETRAZOLE	100MG	152708-24-2
PHR2340-1G	5'-ADENYLIC ACID	1G	61-19-8
PHR2597-100MG	5-HYDROXY- -TRYPTOPHAN	100MG	4350-09-8
PHR9080-100MG	5-METHOXY-1,3-DIHYDROBENZIMIDAZOL-2-ONE	100MG	2080-75-3
PHR8877-100MG	5-METHOXY-2-(4-METHOXY-3,5-DIMETHYLPYRIDIN-2-YL)-1H-BENZIMIDAZOLE	100MG	
PHR3806-100MG	5-METHYLCYTIDINE	100MG	SIGN IN TO VIEW
PHR3800-100MG	5-METHYLURIDINE	100MG	1463-10-1
PHR9196-100MG	5A-ANDROSTAN-17B-OL-3-ONE	100MG	521-18-6
PHR2470-100MG	6-HYDROXYNICOTINIC ACID	100MG	5006-66-6
PHR2218-100MG	6-TERT-BUTYL-3-(4,5-DIHYDRO-1H-IMIDAZOL-2-YLMETHYL)-4-HYDROXY-2,4-DIMETHYL-CYCLOHEX-2,5-DIENONE TRIFLUOROACETATE	100MG	
PHR9144-50MG	7-[(2-AMINOETHYL)AMINO]-1-CYCLOPROPYL-6-FLUORO-1,4-DIHYDRO-4-OXO-3-QUINOLINECARBOXYLIC ACID	50MG	103222-12-4
PHR3142-50MG	7-METHYLYXANTHINE	50MG	552-62-5
PHR1253-500MG	ACARBOSE	500MG	56180-94-0
PHR3166-100MG	ACECLOFENAC	100MG	89796-99-6
PHR1266-500MG	ACESULFAME K	500MG	55589-62-3
PHR3523-100MG	ACETALDEHYDE-2,4-DNPH	100MG	1019-57-4
PHR1086-1G	ACETANILIDE MELTING POINT STANDARD	1G	103-84-4
PHR3614-50MG	ACETYL SIMVASTATIN	50MG	145576-25-6
PHR1383-1G	ADENINE	1G	73-24-5
PHR2341-200MG	ADIPIC ACID	200MG	124-04-9
PHR2478-50MG	ADRENALONE HYDROCHLORIDE	50MG	62-13-5
PHR2342-100MG	ALEURITIC ACID	100MG	533-87-9
PHR1638-1G	ALFUZOSIN HYDROCHLORIDE	1G	81403-68-1
PHR1364-500MG	ALLANTOIN	500MG	97-59-6
PHR2969-1G	ALMOND OIL	1G	8007-69-0
PHR2560-1G	ALPHA CYCLODEXTRIN	1G	10016-20-3
PHR2561-1G	ALPHA LIPOIC ACID	1G	1077-28-7
PHR8912-100MG	ALPHA-APO-OXYTETRACYCLINE	100MG	18695-01-7
PHR2343-1G	ALTRETAMINE	1G	645-05-6
PHR3169-500MG	AMCINONIDE	500MG	51022-69-6
PHR3170-1G	AMILOXATE	1G	71617-10-2
PHR3171-500MG	AMINO BENZOATE POTASSIUM	500MG	138-84-1
PHR3172-500MG	AMINO BENZOATE SODIUM	500MG	555-06-6
PHR1224-1G	AMINOCAPROIC ACID	1G	60-32-2
PHR3173-500MG	AMINOGLUTETHIMIDE	500MG	125-84-8
PHR3174-500MG	AMINOHIPURIC ACID	500MG	61-78-9
PHR1327-1G	AMMONIUM CHLORIDE	1G	12125-02-9
PHR1689-200MG	AMMONIUM GLYCYRRHIZATE	200MG	53956-04-0

Cat. No.	Description	Pack Size	CAS No.
PHR2344-1G	AMODIAQUINE HYDROCHLORIDE HYDRATE	1G	6398-98-7
PHR1662-500MG	AMPHOTERICIN B	500MG	1397-89-3
PHR2345-500MG	AMPROLIUM	500MG	
PHR3175-50MG	ANHYDROTETRACYCLINE HYDROCHLORIDE	50MG	13803-65-1
PHR3176-500MG	ANHYDROUS SODIUM SUCCINATE	500MG	150-90-3
PHR2346-500MG	ANILINE	500MG	62-53-3
PHR1677-1G	ANTHRALIN (DITHRANOL)	1G	1143-38-0
PHR3422-300MG	APREPITANT	300MG	170729-80-3
PHR3530-5G	ARACHIDIC ACID	5G	506-30-9
PHR3402-500MG	ARGATROBAN MONOHYDRATE	500MG	141396-28-3
PHR2813-50MG	ARGININE ASPARTATE	50MG	7675-83-4
PHR2573-300MG	ARTESUNATE	300MG	88495-63-0
PHR2482-200MG	ATRACURIUM BESYLATE	200MG	64228-81-5
PHR3846-100MG	ATROPINE	100MG	51-55-8
PHR1379-1G	ATROPINE SULFATE	1G	5908-99-6
PHR1073-1G	AVOBENZONE	1G	70356-09-1
PHR2484-100MG	AZURE B	100MG	531-55-5
PHR9170-50MG	BACCATIN III	50MG	27548-93-2
PHR1590-1G	BACITRACIN	1G	1405-87-4
PHR1559-1G	BACITRACIN ZINC	1G	1405-89-6
PHR1619-1G	BECLOMETHASONE DIPROPIONATE	1G	5534-09-8
PHR2793-400MG	BEMOTRIZINOL	400MG	187393-00-6
PHR2352-300MG	BENOXINATE HYDROCHLORIDE	300MG	5987-82-6
PHR1203-3X1ML	BENZALDEHYDE	3X1ML	100-52-7
PHR2353-1ML	BENZALDEHYDE PROPYLENE GLYCOL ACETAL	1ML	2568-25-4
PHR2354-300MG	BENZALKONIUM BROMIDE	300MG	
PHR3795-5.2ML	BENZALKONIUM CHLORIDE (4% AQUEOUS SOLUTION)	5.2ML	
PHR1371-5ML	BENZALKONIUM CHLORIDE 10% SOLUTION	5ML	63449-41-2
PHR1681-1ML	BENZALKONIUM CHLORIDE 50% SOLUTION	1ML	63449-41-2
PHR1310-1G	BENZENE	1G	71-43-2
PHR1425-1G	BENZETHONIUM CHLORIDE	1G	121-54-0
PHR9055-100MG	BENZIL	100MG	134-81-6
PHR1158-1G	BENZOCAINE	1G	94-09-7
PHR1050-1G	BENZOIC ACID	1G	65-85-0
PHR1204-1G	BENZOPHENONE	1G	119-61-9
PHR1019-1G	BENZYL ALCOHOL	1G	100-51-6
PHR1213-5G	BENZYL BENZOATE	5G	120-51-4
PHR1370-1ML	BENZYL NICOTINATE	1ML	94-44-0
PHR1502-500MG	BERBERINE CHLORIDE	500MG	633-65-8
PHR1349-1G	BETA ALANINE	1G	107-95-9
PHR8913-100MG	BETA-APO-OXYTETRACYCLINE	100MG	18751-99-0
PHR1439-1G	BETA-CYCLODEXTRIN	1G	7585-39-9
PHR2923-1G	BETADEX SULFOBUTYL ETHER SODIUM	1G	182410-00-0
PHR2355-1G	BETAINE	1G	107-43-7
PHR2356-300MG	BETAINE HYDROCHLORIDE	300MG	590-46-5
PHR2357-300MG	BETHANECHOL CHLORIDE	300MG	590-63-6
PHR1296-1G	BIFONAZOLE	1G	60628-96-8
PHR1233-1G	BIOTIN	1G	58-85-5
PHR2358-1G	BIPHENYL	1G	92-52-4
PHR3434-50MG	BISDESMETHOXYCURCUMIN	50MG	33171-05-0
PHR9197-100MG	BOLDENONE	100MG	846-48-0
PHR2970-1G	BORAGE SEED OIL	1G	84012-16-8
PHR2708-500MG	BOSENTAN	500MG	157212-55-0
PHR3424-400MG	BRINZOLAMIDE	400MG	138890-62-7

Cat. No.	Description	Pack Size	CAS No.
PHR3495-1ML	BROMOFORM	1ML	75-25-2
PHR1304-1G	BUTYLATED HYDROXYANISOLE	1G	25013-16-5
PHR1117-1G	BUTYLATED HYDROXYTOLUENE	1G	128-37-0
PHR1022-10G	BUTYLPARABEN	10G	94-26-8
PHR2730-1.2ML	BUTYRALDEHYDE	1.2ML	123-72-8
PHR2730-5X1.2ML	BUTYRALDEHYDE	5X1.2ML	123-72-8
PHR2359-100MG	CALCIFEDIOL	100MG	63283-36-3
PHR1362-5G	CALCIUM ACETATE	5G	62-54-4
PHR1326-5G	CALCIUM CARBONATE	5G	471-34-1
PHR2669-500MG	CALCIUM GLUCONATE ANHYDROUS	500MG	299-28-5
PHR2670-500MG	CALCIUM GLUCONATE MONHYDRATE	500MG	66905-23-5
PHR1232-500MG	CALCIUM-D-PANTOTHENATE	500MG	137-08-6
PHR2971-1G	CANOLA OIL	1G	120962-03-0
PHR1202-1G	CAPRYLIC ACID (OCTANOIC ACID)	1G	124-07-2
PHR1450-1G	CAPSAICIN	1G	404-86-4
PHR1307-1G	CAPTOPRIL	1G	62571-86-2
PHR1511-500MG	CARBACHOL	500MG	51-83-2
PHR1655-1G	CARBIDOPA	1G	38821-49-7
PHR3417-200MG	CARBOPLATIN	200MG	41575-94-4
PHR2726-2G	CARBOXYMETHYLCELLULOSE SODIUM	2G	9004-32-4
PHR3819-400MG	CARGLUMIC ACID	400MG	1188-38-1
PHR2262-100MG	CARMOISINE	100MG	3567-69-9
PHR2972-3X200MG	CASTOR OIL	3X200MG	8001-79-4
PHR1283-1G	CEFACLOR	1G	70356-03-5
PHR1521-1G	CEFTIOFUR SODIUM	1G	
PHR1382-1G	CEFTRIAXONE SODIUM	1G	104376-79-6
PHR2688-1G	CETRIMONIUM BROMIDE	1G	
PHR1133-1G	CETYL ALCOHOL	1G	36653-82-4
PHR1166-1G	CETYL PALMITATE	1G	540-10-3
PHR1226-1G	CETILPYRIDINIUM CHLORIDE	1G	6004-24-6
PHR1412-1G	CHLORAMPHENICOL	1G	
PHR1545-1G	CHLOROBUTANOL	1G	6001-64-5
PHR1179-500MG	CHLOROTHIAZIDE	500MG	58-94-6
PHR1478-500MG	CHLOROXYLENOL	500MG	88-04-0
PHR1284-1G	CHLORPROPAMIDE	1G	94-20-2
PHR1520-500MG	CHLORTETRACYCLINE HYDROCHLORIDE	500MG	64-72-2
PHR1533-500MG	CHOLESTEROL	500MG	57-88-5
PHR1460-500MG	CHOLINE BITARTRATE	500MG	87-67-2
PHR1251-500MG	CHOLINE CHLORIDE	500MG	67-48-1
PHR3284-3X200MG	CINNAMALDEHYDE	3X200MG	14371-10-9
PHR2694-200MG	CINNAMIC ACID	200MG	140-10-3
PHR2811-200MG	CINNARIZINE	200MG	298-57-7
PHR3506-1ML	CIS-DICHLOROETHYLENE	1ML	156-59-2
PHR1624-200MG	CISPLATIN	200MG	15663-27-1
PHR2903-200MG	CITRIC ACID MONOHYDRATE	200MG	5949-29-1
PHR1071-1G	CITRIC ACID, ANHYDROUS	1G	77-92-9
PHR1597-1ML	CMIT/MIT	1ML	
PHR3520-1G	COCONUT OIL, 1000MG, NEAT	1G	
PHR2974-1G	COD LIVER OIL	1G	8001-69-2
PHR1477-1G	COPPER SULFATE	1G	7758-98-7
PHR2897-5X400MG	CORN OIL	5X400MG	8001-30-7
PHR2976-1G	COTTONSEED OIL	1G	8001-29-4
PHR1462-1G	CREATININE	1G	60-27-5
PHR2695-50MG	CYANIDIN CHLORIDE	50MG	528-58-5"

Cat. No.	Description	Pack Size	CAS No.
PHR1234-1G	CYANOCOBALAMIN	1G	68-19-9
PHR2898-3X1.2ML	CYCLOHEXYLMETHANOL	3X1.2ML	100-49-2
PHR1092-500MG	CYCLOSPORINE	500MG	59865-13-3
PHR2809-50MG	CYPROTERONE ACETATE	50MG	427-51-0
PHR2899-500MG	CYROMAZINE	500MG	66215-27-8
PHR3802-100MG	CYTIDINE	100MG	65-46-3
PHR1350-500MG	CYTOSINE	500MG	71-30-7
PHR9060-25MG	-(-)-A-PHENYLGLYCINE	25MG	875-74-1
PHR2623-50MG	DAIDZIN	50MG	552-66-9
PHR1464-1G	DAPSONE	1G	80-08-0
PHR3433-400MG	DAUNORUBICIN HYDROCHLORIDE	400MG	23541-50-6
PHR3411-500MG	DEFEROXAMINE MESYLATE	500MG	138-14-7
PHR1070-5X1.2ML	ETHANOL (DEHYDRATED ALCOHOL)	5X1.2ML	64-17-5
PHR1582-1G	DEHYDROACETIC ACID	1G	520-45-6
PHR3798-100MG	DEOXYADENOSINE MONOHYDRATE	100MG	958-09-8
PHR3803-50MG	DEOXYCYTIDINE	50MG	SIGN IN TO VIEW
PHR3799-50MG	DEOXYGUANOSINE MONOHYDRATE	50MG	961-07-9
PHR3801-100MG	DEOXYURIDINE	100MG	
PHR1300-500MG	DEQUALINIUM CHLORIDE	500MG	522-51-0
PHR2701-500MG	DEXMETETOMIDINE HYDROCHLORIDE	500MG	145108-58-3
PHR1228-500MG	DEXPANTHENOL	500MG	81-13-0
PHR9198-100MG	DI-2-THIENYL KETONE	100MG	704-38-1
PHR9201-100MG	DI-2-THIENYLGLYCOLIC ACID METHYL ESTER	100MG	
PHR2900-5G	DIBASIC POTASSIUM PHOSPHATE	5G	7758-11-4
PHR3830-500MG	DIBUCAINE HYDROCHLORIDE	500MG	61-12-1
PHR1366-1G	DIBUTYL PHTHALATE	1G	84-74-2
PHR1216-3X1.2ML	DIBUTYL SEBACATE	3X1.2ML	109-43-3
PHR1476-1G	DIETHYL PHTHALATE	1G	84-66-2
PHR1045-1G	DIETHYLENE GLYCOL	1G	111-46-6
PHR2733-200MG	DIGLYCINE (GLYCYLGLYCINE)	200MG	556-50-3
PHR2627-50MG	DIHYDROCPSAICIN	50MG	19408-84-5
PHR3589-500MG	DIHYDROERGOTAMINE MESYLATE	500MG	6190-39-2
PHR1430-1G	DIHYDROXYACETONE	1G	96-26-4
PHR8891-100MG	DIISOPROPYLHOMOVERATRYLCYANIDE	100MG	
PHR1309-3G	DIMETHYL SULFOXIDE	3G	67-68-5
PHR2799-500MG	DIOSMIN	500MG	520-27-4
PHR8919-100MG	DIOXINDOLYL-L-ALANINE	100MG	184955-21-3
PHR9184-50MG	DIPHENYL ETHER	50MG	101-84-8
PHR3772-400MG	DIPYRIDAMOLE	400MG	58-32-2
PHR1514-1G	DISODIUM GUANYLATE	1G	5550-12-9
PHR1475-1G	DISODIUM INOSINATE	1G	352195-40-5
PHR9051-100MG	-KYNURENINE	100MG	343-65-7
PHR1358-1ML	DMDM HYDANTOIN	1ML	
PHR3842-400MG	DOBUTAMINE HYDROCHLORIDE	400MG	49745-95-1
PHR2264-1G	DOCOSANOL (N-DOCOSANOL)	1G	661-19-8
PHR2728-1G	DOCUSATE SODIUM	1G	577-11-7
PHR3306-500MG	DOFETILIDE	500MG	115256-11-6
PHR1145-1G	DOXYCYCLINE HYCLATE	1G	24390-14-5
PHR1420-1G	DOXYLAMINE SUCCINATE	1G	562-10-7
PHR2585-500MG	DROPERIDOL	500MG	548-73-2
PHR3314-500MG	DROSPIRENONE	500MG	67392-87-4
PHR2586-1G	EDETIC ACID	1G	60-00-4
PHR1289-1G	ENALAPRIL MALEATE	1G	76095-16-4
PHR3440-500MG	ENALAPRILAT	500MG	84680-54-6

Cat. No.	Description	Pack Size	CAS No.
PHR1513-500MG	ENROFLOXACIN	500MG	93106-60-6
PHR1301-1G	ENSULIZOLE(2-PHENYLBENZIMIDAZOLE 5-SULFONIC ACID)	1G	27503-81-7
PHR3445-500MG	ENTECAVIR MONOHYDRATE	500MG	209216-23-9
PHR1333-1G	EPIGALLOCATECHIN GALLATE	1G	989-51-5
PHR1509-500MG	EPINEPHRINE BITARTRATE	500MG	51-42-3
PHR3815-250MG	EPLERENONE	250MG	107724-20-9
PHR1238-500MG	ERGOCALCIFEROL (VITAMIN D2)	500MG	50-14-6
PHR1512-200MG	ERGOSTEROL	200MG	57-87-4
PHR1479-1G	ERYTHRITOL	1G	149-32-6
PHR2901-500MG	ESCIN	500MG	11072-93-8
PHR1535-500MG	ESTRONE	500MG	53-16-7
PHR3179-1G	ESTROPIPATE	1G	7280-37-7
PHR3513-500MG	ETHOXYQUIN	500MG	91-53-2
PHR9200-100MG	ETHYL 2,2-DIPHENYLACETATE	100MG	3468-99-3
PHR9174-100MG	ETHYL 2-AMINO BENZOATE	100MG	87-25-2
PHR9173-100MG	ETHYL 3-AMINO BENZOATE	100MG	582-33-2
PHR9212-100MG	ETHYL 4-(4,8-DICHLORO-5,6-DIHYDRO-11H-BENZO[5,6]CYCLOHEPTA[1,2-B]PYRIDIN-11-YLIDENE)-1-PIPERIDINECARBOXYLATE	100MG	165739-83-3
PHR2225-100MG	ETHYL 4-NITROBENZOATE	100MG	99-77-4
PHR9082-100MG	ETHYL VALPROATE	100MG	17022-31-0
PHR1268-1G	ETHYL VANILLIN	1G	121-32-4
PHR1046-1G	ETHYLENE GLYCOL	1G	107-21-1
PHR2093-5ML	ETHYLENE OXIDE AND ACETALDEHYDE STANDARD SOLUTION	5ML	
PHR2092-5ML	ETHYLENE OXIDE AND DIOXANE STANDARD SOLUTION	5ML	
PHR1295-1.2ML	ETHYLENE OXIDE STOCK SOLUTION	1.2ML	
PHR2801-400MG	ETHYLHEXYL TRIAZONE	400MG	88122-99-0
PHR1011-10G	ETHYL PARABEN	10G	120-47-8
PHR2668-20MG	ETHYL PARABEN SODIUM	20MG	35285-68-8
PHR2453-200MG	ETONOGESTREL	200MG	54048-10-1
PHR2977-1G	EUCALYPTUS OIL	1G	8000-48-4
PHR1269-1G	EUGENOL	1G	97-53-0
PHR2978-1G	EVENING PRIMROSE OIL	1G	308064-97-3
PHR2263-100MG	FD&C YELLOW NO. 6	100MG	2783-94-0
PHR2790-500MG	FINGOLIMOD HYDROCHLORIDE	500MG	162359-56-0
PHR2980-1G	FLAX SEED OIL	1G	8001-26-1
PHR2921-500MG	FLECAINIDE ACETATE	500MG	54143-56-5
PHR2589-500MG	FLOXURIDINE	500MG	50-91-9
PHR1659-1G	FLUCYTOSINE	1G	2022-85-7
PHR1397-500MG	FLUCINOLONE ACETONIDE	500MG	67-73-2
PHR3405-400MG	FLUROMETHOLONE	400MG	426-13-1
PHR3490-400MG	FLUROMETHOLONE ACETATE	400MG	3801-06-7
PHR1620-1G	FLUVASTATIN SODIUM	1G	93957-55-2
PHR2703-500MG	FORMOTEROL FUMARATE	500MG	183814-30-4
PHR1002-1G	FRUCTOSE	1G	57-48-7
PHR1270-1G	FUMARIC ACID	1G	110-17-8
PHR2598-500MG	FURAZOLIDONE	500MG	67-45-8
PHR2810-50MG	FUSIDIC ACID	50MG	6990-06-3
PHR1206-500MG	GALACTOSE	500MG	59-23-4
PHR2582-50MG	GEMCITABINE HYDROCHLORIDE	50MG	122111-03-9
PHR1286-1G	GEMFIBROZIL	1G	25812-30-0
PHR1077-1G	GENTAMICIN SULFATE	1G	1405-41-0
PHR2200-50MG	GINSENOSIDE RG1	50MG	22427-39-0
PHR1288-1G	GLICLAZIDE	1G	21187-98-4
PHR1020-5G	GLYCERIN	5G	56-81-5

Cat. No.	Description	Pack Size	CAS No.
PHR3180-1G	GLYCERYL DISTEARATE	1G	1323-83-7
PHR1111-1G	GLYCINE	1G	56-40-6
PHR1427-1G	GLYCOLIC ACID	1G	79-14-1
PHR3435-400MG	GLCYL- -TYROSINE DIHYDRATE	400MG	39630-46-1
PHR3009-100MG	GONADORELIN ACETATE	100MG	52699-48-6
PHR3776-100MG	GOSERELIN ACETATE	100MG	145781-92-6
PHR1534-1G	GRISEOFULVIN	1G	126-07-8
PHR1243-500MG	GUANINE	500MG	73-40-5
PHR1428-1G	HEPES	1G	7365-45-9
PHR3181-250MG	HEXYLENE GLYCOL	250MG	107-41-5
PHR1357-500MG	HISTAMINE DIHYDROCHLORIDE	500MG	
PHR3182-500MG	HOMATROPINE HYDROBROMIDE	500MG	51-56-9
PHR3183-500MG	HOMATROPINE METHYLBROMIDE	500MG	80-49-9
PHR1085-1G	HOMOSALATE	1G	118-56-9
PHR1085-10G	HOMOSALATE	10G	118-56-9
PHR8925-100MG	HUMAN INSULIN	100MG	11061-68-0
PHR3184-500MG	HYDRALAZINE HYDROCHLORIDE	500MG	304-20-1
PHR1032-1G	HYDROCHLOROTHIAZIDE	1G	58-93-5
PHR2981-500MG	HYDROGENATED CASTOR OIL	500MG	
PHR1469-1G	HYDROQUINONE	1G	123-31-9
PHR3186-250MG	HYDROXOCOBALAMIN ACETATE	250MG	22465-48-1
PHR1485-1G	HYDROXYPROPYL CELLULOSE	1G	9004-64-2
PHR2629-300MG	HYDROXYUREA	300MG	127-07-1
PHR2691-100MG	HYPEROSIDE	100MG	482-36-0
PHR2704-75MG	IDARUBICIN HYDROCHLORIDE	75MG	57852-57-0
PHR2884-500MG	IDOXURIDINE	500MG	54-42-2
PHR3187-1G	IFOSFAMIDE	1G	3778-73-2
PHR1180-500MG	IMIDAZOLE	500MG	288-32-4
PHR1500-1G	IMIDUREA	1G	39236-46-9
PHR3524-1G	INDOLE	1G	120-72-9
PHR1351-1G	INOSITOL	1G	87-89-8
PHR3188-500MG	IODIPAMIDE	500MG	606-17-7
PHR3015-500MG	IODIXANOL	500MG	92339-11-2
PHR3189-1.2ML	ISOBUTYL ACETATE	1.2ML	110-19-0
PHR3818-250MG	ISOFLUPREDONE ACETATE	250MG	338-98-7
PHR3508-100MG	ISOMALTOTRIOSE, DP3	100MG	3371-50-4
PHR1661-2ML	ISOPENTANE SOLUTION	2ML	78-78-4
PHR2129-50MG	ISOPROPYL 2-(4-ISOBUTYLPHENYL)PROPANOATE	50MG	
PHR1123-1G	ISOPROPYL MYRISTATE	1G	110-27-0
PHR1137-1G	ISOPROPYL PALMITATE	1G	142-91-6
PHR2612-50MG	ISOPROPYL VALPROATE	50MG	60498-67-1
PHR2722-250MG	ISOPROTERENOL HYDROCHLORIDE	250MG	51-30-9
PHR1486-1G	ISOSORBIDE SOLUTION	1G	
PHR1188-3X100MG	ISOTRETINOIN	3X100MG	4759-48-2
PHR3190-500MG	ISOXSUPRINE HYDROCHLORIDE	500MG	579-56-6
PHR1380-1G	IVERMECTIN	1G	70288-86-7
PHR2705-30MG	KAEMPFEROL	30MG	520-18-3
PHR1487-500MG	KANAMYCIN SULFATE	500MG	25389-94-0
PHR1385-1G	KETOCONAZOLE	1G	65277-42-1
PHR1335-1G	LABETALOL HYDROCHLORIDE	1G	32780-64-6
PHR1215-3X1.5ML	LACTIC ACID	3X1.5ML	50-21-5
PHR3192-1G	LACTITOL	1G	81025-04-9
PHR1024-1G	LACTOSE (MONOHYDRATE)	1G	64044-51-5
PHR1025-1G	LACTOSE, ANHYDROUS	1G	63-42-3

Cat. No.	Description	Pack Size	CAS No.
PHR1608-1G	LACTULOSE	1G	4618-18-2
PHR1110-1G	L-ALANINE	1G	56-41-7
PHR2485-300MG	L-ALANYL-L-GLUTAMINE	300MG	39537-23-0
PHR1489-1G	L-ARABINITOL	1G	7643-75-6
PHR1106-1G	L-ARGININE	1G	74-79-3
PHR1104-1G	L-ASPARTIC ACID	1G	56-84-8
PHR1580-1G	LAURIC ACID	1G	143-07-7
PHR2599-1G	LAURYL ALCOHOL	1G	112-53-8
PHR3510-500MG	LAURYL GALLATE	500MG	1166-52-5
PHR3191-500MG	L-CITRULLINE	500MG	372-75-8
PHR1102-1G	L-CYSTEINE HYDROCHLORIDE MONOHYDRATE	1G	7048-04-6
PHR1323-500MG	L-CYSTINE	500MG	56-89-3
PHR1540-1G	LETROZOLE	1G	112809-51-5
PHR1541-1G	LEUCOVORIN CALCIUM	1G	1492-18-8
PHR3193-5MG	LEUPROLIDE ACETATE	5MG	74381-53-6
PHR2630-50MG	LEVOMEDETOMIDINE	50MG	190000-46-5
PHR3196-200MG	LEVOMENOL	200MG	23089-26-1
PHR1613-1G	LEVOTHYROXINE	1G	51-48-9
PHR2880-100MG	LEVOTHYROXINE SODIUM HYDRATE	100MG	25416-65-3
PHR1107-1G	L-GLUTAMIC ACID	1G	56-86-0
PHR1125-1G	L-GLUTAMINE	1G	56-85-9
PHR1108-1G	L-HISTIDINE	1G	71-00-1
PHR3294-50MG	LICARBAZEPINE	50MG	29331-92-8
PHR1657-1G	LINCOMYCIN HYDROCHLORIDE	1G	7179-49-9
PHR2723-5X200MG	LINOLEYL ALCOHOL	5X200MG	506-43-4
PHR1099-1G	L-ISOLEUCINE	1G	73-32-5
PHR1105-1G	L-LEUCINE	1G	61-90-5
PHR1096-1G	L-LYSINE ACETATE	1G	57282-49-2
PHR1101-1G	L-LYSINE MONOHYDROCHLORIDE	1G	657-27-2
PHR1116-1G	L-MENTHOL	1G	2216-51-5
PHR1241-1G	L-METHIONINE	1G	63-68-3
PHR1100-1G	L-PHENYLALANINE	1G	63-91-2
PHR1332-500MG	L-PROLINE	500MG	147-85-3
PHR1103-1G	L-SERINE	1G	56-45-1
PHR1242-500MG	L-THREONINE	500MG	72-19-5
PHR2600-100MG	L-TRYPTOPHAN ETHYL ESTER HYDROCHLORIDE	100MG	2899-28-7
PHR8909-100MG	L-TRYPTOPHAN ISOPROPYL ESTER HYDROCHLORIDE	100MG	71959-71-2
PHR2601-100MG	L-TRYPTOPHAN METHYL ESTER HYDROCHLORIDE	100MG	7524-52-9
PHR1097-1G	L-TYROSINE	1G	60-18-4
PHR1172-1G	L-VALINE	1G	72-18-4
PHR2631-250MG	LYSOPHOSPHATIDYLCHOLINE	250MG	9008-30-4
PHR1328-5G	MAGNESIUM CARBONATE	5G	23389-33-5
PHR2486-1G	MAGNESIUM CHLORIDE	1G	7791-18-6
PHR2604-2G	MAGNESIUM SULFATE	2G	10034-99-8
PHR1272-500MG	MALEIC ACID	500MG	110-16-7
PHR1273-1G	MALIC ACID	1G	6915-15-7
PHR1248-500MG	MALTITOL	500MG	585-88-6
PHR3198-1G	MALTOL	1G	118-71-8
PHR1497-1G	MALTOSE MONOHYDRATE	1G	6363-53-7
PHR3199-1G	MANDELIC ACID	1G	90-64-2
PHR1007-1G	MANNITOL	1G	69-65-8
PHR3033-500MG	MAPROTILINE HYDROCHLORIDE	500MG	10347-81-6
PHR2448-300MG	MEBENDAZOLE	300MG	31431-39-7
PHR3200-500MG	MECLOCYCLINE SULFOSALICYLATE	500MG	73816-42-9

Cat. No.	Description	Pack Size	CAS No.
PHR3201-1G	MECLOFENAMATE SODIUM	1G	6385-02-0
PHR1705-1G	MEFLOQUINE HYDROCHLORIDE	1G	51773-92-3
PHR1324-1G	MEGLUMINE	1G	6284-40-8
PHR3035-250MG	MELENGESTROL ACETATE	250MG	2919-66-6
PHR1411-500MG	MENADIONE	500MG	58-27-5
PHR2363-3X75MG	MENAQUINONE -7	3X75MG	2124-57-4
PHR2271-100MG	MENAQUINONE-4	100MG	863-61-6
PHR1611-1G	MERCAPTOPYRINE	1G	6112-76-1
PHR3426-500MG	METARAMINOL BITARTRATE	500MG	33402-03-8
PHR3735-1G	METHAPYRILENE HCL	1G	135-23-9
PHR3202-1G	METHAZOLAMIDE	1G	554-57-4
PHR1395-1G	METHOCARBAMOL	1G	532-03-6
PHR3203-500MG	METHSCOPOLAMINE BROMIDE	500MG	155-41-9
PHR9168-50MG	METHYL 3-[[(2RS)-2-(BUTYLAMINO)PROPANOYL]AMINO]-4-METHYLTHIOPHENE-2- CARBOXYLATE HYDROCHLORIDE	50MG	23964-59-2
PHR9159-100MG	METHYL 3-AMINO-4-METHYLTHIOPHENE-2-CARBOXYLATE	100MG	85006-31-1
PHR2508-100MG	METHYL 4-AMINO-2-HYDROXYBENZOATE	100MG	
PHR9177-100MG	METHYL 4-AMINO BENZOATE	100MG	619-45-4
PHR2926-3X1.2ML	METHYL ACETATE	3X1.2ML	79-20-9
PHR3539-1G	METHYL ARACHIDATE	1G	1120-28-1
PHR3540-1G	METHYL BEHENATE	1G	929-77-1
PHR3204-3X200MG	METHYL CAPRATE	3X200MG	110-42-9
PHR3535-1G	METHYL CAPRATE	1G	110-42-9
PHR3533-1G	METHYL HEXANOATE	1G	106-70-7
PHR2605-200MG	METHYL LAURATE	200MG	111-82-0
PHR2606-5X100MG	METHYL LINOLEATE	5X100MG	112-63-0
PHR3734-5X100MG	METHYL LINOLENATE	5X100MG	301-00-8
PHR3537-1G	METHYL MYRISTATE	1G	124-10-7
PHR3534-1G	METHYL OCTANOATE	1G	111-11-5
PHR3538-1G	METHYL PALMITATE	1G	112-39-0
PHR3733-5X100MG	METHYL PALMITOLEATE	5X100MG	1120-25-8
PHR3837-200MG	METHYL RICINOLEATE	200MG	141-24-2
PHR3541-100MG	METHYL TETRACOSANOATE	100MG	2442-49-1
PHR3835-250MG	METHYL TRICOSANOATE	250MG	2433-97-8
PHR9099-100MG	METHYL VALPROATE	100MG	22632-59-3
PHR2332-40MG	METHYL 5-[2-[BENZYL(TERT-BUTYL)AMINO]ACETYL]-2-HYDROXYBENZOATE HYDROCHLORIDE	40MG	27475-26-9
PHR3410-2X200MG	METHYLCOBALAMIN	2X200MG	13422-55-4
PHR3838-1G	METHYLENE BLUE	1G	122965-43-9
PHR3575-50MG	METHYLENE SIMVASTATIN	50MG	1449248-72-9
PHR1012-1G	METHYLPARABEN	1G	99-76-3
PHR1012-10G	METHYLPARABEN	10G	99-76-3
PHR1346-1G	METHYLSULFONYLMETHANE	1G	67-71-0
PHR3206-500MG	MEXILETINE HYDROCHLORIDE	500MG	5370-01-4
PHR2633-500MG	MILBEMYCIN OXIME	500MG	129496-10-2
PHR2881-1.5ML	MINERAL OIL	1.5ML	8042-47-5
PHR3207-1G	MITOTANE	1G	53-19-0
PHR2696-200MG	MOGROSIDE V	200MG	88901-36-4
PHR3432-1G	MOLINDONE HYDROCHLORIDE	1G	15622-65-8
PHR1400-500MG	MOMETASONE FUROATE	500MG	83919-23-7
PHR2775-300MG	MONENSIN SODIUM	300MG	22373-78-0
PHR1330-5G	MONOBASIC POTASSIUM PHOSPHATE	5G	7778-77-0
PHR3208-500MG	MONOBENZONE	500MG	103-16-2
PHR2634-1G	MONOSODIUM GLUTAMATE	1G	
PHR2635-200MG	MONOSTEARYL MALEATE	200MG	2424-62-6

Cat. No.	Description	Pack Size	CAS No.
PHR3777-300MG	MONTELUKAST DICYCLOHEXYLAMINE	300MG	577953-88-9
PHR1603-1G	MONTELUKAST SODIUM	1G	151767-02-1
PHR3209-200MG	MORANTEL TARTRATE	200MG	26155-31-7
PHR2806-100MG	MUPIROCIN	100MG	12650-69-0
PHR3210-400MG	MUPIROCIN CALCIUM	400MG	115074-43-6
PHR3421-200MG	MUPIROCIN LITHIUM	200MG	73346-79-9
PHR1313-1G	M-XYLENE	1G	108-38-3
PHR3211-400MG	MYCOPHENOLATE SODIUM	400MG	37415-62-6
PHR1124-1G	MYRISTIC ACID	1G	544-63-8
PHR1135-1G	MYRISTYL ALCOHOL	1G	112-72-1
PHR9147-50MG	N-(2,6-DIMETHYLPHENYL)PYRIDINE-2-CARBOXAMIDE	50MG	39627-98-0
PHR9202-75MG	N,N-DIMETHYL-2-[1-(METHYLSULFONYL)-5-(1H-1,2,4-TRIAZOL-1-YLMETHYL)-1H-INDOL-3-YL]ETHANAMINE	75MG	
PHR9151-50MG	N,N-DIMETHYL-2-[5-(1H-1,2,4-TRIAZOL-1-YLMETHYL)-1H-INDOL-3-YL]ETHANAMINE N-OXIDE	50MG	260435-42-5
PHR3056-500MG	NABUMETONE	500MG	42924-53-8
PHR2636-200MG	N-ACETYL-D,L-METHIONINE	200MG	1115-47-5
PHR1098-1G	N-ACETYL- -CYSTEINE	1G	616-91-1
PHR1177-500MG	N-ACETYL- -TRYPTOPHAN	500MG	1218-34-4
PHR1173-1G	N-ACETYL- -TYROSINE	1G	537-55-3
PHR3212-400MG	NADOLOL	400MG	42200-33-9
PHR1275-1G	NAPHTHALENE	1G	91-20-3
PHR1703-1G	NATAMYCIN	1G	7681-93-8
PHR3213-50MG	NEOHESPERIDIN	50MG	13241-33-3
PHR1491-1G	NEOMYCIN SULFATE	1G	1405-10-3
PHR8908-100MG	N-FORMYLKYNURENINE	100MG	3978-11-8
PHR8748-25MG	N-FORMYL- -LEUCINE (3S,4R,6S)-3-HEXYL-2-OXO-6-UNDECYLTETRAHYDRO-2H-PYRAN-4-YL ESTER	25MG	130793-27-0
PHR3275-3X1.2ML	N-HEXANOL	3X1.2ML	111-27-3
PHR1276-1G	NIACIN (NICOTINIC ACID)	1G	59-67-6
PHR1033-1G	NIACINAMIDE	1G	98-92-0
PHR2531-500MG	NICARDIPINE HYDROCHLORIDE	500MG	54527-84-3
PHR1290-1G	NIFEDIPINE	1G	21829-25-4
PHR3778-50MG	NIFEDIPINE NITROPHENYLPYRIDINE ANALOG	50MG	67035-22-7
PHR3214-500MG	NILUTAMIDE	500MG	63612-50-0
PHR1336-1G	NIMESULIDE	1G	51803-78-2
PHR1191-1G	NITROFURANTOIN	1G	67-20-9
PHR2456-200MG	NITROFURFURAL DIACETATE	200MG	92-55-7
PHR3215-500MG	NIZATIDINE	500MG	76963-41-2
PHR9100-100MG	N-METHYL-5-(METHYLAMINO)-1H-IMIDAZOLE-4-CARBOXAMIDE	100MG	6736-40-9
PHR1352-2G	N-METHYLPYRROLIDONE	2G	872-50-4
PHR9227-25MG	N-NITROSO-BISOPROLOL	25MG	2820170-76-9
PHR3608-1ML	N-NITROSODIBUTYLAMINE (NDBA)	1ML	924-16-3
PHR2408-1G	N-NITROSODIETHYLAMINE (NDEA)	1G	55-18-5
PHR3607-1ML	N-NITROSODIISOPROPYLAMINE (NDIPA)	1ML	601-77-4
PHR3609-1ML	N-NITROSOETHYLISOPROPYLAMINE (NEIPA / NIPEA)	1ML	16339-04-1
PHR3610-1ML	N-NITROSOMETHYLAMINO BUTYRIC ACID (NMBA)	1ML	61445-55-4
PHR3611-1ML	N-NITROSOMETHYLPHENYLAMINE (NMPA)	1ML	614-00-6
PHR2617-100MG	N-NORVERAPAMIL	100MG	67812-42-4
PHR2731-1G	NONOXYNOL 9	1G	
PHR3509-500MG	NORDIHYDROGUAIARETIC ACID	500MG	500-38-9
PHR3787-400MG	NORELGESTROMIN	400MG	53016-31-2
PHR3820-250MG	NOREPINEPHRINE BITARTRATE	250MG	108341-18-0
PHR1536-200MG	NORPHENYLEPHRINE HYDROCHLORIDE	200MG	4779-94-6
PHR8915-30MG	NORSCOPOLAMINE	30MG	4684-28-0

Cat. No.	Description	Pack Size	CAS No.
PHR3216-500MG	NORTRIPTYLINE HYDROCHLORIDE	500MG	894-71-3
PHR1080-1G	OCTINOXATE	1G	5466-77-3
PHR1080-10G	OCTINOXATE	10G	5466-77-3
PHR1081-10G	OCTISALATE	10G	118-60-5
PHR1083-1G	OCTOCRYLENE	1G	6197-30-4
PHR1083-10G	OCTOCRYLENE	10G	6197-30-4
PHR3511-500MG	OCTYL GALLATE	500MG	1034-01-1
PHR1155-1G	OCTYLDODECANOL	1G	5333-42-6
PHR3764-50MG	O-FLUTAMIDE	50MG	151262-93-0
PHR1586-1G	OLEIC ACID	1G	112-80-1
PHR2637-5X100MG	OLEYL ALCOHOL	5X100MG	143-28-2
PHR2902-5X400MG	OLIVE OIL	5X400MG	8001-25-0
PHR1473-500MG	O-TOLUENESULFONAMIDE	500MG	88-19-7
PHR3291-500MG	OXALIC ACID	500MG	144-62-7
PHR3217-500MG	OXAPROZIN	500MG	21256-18-8
PHR1074-1G	OXYBENZONE	1G	131-57-7
PHR1312-1G	O-XYLENE	1G	95-47-6
PHR1537-1G	OXYTETRACYCLINE	1G	6153-64-6
PHR1325-500MG	OXYTETRACYCLINE HYDROCHLORIDE	500MG	2058-46-0
PHR1355-10G	PADIMATE O	10G	21245-02-3
PHR1120-1G	PALMITIC ACID	1G	57-10-3
PHR2552-200MG	PALONOSETRON HYDROCHLORIDE	200MG	135729-62-3
PHR1182-500MG	PAPAVERINE HYDROCHLORIDE	500MG	61-25-6
PHR2596-1G	PARACHLOROPHENOL	1G	106-48-9
PHR2985-1G	PEANUT OIL	1G	8002-03-7
PHR1596-500MG	PEMETREXED	500MG	357166-29-1
PHR3675-50MG	PEMETREXED ENANTIOMER	50MG	182009-04-7
PHR2638-500MG	PENICILLAMINE	500MG	52-67-5
PHR2639-500MG	PENICILLAMINE DISULFIDE	500MG	20902-45-8
PHR3156-50MG	PENTANAMIDE (VALERAMIDE)	50MG	626-97-1
PHR3157-50MG	PENTANOIC ACID (VALERIC ACID)	50MG	109-52-4
PHR2645-200MG	PENTETIC ACID	200MG	67-43-6
PHR1249-1G	PHENACETIN	1G	62-44-2
PHR1094-1G	PHENACETIN MELTING POINT STANDARD	1G	62-44-2
PHR3220-500MG	PHENAZOPYRIDINE HYDROCHLORIDE	500MG	136-40-3
PHR1573-500MG	PHENFORMIN HYDROCHLORIDE	500MG	834-28-6
PHR1047-1G	PHENOL	1G	108-95-2
PHR1402-1G	PHENOXYBENZAMINE HYDROCHLORIDE	1G	63-92-3
PHR1403-1G	PHENTOLAMINE MESYLATE	1G	65-28-1
PHR1152-1G	PHENYL SALICYLATE MELTING POINT STANDARD	1G	118-55-8
PHR2607-500MG	PHENYLBUTAZONE	500MG	50-33-9
PHR3221-50MG	PHENYLEPHRINE BITARTRATE	50MG	17162-39-9
PHR1122-1.5G	PHENYLETHYL ALCOHOL	1.5G	60-12-8
PHR1122-2G	PHENYLETHYL ALCOHOL	2G	60-12-8
PHR3083-200MG	PHENYLTOLOXAMINE CITRATE	200MG	1176-08-5
PHR2649-250MG	PHOSPHATIDYLCHOLINE (SOY)	250MG	97281-47-5
PHR3222-5X1.5ML	PHOSPHORIC ACID	5X1.5ML	7664-38-2
PHR1078-1G	PHYTONADIONE	1G	84-80-0
PHR1493-500MG	PILOCARPINE HYDROCHLORIDE	500MG	54-71-7
PHR1494-500MG	PILOCARPINE NITRATE	500MG	148-72-1
PHR3084-200MG	PIMOBENDAN	200MG	74150-27-9
PHR3223-500MG	PIMOZIDE	500MG	2062-78-4
PHR3224-500MG	PINDOLOL	500MG	13523-86-9
PHR2211-50MG	PIPERINE	50MG	94-62-2

Cat. No.	Description	Pack Size	CAS No.
PHR2675-500MG	PIRACETAM	500MG	7491-74-9
PHR1595-500MG	POLYMYXIN B SULFATE	500MG	1405-20-5
PHR3789-5X1G	POLYSORBATE 20	5X1G	9005-64-5
PHR3797-5X1G	POLYSORBATE 60	5X1G	9005-67-8
PHR3790-5X1G	POLYSORBATE 80	5X1G	9005-65-6
PHR1329-5G	POTASSIUM CHLORIDE	5G	7447-40-7
PHR2214-1G	POTASSIUM CYANIDE	1G	151-50-8
PHR1130-1G	POTASSIUM GLUCONATE	1G	299-27-4
PHR1360-5G	POTASSIUM IODIDE	5G	7681-11-0
PHR1278-1G	POTASSIUM SORBATE	1G	24634-61-5
PHR2792-500MG	POTASSIUM SUCROSE OCTASULFATE	500MG	73264-44-5
PHR1250-500MG	POVIDONE	500MG	9003-39-8
PHR3225-500MG	PRALIDOXIME CHLORIDE	500MG	14018-50-9
PHR1302-1G	PRAMOXINE HYDROCHLORIDE	1G	637-58-1
PHR3226-300MG	PRASUGREL HYDROCHLORIDE	300MG	389574-19-0
PHR2564-100MG	PREGNENOLONE	100MG	145-13-1
PHR3108-500MG	PRIMAQUINE PHOSPHATE	500MG	63-45-6
PHR2608-500MG	PROBENECID	500MG	57-66-9
PHR1252-500MG	PROCAINAMIDE HYDROCHLORIDE	500MG	614-39-1
PHR3228-500MG	PROMAZINE HYDROCHLORIDE	500MG	53-60-1
PHR3229-500MG	PROPARACAINE HYDROCHLORIDE	500MG	5875-06-9
PHR1051-1.5G	PROPYLENE GLYCOL	1.5G	57-55-6
PHR1010-10G	PROPYLPARABEN	10G	94-13-3
PHR2609-500MG	PROPYLTHIOURACIL	500MG	51-52-5
PHR1389-1G	PROPYPHENAZONE	1G	479-92-5
PHR2610-25MG	PROTAMINE SULFATE	25MG	9009-65-8
PHR1474-500MG	P-TOLUENESULFONAMIDE	500MG	70-55-3
PHR1314-1G	P-XYLENE	1G	106-42-3
PHR1576-500MG	PYRAZINAMIDE	500MG	98-96-4
PHR3120-500MG	PYRIDOSTIGMINE BROMIDE	500MG	101-26-8
PHR1340-500MG	PYRILAMINE MALEATE	500MG	59-33-6
PHR3231-500MG	PYRIMETHAMINE	500MG	58-14-0
PHR2817-100MG	PYROGLUTAMIC ACID	100MG	98-79-3
PHR1488-1G	QUERCETIN	1G	849061-97-8
PHR1468-1G	QUINIC ACID	1G	77-95-2
PHR3232-500MG	QUINIDINE GLUCONATE	500MG	7054-25-3
PHR3233-1G	QUINIDINE SULFATE	1G	6591-63-5
PHR2653-500MG	RACEMETHIONINE	500MG	59-51-8
PHR3234-50MG	RACEPINEPHRINE HYDROCHLORIDE	50MG	329-63-5
PHR3441-400MG	RACTOPAMINE HYDROCHLORIDE	400MG	90274-24-1
PHR2697-1G	REBAUDIOSIDE A	1G	58543-16-1
PHR1201-500MG	RESORCINOL	500MG	108-46-3
PHR3526-500MG	RIBITOL	500MG	
PHR2611-100MG	RIFABUTIN	100MG	72559-06-9
PHR1664-500MG	RIFAMYCIN SODIUM	500MG	14897-39-3
PHR3267-500MG	RIMANTADINE HYDROCHLORIDE	500MG	1501-84-4
PHR1631-500MG	RISPERIDONE	500MG	106266-06-2
PHR1615-1G	RITODRINE HYDROCHLORIDE	1G	23239-51-2
PHR2397-500MG	ROCURONIUM BROMIDE	500MG	119302-91-9
PHR2932-300MG	ROSIGLITAZONE	300MG	122320-73-4
PHR2513-300MG	ROSIGLITAZONE MALEATE	300MG	155141-29-0
PHR3235-50MG	ROSMARINIC ACID	50MG	20283-92-5
PHR3236-500MG	ROXARSONE	500MG	121-19-7
PHR1341-1G	SACCHARIN	1G	81-07-2

Cat. No.	Description	Pack Size	CAS No.
PHR1348-500MG	SACCHARIN SODIUM	500MG	128-44-9
PHR2987-1G	SAFFLOWER OIL	1G	8001-23-8
PHR3782-50MG	SALIDROSIDE	50MG	10338-51-9
PHR1470-500MG	SCOPOLAMINE HYDROBROMIDE	500MG	6533-68-2
PHR3134-500MG	SELEGILINE HYDROCHLORIDE	500MG	14611-52-0
PHR2613-200MG	SELENOMETHIONINE	200MG	1464-42-2
PHR2654-500MG	SILVER SULFADIAZINE	500MG	22199-08-2
PHR2655-100MG	SILYBIN	100MG	22888-70-6
PHR1438-1G	SIMVASTATIN	1G	79902-63-9
PHR2615-2G	SODIUM ACETATE	2G	127-09-3
PHR2618-2G	SODIUM ACETATE TRIHYDRATE	2G	6131-90-4
PHR1279-1G	SODIUM ASCORBATE	1G	134-03-2
PHR1231-1G	SODIUM BENZOATE	1G	532-32-1
PHR1321-5G	SODIUM CHLORIDE	5G	7647-14-5
PHR1416-1G	SODIUM CITRATE	1G	6132-04-3
PHR3740-500MG	SODIUM CYCLAMATE	500MG	139-05-9
PHR1609-10G	SODIUM DEHYDROACETIC ACID	10G	4418-26-2
PHR1305-1G	SODIUM DODECYLBENZENESULFONATE	1G	25155-30-0
PHR2657-200MG	SODIUM -PANTOATE	200MG	60979-68-2
PHR1408-1G	SODIUM FLUORIDE	1G	7681-49-4
PHR1113-1G	SODIUM -LACTATE	1G	867-56-1
PHR1434-1G	SODIUM METABISULFITE	1G	7681-57-4
PHR2787-1G	SODIUM MONOFLUOROPHOSPHATE	1G	10163-15-2
PHR1423-1G	SODIUM NITROPRUSSIDE	1G	13755-38-9
PHR2540-100MG	SODIUM PICOSULFATE	100MG	10040-45-6
PHR3238-1G	SODIUM SALICYLATE	1G	54-21-7
PHR2658-1G	SODIUM SULFATE ANHYDROUS	1G	
PHR1409-4X1G	SODIUM TARTRATE DIHYDRATE	4X1G	6106-24-7
PHR2690-2G	SODIUM THIOSULFATE	2G	7772-98-7
PHR2204-1G	SODIUM TRIMETAPHOSPHATE	1G	7785-84-4
PHR1367-1G	SORBIC ACID	1G	110-44-1
PHR1006-1G	SORBITOL	1G	50-70-4
PHR2882-1G	SOYBEAN OIL	1G	8001-22-7
PHR1426-500MG	SPECTINOMYCIN HYDROCHLORIDE	500MG	
PHR1441-1G	SPECTINOMYCIN SULFATE	1G	64058-48-6
PHR1417-1G	SQUALANE	1G	111-01-3
PHR3519-1G	SQUALENE	1G	111-02-4
PHR3398-200MG	S-SELENOMETHIONINE	200MG	3211-76-5
PHR2988-1G	STAR ANISE OIL	1G	68952-43-2
PHR3447-500MG	STAVUDINE	500MG	3056-17-5
PHR1114-1G	STEARIC ACID	1G	57-11-4
PHR1115-1G	STEARYL ALCOHOL	1G	112-92-5
PHR2687-500MG	STREPTOMYCIN SULFATE	500MG	3810-74-0
PHR3500-1ML	STYRENE	1ML	100-42-5
PHR1418-1G	SUCCINIC ACID	1G	110-15-6
PHR2659-1G	SUCCINYLCHOLINE CHLORIDE	1G	6101-15-1
PHR2660-250MG	SUCCINYLMONCHOLINE CHLORIDE	250MG	5297-17-6
PHR1342-500MG	SUCRALOSE	500MG	56038-13-2
PHR1519-1G	SULFACETAMIDE SODIUM	1G	6209-17-2
PHR3240-500MG	SULFACHLORPYRIDAZINE	500MG	80-32-0
PHR1343-500MG	SULFADIAZINE	500MG	68-35-9
PHR1448-500MG	SULFADIMETHOXINE	500MG	122-11-2
PHR3241-500MG	SULFADOXINE	500MG	2447-57-6
PHR3136-1G	SULFAMETHAZINE	1G	57-68-1

Cat. No.	Description	Pack Size	CAS No.
PHR3242-500MG	SULFAMETHIZOLE	500MG	144-82-1
PHR1184-5G	SULFANILAMIDE	5G	
PHR1093-1G	SULFANILAMIDE MELTING POINT STANDARD	1G	63-74-1
PHR1419-500MG	SULFANILIC ACID	500MG	121-57-3
PHR1087-1G	SULFAPYRIDINE MELTING POINT STANDARD	1G	144-83-2
PHR3243-250MG	SULFASALAZINE	250MG	599-79-1
PHR3244-500MG	SULFINPYRAZONE	500MG	57-96-5
PHR3359-500MG	SULINDAC	500MG	38194-50-2
PHR2989-1G	SUNFLOWER OIL	1G	8001-21-6
PHR3245-500MG	TAGATOSE	500MG	87-81-0
PHR2706-300MG	TAMOXIFEN CITRATE	300MG	54965-24-1
PHR1524-200MG	TAMSULOSIN HYDROCHLORIDE	200MG	106463-17-6
PHR1667-3X1.2ML	T-AMYL ALCOHOL	3X1.2ML	75-85-4
PHR1472-1G	TARTARIC ACID	1G	87-69-4
PHR1109-1G	TAURINE	1G	107-35-7
PHR1686-1G	TAZOBACTAM	1G	89786-04-9
PHR1437-1G	TEMOZOLOMIDE	1G	85622-93-1
PHR3247-500MG	TERCONAZOLE	500MG	
PHR3248-1G	TERPIN HYDRATE	1G	2451-01-6
PHR3512-500MG	TERT-BUTYLHYDROQUINONE	500MG	1948-33-0
PHR3499-1ML	TETRACHLOROETHYLENE	1ML	127-18-4
PHR1041-500MG	TETRACYCLINE HYDROCHLORIDE	500MG	64-75-5
PHR3249-500MG	THALIDOMIDE	500MG	50-35-1
PHR1259-500MG	THIABENDAZOLE	500MG	148-79-8
PHR2265-250MG	THIO-ACID PROPIONATE	250MG	
PHR3250-500MG	THIOGUANINE	500MG	154-42-7
PHR3143-500MG	THIORIDAZINE HYDROCHLORIDE	500MG	130-61-0
PHR1758-200MG	THIOUREA	200MG	62-56-6
PHR1134-1G	THYMOL	1G	89-83-8
PHR2814-200MG	TIBOLONE	200MG	5630-53-5
PHR2788-400MG	TICAGRELOR	400MG	274693-27-5
PHR3367-300MG	TICLOPIDINE HYDROCHLORIDE	300MG	53885-35-1
PHR1079-1G	TOBRAMYCIN	1G	32986-56-4
PHR3251-500MG	TOLBUTAMIDE	500MG	64-77-7
PHR1311-1G	TOLUENE	1G	108-88-3
PHR3493-1ML	TRANS-1,2-DICHLOROETHYLENE	1ML	156-60-5
PHR1218-2.2ML	TRANS-ANETHOLE	2.2ML	4180-23-8
PHR3252-50MG	TRANSPLATIN	50MG	14913-33-8
PHR1344-500MG	TREHALOSE	500MG	6138-23-4
PHR1187-3X100MG	TRETINOIN	3X100MG	302-79-4
PHR1410-1G	TRIACETIN	1G	102-76-1
PHR1156-1G	TRIBUTYL 2-ACETYLCITRATE	1G	77-90-7
PHR1205-0.5ML	TRIBUTYL PHOSPHATE	0.5ML	126-73-8
PHR1562-3X1.2ML	TRICHLOROETHYLENE	3X1.2ML	79-01-6
PHR1338-1G	TRICLOSAN	1G	3380-34-5
PHR1157-1G	TRIETHYL 2-ACETYLCITRATE	1G	77-89-4
PHR1153-1G	TRIETHYL CITRATE	1G	77-93-0
PHR2883-5X2ML	TRIETHYLAMINE	5X2ML	121-44-8
PHR3826-1ML	TRIETHYLENE GLYCOL	1ML	112-27-6
PHR3150-200MG	TRIFLURIDINE	200MG	70-00-8
PHR2661-100MG	TRIGLYCINE	100MG	556-33-2
PHR2804-50MG	TRIMEBUTINE MALEATE	50MG	34140-59-5
PHR3257-500MG	TRIOXSALEN	500MG	3902-71-4
PHR3258-500MG	TRIPLENNAMINE HYDROCHLORIDE	500MG	154-69-8

Cat. No.	Description	Pack Size	CAS No.
PHR2662-500MG	TRIPROLIDINE	500MG	6138-79-0
PHR2928-100MG	TRIPROLIDINE HYDROCHLORIDE Z-ISOMER	100MG	51657-91-1
PHR1347-1G	TROMETHAMINE	1G	77-86-1
PHR2815-50MG	TROXERUTIN	50MG	7085-55-4
PHR2652-500MG	TYLOSIN	500MG	1401-69-0
PHR1508-500MG	TYLOSIN TARTRATE	500MG	1405-54-5
PHR1429-1G	UNDECYLENOYL PHENYLALANINE (SEPIWHITE)	1G	
PHR1581-100MG	URACIL	100MG	66-22-8
PHR3260-100MG	URACIL ARABINOSIDE	100MG	3083-77-0
PHR1406-1G	UREA	1G	57-13-6
PHR1579-500MG	URSODIOL	500MG	128-13-2
PHR1626-1G	VALGANCICLOVIR HYDROCHLORIDE	1G	175865-59-5
PHR9137-100MG	VALPROMIDE	100MG	2430-27-5
PHR1245-1G	VANILLIN	1G	121-33-5
PHR1091-1G	VANILLIN MELTING POINT STANDARD	1G	121-33-5
PHR2720-500MG	VARDENAFIL HYDROCHLORIDE	500MG	330808-88-3
PHR1627-100MG	VECURONIUM BROMIDE	100MG	50700-72-6
PHR2698-100MG	VINBLASTINE SULFATE	100MG	143-67-9
PHR3160-500MG	VINORELBINE TARTRATE	500MG	125317-39-7
PHR3501-1ML	VINYL ACETATE	1ML	108-05-4
PHR2692-50MG	VITEXIN	50MG	3681-93-4
PHR3261-200MG	XANTHANOIC ACID	200MG	82-07-5
PHR3262-200MG	XANTHONE	200MG	90-47-1
PHR3263-500MG	XYLAZINE	500MG	7361-61-7
PHR3264-500MG	XYLAZINE HYDROCHLORIDE	500MG	23076-35-9
PHR1280-1G	XYLITOL	1G	87-99-0
PHR3265-500MG	YOHIMBINE HYDROCHLORIDE	500MG	65-19-0
PHR3266-500MG	ZALCITABINE	500MG	7481-89-2
PHR1401-5G	ZINC PYRITHIONE	5G	
PHR2917-1G	ZINC SULFATE	1G	7446-19-7
PHR9043-50MG	A,A-DIPHENYL-4-PIPERIDINEMETHANOL	50MG	115-46-8
PHR9209-50MG	A,A-DIPHENYL-Γ-BUTYROLACTONE	50MG	956-89-8
PHR1239-1G	B-CAROTENE	1G	7235-40-7
PHR1239-3X100MG	B-CAROTENE	3X100MG	7235-40-7
PHR2693-200MG	B-SITOSTEROL	200MG	83-46-5


Residual Solvents

- Residual solvents are defined as Organic Volatile Chemicals
- Used or produced in Drugs/Excipients manufacturing
- Trace amounts can impact Drug Safety & Efficacy
- Evaluation is quite crucial due to possible risk to human health
- Divided in 3 categories (Class-I, Class-II & Class-III)
- Strict guidelines to be followed by all Pharma companies



Cat. No.	Description	Pack Size	CAS No.
Residual Solvents			
PHR1549-3X1.2ML	1,4-DIOXANE	3X1.2ML	123-91-1
RTC000093-20ML	1,4-DIOXANE	20ML	123-91-1
PHR1209-3X1.2ML	1-BUTANOL	3X1.2ML	71-36-3
PHR1217-3X1.2ML	1-PENTANOL	3X1.2ML	71-41-0
PHR1208-3X1.2ML	1-PROPANOL	3X1.2ML	71-23-8
PHR1563-3X1.2ML	2-BUTANOL	3X1.2ML	78-92-2
PHR1550-3X1.2ML	2-ETHOXYETHANOL	3X1.2ML	110-80-5
PHR1544-3X1.2ML	2-METHOXY ETHANOL	3X1.2ML	109-86-4
PHR2729-5X1.2ML	2-METHYL-1-PROPANOL	5X1.2ML	78-83-1
PHR1072-3X1.2ML	2-PROPANOL	3X1.2ML	67-63-0
RTC000079-20ML	2-PROPANOL	20ML	67-63-0
PHR2337-5X1.2ML	3-METHYL-1-BUTANOL	5X1.2ML	123-51-3
PHR2107-3X1.2ML	ACETALDEHYDE	3X1.2ML	75-07-0
RTC000075-20ML	ACETALDEHYDE	20ML	75-07-0
PHR1219-3X1.2ML	ACETONE	3X1.2ML	67-64-1
RTC000078-20ML	ACETONE	20ML	67-64-1
PHR1551-3X1.2ML	ACETONITRILE	3X1.2ML	75-05-8
RTC000080-20ML	ACETONITRILE	20ML	75-05-8
PHR1212-3X1.2ML	ANISOLE	3X1.2ML	100-66-3
PHR1211-3X1.2ML	BUTYL ACETATE	3X1.2ML	123-86-4
PHR1665-3X1.2ML	CHLOROBENZENE	3X1.2ML	108-90-7
RTC000099-20ML	CHLOROBENZENE	20ML	108-90-7
PHR1552-3X1.2ML	CHLOROFORM	3X1.2ML	67-66-3
PHR1210-3X1.2ML	CUMENE	3X1.2ML	98-82-8
PHR1687-3X1.2ML	CYCLOHEXANE	3X1.2ML	110-82-7
RTC000089-20ML	CYCLOHEXANE	20ML	110-82-7
RTC000077-20ML	DIETHYL ETHER	20ML	60-29-7
PHR2183-3X1.2ML	DIISOPROPYL ETHER	3X1.2ML	108-20-3
RTC000085-20ML	DIISOPROPYL ETHER	20ML	108-20-3
RTC000103-20ML	DIMETHYL SULFOXIDE	20ML	67-68-5
PHR1553-3X1.2ML	DIMETHYLFORMAMIDE	3X1.2ML	68-12-2
RTC000098-20ML	DIMETHYLFORMAMIDE	20ML	68-12-2
PHR1481-3X1.2ML	ETHYL ACETATE	3X1.2ML	141-78-6
RTC000087-20ML	ETHYL ACETATE	20ML	141-78-6

Cat. No.	Description	Pack Size	CAS No.
PHR1229-3X1.2ML	ETHYLBENZENE	3X1.2ML	100-41-4
RTC000100-20ML	ETHYLBENZENE	20ML	100-41-4
PHR1554-3X1.2ML	HEPTANE	3X1.2ML	142-82-5
RTC000091-20ML	HEPTANE	20ML	142-82-5
PHR1555-3X1.2ML	HEXANE	3X1.2ML	110-54-3
RTC000084-20ML	HEXANE	20ML	110-54-3
PHR1915-3X1.2ML	ISOOCTANE	3X1.2ML	540-84-1
PHR2109-3X1.2ML	ISOPROPYL ACETATE	3X1.2ML	108-21-4
RTC000090-20ML	ISOPROPYL ACETATE	20ML	108-21-4
PHR1547-3X1.2ML	MESITYL OXIDE	3X1.2ML	141-79-7
RTC000097-20ML	MESITYL OXIDE	20ML	141-79-7
PHR1372-3X1.5ML	METHANOL	3X1.5ML	67-56-1
RTC000076-20ML	METHANOL	20ML	67-56-1
PHR1556-3X1.2ML	METHYL ISOBUTYL KETONE	3X1.2ML	108-10-1
RTC000095-20ML	METHYL ISOBUTYL KETONE	20ML	108-10-1
PHR1557-3X1.2ML	METHYLENECHLORIDE	3X1.2ML	75-09-2
RTC000101-20ML	M-XYLENE	20ML	108-38-3
PHR2110-3X1.2ML	N,N-DIMETHYLACETAMIDE	3X1.2ML	127-19-5
RTC000104-20ML	N,N-DIMETHYLACETAMIDE	20ML	127-19-5
RTC000102-20ML	O-XYLENE	20ML	95-47-6
RTC000094-20ML	PROPYL ACETATE	20ML	109-60-4
PHR1558-3X1.2ML	PYRIDINE	3X1.2ML	110-86-1
PHR1320-3X1.2ML	RESIDUAL SOLVENT - ACETONITRILE	3X1.2ML	75-05-8"
PHR1318-3X1.2ML	RESIDUAL SOLVENT - BENZENE	3X1.2ML	—
PHR1316-3X1.2ML	RESIDUAL SOLVENT - CHLOROFORM	3X1.2ML	67-66-3"
PHR1170-3X1.2ML	RESIDUAL SOLVENT - METHANOL	3X1.2ML	—
PHR1171-3X1.2ML	RESIDUAL SOLVENT - METHYLENE CHLORIDE	3X1.2ML	—
PHR1063-3X1.2ML	RESIDUAL SOLVENTS MIXTURE - CLASS I	3X1.2ML	—
PHR1064-3X1.2ML	RESIDUAL SOLVENTS MIXTURE - CLASS IIA	3X1.2ML	—
PHR1065-3X1.2ML	RESIDUAL SOLVENTS MIXTURE - CLASS IIB	3X1.2ML	—
PHR1066-3X1.2ML	RESIDUAL SOLVENTS MIXTURE - CLASS IIC	3X1.2ML	—
PHR2184-3X1.2ML	TERT-BUTANOL	3X1.2ML	75-65-0
RTC000082-20ML	TERT-BUTANOL	20ML	75-65-0
PHR1560-3X1.2ML	TERT-BUTYL METHYL ETHER	3X1.2ML	1634-04-4
RTC000083-20ML	TERT-BUTYL METHYL ETHER	20ML	1634-04-4
PHR1561-3X1.2ML	TETRAHYDROFURAN	3X1.2ML	109-99-9
RTC000088-20ML	TETRAHYDROFURAN	20ML	109-99-9
PHR1317-3X1.2ML	TOLUENE	3X1.2ML	108-88-3
RTC000096-20ML	TOLUENE	20ML	108-88-3



To know more contact:

Merck Life Science Private Limited

8th Floor, Godrej One, Pirojsha Nagar, Eastern Express Highway, Vikhroli (E), Mumbai – 400079

Website : www.sigmaaldrich.com